
# Spinocan® / Pencan®

---

NRFit®

---

**B | BRAUN**

LLDorder 5701 - grunewald 211525

<b>GB</b>	Spinal Needle for spinal and diagnostic puncture	<b>8</b>
<b>DE</b>	Spinalkanüle für Spinal- und diagnostische Punktion	<b>10</b>
<b>BG</b>	Спинална игла за спинална и диагностична пункция	<b>12</b>
<b>CZ</b>	Mišní jehla pro mišní a diagnostickou punkci	<b>14</b>
<b>DK</b>	Lumbalkanyle til lumbalpunktur og diagnostisk punktur	<b>16</b>
<b>EE</b>	Spinaalnõel seljaaju- ja diagnostilise punktsiooni jaoks	<b>17</b>
<b>ES</b>	Aguja espinal para punción espinal y de diagnóstico	<b>19</b>
<b>FI</b>	Selkäydinneula selkäydin- ja diagnostiikkapistoihin	<b>21</b>
<b>FR</b>	Aiguille spinale pour ponction spinale et de diagnostic	<b>23</b>
<b>GR</b>	Βελόνα οσφυονωτιαίας παρακέντησης για οσφυονωτιαία και διαγνωστική παρακέντηση	<b>25</b>
<b>HR</b>	Spinalna igla za spinalnu i dijagnostičku punkciju	<b>27</b>
<b>HU</b>	Spinális tű spinális és diagnosztikai punkcióhoz	<b>29</b>
<b>IT</b>	Ago spinale per puntura spinale e diagnostica	<b>31</b>
<b>LT</b>	Spinalinė adata, skirta atlikti spinalinę bei diagnostinę punkciją	<b>33</b>
<b>LV</b>	Spinālā adata spinālās un diagnostiskās punkcijas veikšanai	<b>35</b>
<b>NL</b>	Spinale naald voor spinale en diagnostische punctie	<b>37</b>
<b>NO</b>	Spinalnål for spinal- og diagnostisk punksjon	<b>39</b>
<b>PL</b>	Igła do nakłuc rdzeniowych oraz diagnostycznych	<b>41</b>
<b>PT</b>	Aguilha Espinal para punção lombar e de diagnóstico	<b>43</b>
<b>RO</b>	Ac spinal pentru punctie spinală și de diagnostic	<b>45</b>
<b>SE</b>	Spinalnål för spinalpunktion och diagnostisk punktion	<b>47</b>

**SI** Spinalna igla za spinalno in diagnostično punkcijo

**49**

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### **Content:**

Spinal needle with NRFit® connector (either Spinocan® or Pencan®/ Pencan® Paed)

Optional: Spinal Introducer needle (depending on size of Spinal needle)

### **Important!**

Spinal needles with NRFit® connector are not compatible with Luer devices. They are only suitable for products with NRFit® connector!

### **Materials used**

PC, PE, PP, stainless steel, UV-adhesive

### **Indications**

- Spinal anesthesia (Injection of local anesthetics into the subarachnoid space for diagnosis, operative procedures and pain management).
- Lumbar puncture (diagnosis)

### **Contraindications**

The recognized contraindications for spinal anesthesia and lumbar puncture must be observed.

Be aware of:

- blood coagulation disorders
- skin infection at or near the puncture / injection site
- sepsis
- History of hypersensitivity to local anaesthetic or any of the materials employed is known
- Patient refusal

- Severe untreated hypovolaemia
- Increased intracranial pressure

Spinal anesthesia should not be performed in patients with the following conditions:

- severe decompensated hypovolemia
- shock
- acute cerebral or spinal cord disease
- increased intracranial pressure
- existing infection at the injection site
- endogenous or iatrogenic coagulopathy
- Anatomical deformities of the patient's back
- Amniocentesis

### **Risks**

The known risks of spinal anesthesia and lumbar puncture are especially:

- Post dural puncture headache
- Hypotension and bradycardia
- Cardiac arrest
- Hypothermia
- Nausea and vomiting
- Total spinal anaesthesia
- Low cerebrospinal fluid pressure syndrome
- Apnoea
- Acute toxicity to the local anesthetic

Neurological disorders such as a consequence of:

- Epidural haematoma or abscess formation
- Anterior spinal artery syndrome

- Cauda equina syndrome
- Urinary retention
- Herniation
- Pain
- Infections
- Mispositioning

Spinal anesthesia and lumbar puncture should only be performed with the appropriate equipment and trained staff to manage such complications.

Very rare complications:

- hearing loss
- epidermoid tumor
- retroperitoneal abscess

### **Caution**

Products with NRFit® connector are not compatible with Luer devices. They are only suitable for products with NRFit® connector!

Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient. Do not use if packaging is damaged. Do not resterilize. Avoid applying excessive force during advancement as the spinal needle might bend and eventually break. The tip of the spinal needle may become blunt or deform after bone contact. A needle with a damaged tip increases the risk of post dural puncture headache due to extended perforation of the dura mater.

### **Duration of use**

Spinal needles are used for a few minutes for diagnostic punc-

ture (lumbar puncture) or for the injection of drugs into the subarachnoid space (spinal anesthesia) and are withdrawn from the patient immediately after completion of the procedure.

#### Instructions for use

When performing spinal anesthesia /Lumbar puncture, aseptic conditions are required.

1. Position the patient according to his/her condition and/or the planned puncture/injection site. Choose between lateral decubitus or sitting (prone) position. Locate the midline and follow it during needle placement.
2. Select a suitable spinal needle. It is recommended to use an introducer needle for spinal needles that are 25G (0.5 mm) or less in diameter. Check emergency equipment (intubation kit, drugs) and ensure continuous monitoring of the patient.
3. Disinfect puncture site.
4. Inject 1-1.5 cc of local anesthetic into the skin at the selected puncture site of the spinal needle to anesthetize the tissue and the supraspinous and intraspinal ligaments.
5. The introducer needle (for needles 25G or smaller) is inserted at the puncture site until reaching the interspinous ligament. The introducer is fixed with thumb and index finger of the left hand

(for right-hand user).

6. The spinal needle including the stylet is advanced through the introducer needle. The typical "dural click" occurs when the ligamentum flavum has been passed and the subarachnoid space has been reached.

#### Caution

**If any resistance is felt during advancement of the spinal needle carefully correct the orientation of the needle but never apply strong forces in order to overcome the obstacles.**

7. Remove the stylet from the spinal needle. Carefully check for the occurrence of cerebrospinal fluid (CSF) which may be slightly contaminated by blood but must never be pure blood.
  - Lumbar puncture procedure: take CSF for diagnosis, withdraw spinal needle, cover the puncture site with sterile plaster/dressing.
  - Spinal anesthesia: before injection of the local anesthetic, withdraw CSF from the subarachnoid space to ensure proper needle placement. Inject the drug according to the manufacturers instructions.

#### Note

- **The local anesthetic must never be injected without the evidence of cerebrospinal fluid. If pure blood is flowing instead of cerebro-**

**spinal fluid withdraw the needle and retry at a different intervertebral space.**

- **If no cerebrospinal fluid is flowing back rotate the needle to all four quadrants and carefully aspirate until cerebrospinal fluid is visible. If this procedure does not show any backflow of cerebrospinal fluid, repeat the spinal puncture using a different needle direction.**
- **If paresthesia occurs after reaching the subarachnoid space, the needle must be slightly withdrawn.**
- **If paresthesia occurs during injection, the needle must be repositioned before continuing the injection.**
- 8. After complete drug injection withdraw the needle and cover the puncture site with sterile plaster/dressing.
- 9. Spreading of the anesthetic drug is achieved by turning the patient into a favourable position depending on the type of drug used (hyperbaric, isobaric or hypobaric).
- 10. The level and spread of local anesthesia needs to be carefully monitored. After reaching a complete block of the motoric and/or sensoric and/or sympathetic nerves in the desired dermatome, surgery can start.

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B. Braun Australia Pty. Ltd.
B. Braun Melsungen AG