



Edwards

AORTIC PERFUSION CANNULA
Instructions for Use
Directory

Table with 2 columns: Language and Page number. Includes English (EN), Français (FR), Deutsch (DE), Español (ES), Italiano (IT), Nederlands (NL), Dansk (DA), Svenska (SV), Ελληνικά (EL), Português (PT), Česky (CS), Magyar (HU), Polski (PL), Slovenský (SK), Norsk (NO), Български (BG), Română (RO), Eesti (ET), Lietuvių (LT), Latvīšu (LV), Türkçe (TR), Pycckий (RU), Српски (RS).

AORTIC PERFUSION CANNULA
Instructions for Use

CAUTION: Federal Law (U.S.A.) restricts this device to sale by, or on the order of, a physician.

DESCRIPTION

Edwards Lifesciences aortic perfusion cannulae terminate in a connector acceptance of 3/8" (9.5 mm) or 1/4" (6.3 mm). The catalog identifies models that are supplied with a 3/8" (9.5 mm), 1/4" (6.3 mm) connector, or vent plug.

INDICATIONS FOR USE

Aortic perfusion cannulae are intended for perfusion of the aorta during short-term (<=6 hours) cardiopulmonary bypass procedures.

CONTRAINDICATIONS FOR USE

This device is not intended for use other than as indicated and should not be used when any physical impairment would contraindicate its use.

WARNINGS AND PRECAUTIONS

Supplied sterile and non-pyrogenic in undamaged package. Do not use if device shows signs of damage (i.e., cuts, kinks, crushed areas, leakage), or if package is damaged or open as this may indicate compromised sterility and/or product damage.

FOR SINGLE USE ONLY. This device is designed, intended, and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE THIS DEVICE.

Products are known to contain phthalates, which may be found in device materials containing plasticizers such as DEHP and BPP. High exposure to such phthalates during medical treatments in children and pregnant or nursing women may result in exposure.

This device is intended for short-term use only (<=6 hours). There are no data to support the function and performance of the device beyond six hours of use.

DO NOT PRECOOL the EZ Glide cannula in an ice bath prior to use. Precooling can cause the tip to become brittle and prone to damage.

Improper orientation may lead to patient harm or excessive line pressures.

pressures. Cannula with a curved tip will have a printed orientation line on the cannula body; do not use if line is missing or illegible. Ensure that the vessel used is of adequate size to allow sufficient perfusion distal to the cannula after insertion.

DIRECTIONS FOR USE

Note: If suture ring assembly is required prior to aortic cannulation, remove suture ring from pouch and slide over cannula tip to desired location.

- 1. Using standard surgical technique, place cannula so that proper perfusion can take place.
2. For cannulae marked with an orientation line, tip direction is indicated by referencing the orientation mark on the cannula body in relationship to the direction of flow as it exits the cannula tip.
3. If supplied, the white porous vent plug is designed to fit securely into the back of the lumen of aortic perfusion cannulae without a pre-attached connector.
4. Products that state "quick venting cap" on the label (EZQ) contain a vent cap designed to allow venting while wet.
5. Secure cannula. Connect to circuit tubing and securely tie-band all line connections.
6. Upon completion of procedure, remove cannula and repair aortotomy.

Edwards, Edwards Lifesciences, the stylized E logo, and EZ Glide are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

Français (FR) French

CANULE DE PERFUSION AORTIQUE
Mode d'emploi

DESCRIPTION

L'extrémité des canules de perfusion aortique Edwards Lifesciences accepte des raccords de 9,5 mm (3/8 po) ou 6,3 mm (1/4 po). Le catalogue indique les modèles fournis avec un bouchon de ventilation ou un raccord de 9,5 mm (3/8 po) ou 6,3 mm (1/4 po).

INDICATIONS

Les canules de perfusion aortique sont indiquées pour la perfusion de l'aorte au cours d'interventions rapides d'une durée inférieure ou égale à 6 heures de circulation extracorporelle.




CHART LEGEND - LEGENDE DU GRAPHIQUE - ZEICHNERLEGENDE - SIGNIFICADO DEL CUADRO - LEGENDA DEL GRAFICO - LEGENDA TABEL - DIAGRAMFORKLÄRUNG - DIAGRAMFÖRKLÄRUNG - YEDIMNİMA (FİYAHMAT) - LEGENDA DA TABELA - LEGENDA GRAFICU - ABRA SOZVEGE - LEGENDA MYNYSSEU - YFYSYTYKYY K TABUCCKE - SIMBOLFORKLÄRUNG - KÄÄNNÖN SELITYKSET - ZEPENJA NA DIZAPKANATA - LEGENDA GRAFIC - TABELDUSE SELITYTUS - DIAGRAMOS PAUKSKINIMAS - DIAGRAMMAS SKAIDROJIMAS - GRAFİK AÇIKLAMALARI - YCDÖRMEH ÖZGÖNMEYENH DUKATPAMBA - LEGENDA GRAFICOMA

Table with 3 columns: Language, Description, and Flow Rate (L/min), H2O at Room Temperature. Includes EN, FR, DE, ES, IT, NL, DA, SV, PT, CS, HU, PL, SK, NO, FI, BG, RO, ET, LT, LV, TR, RU, RS.