
# SORIN

# CARDIOACCESSORIES

*Intracardiac Suckers - Aspiratori intracardiaci - Ventouse intracardiaque – Intrakardialsauger - Succionador intracardiaco - Aspirador Intracardiaco - Ενδοκαρδιακή αναρρόφηση - Intracardiale Zuigbuis – Intrakardialsugare - Intrakardial suger - Intrakardiaali-imu - Intrakardiální odsávačka - Ssak wewnątrzsercowy - Intrakardiálny nasávací ventil - Intrakardiális szívó - Intrakardinis siurbtukas - Súdamesisene imitoru - Aspiratoare intracardiace - Intrakardiyak Emiciler*

en - ENGLISH - INSTRUCTIONS FOR USE.....	3
it - ITALIANO - ISTRUZIONI PER L'USO .....	5
fr - FRANCAIS - MODE D'EMPLOI.....	7
de - DEUTSCH - GEBRAUCHSANWEISUNG .....	9
es - ESPAÑOL - INSTRUCCIONES DE USO .....	11
pt - PORTUGUÊS - INSTRUÇÕES DE UTILIZAÇÃO.....	13
el - ΕΛΛΗΝΙΚΑ - ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ .....	15
nl - NEDERLANDS - GEBRUIKSAANWIJZINGEN .....	17
sv - SVENSKA - BRUKSANVISNING.....	19
da - DANSK - BRUGSANVISNING .....	21

fi - SUOMI - KÄYTTÖOHJEET .....	23
cs - ČEŠTINA - NÁVOD K POUŽITÍ.....	25
pl - POLSKI - INSTRUKCJE UŻYTKOWANIA .....	27
sk - SLOVENČINA - NÁVOD NA POUŽITIE .....	29
hu - MAGYAR - KEZELÉSI ÚTMUTATÓ .....	31
lt - LIETUVIŠKAI - NAUDOJIMO INSTRUKCIJOS .....	33
et - EESTI - KASUTUSJUHEND.....	35
ro - ROMÂNĂ - INSTRUCȚIUNI DE UTILIZARE.....	37
tr - TÜRKÇE - KULLANMA TALİMATI.....	39

**I. CONTENTS**

- A. Description
- B. Intended use
- C. Contraindications
- D. Safety information
- E. Set up
- F. Use
- G. Device replacement
- H. Medical devices for use with the Intracardiac suckers
- I. Return of used products
- J. Limited Warranty

**A. DESCRIPTION**

The intracardiac suckers consist of a PVC, ABS or stainless steel body. These suckers are available in several styles and sizes:

- S331 series:  
fluted tip, flexible, with mounting clip  
Product designation: S331-70
- S351 series:  
fluted tip with ergonomic grip  
Product designation: S351-70
- S352 series:  
smooth, rounded, non-block tip with ergonomic grip  
Product designation: S352-80
- S353 series:  
pediatric, with fluted tip and ergonomic grip  
Product designation: S353-50
- S354 series:  
soft silicone tip with ergonomic grip  
Product designation: S354-25
- S361 series:  
pediatric, with fluted tip and flexible tubing  
Product designation: S361-50
- Sucker "Wessex" series  
Product designation: 072600104
- S900 series:  
Customised suckers  
Product designation: S900-xx

"B" preceding the part number on the product label indicates that the device is available with phospholipid PC coating. Devices coated with PC are used when a coated blood path is desired. The PC coating improves the blood compatibility of the device by reducing platelet adhesion on the coated surfaces.

The device is single use, non-toxic, non pyrogenic, supplied STERILE and packaged singly. Sterilised by ethylene oxide. The level of ethylene oxide residuals in the device is within the limits established by national regulations in the country of use.

**B. INTENDED USE**

The intracardiac suckers are intended to be used to remove excess fluid from the surgical field during cardiopulmonary surgery for periods of up to six hours.

All of the suckers should be used in combination with medical devices listed in section H (Medical devices for use with the Cannulae).

**C. CONTRAINDICATIONS**

The suckers are for use only as indicated. Not recommended r use in diseased or abnormal cannulation sites.

"B" preceding the part number on the product label indicates that the device is PC coated, currently Sorin Group Italia is not aware of any contraindications to the use of this coated device.

**D. SAFETY INFORMATION**

Information intended to attract the attention of the user to potentially dangerous situations and to ensure correct and safe use of the device is indicated in the text in the following way:

**WARNING**

WARNING indicates serious adverse reactions and potential safety hazards for practitioner and/or patient that can occur in the proper use or misuse of the device and also the limitations of use and the measures to be adopted in such cases.

**CAUTION**

CAUTION indicates any special care to be exercised by a practitioner for the safe and effective use of the device.

EXPLANATION OF THE SYMBOLS USED ON THE LABELS

-  For single use only (Do not reuse)
-  Batch code (number)  
(reference for product traceability)
-  Use by (Expiry date)

Manufactured by:  
Sterile - Ethylene oxide sterilised

Non Pyrogenic

Contains PHTHALATE

Warning: Do not resterilize.

Contents sterile only if package is not opened, damaged or broken

Catalogue (code) number  
Units

Attention, see instruction for use

Consult the instructions for use on the website at www.sorinmanuals.com

This way up

Fragile; handle with care

Keep away from heat

Keep dry

Diameter

Length

Shank length

Connector

Tip style / Angle

The following is general safety information with the aim of advising the operator in preparing to use the device.

Also, specific safety information is given in the instructions for use at locations in the text where that information is relevant for correct operation.

**WARNING**

- The User should carefully check the device during set-up and priming for leaks. Do not use if any leak is detected.
- The device must be used in accordance with the instructions for use provided in this manual.
- For use by professionally trained personnel only
- Sorin Group Italia is not responsible for problems arising from inexperience or improper use.
- FRAGILE, handle with care.
- Keep dry. Store at room temperature.
- Always give and maintain a correct dose and accurate monitoring of the anticoagulant before, during and after the bypass.
- For single use and for single-patient use only. During use the device is in contact with human blood, body fluids, liquids or gases for the purpose of eventual infusion, administration or introduction into the body, and due to its specific design it cannot be fully cleaned and disinfected after use. Therefore, reuse on other patients might cause cross-contamination, infection and sepsis. In addition, the reuse increases the probability of product failure (integrity, functionality and clinical effectiveness).
- If on the label is present the symbol  then the device contains Phthalates. Considering the nature of contact with the body, the limited duration of contact and the number of treatments per patient, the amount of phthalates which might be released from the device does not raise specific concerns about residual risks. Further information is available on request from Sorin Group Italia
- The device must not undergo any further processing.
- Do not resterilize.
- After use, dispose of the device in accordance with applicable regulations in force in the country of use.
- The device must only be used if sterile.- For further information and/or in case of complaint contact SORIN GROUP ITALIA or the authorised local representative.

- Do not use if sterile packaging is damaged, unsealed, or has been exposed to moisture or other conditions that would compromise the sterility of the device.
- Check the expiry date on the label attached. Do not use the device after the date shown.
- The device must be used immediately after opening the sterile packaging.
- The device must be handled aseptically.
- Carry out a visual inspection and carefully check the device before use. Transport and/or storage conditions other than those prescribed may have caused damage to the device.
- Do not use solvents such as alcohol, ether, acetone, etc.: as contact may cause damage to the device.

#### **CAUTION**

- Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

### E. SET UP

- Using care not to contaminate the sterile product, remove the sucker from its sterile pouch and place it into the sterile field.
- Check the sucker for patency.
- Attach the sucker to the suction tubing line.

Note: The intracardiac suckers with flexible tubing may be shaped to desired configuration.

#### **WARNING**

- Attachment of the suckers to tubing must be secure to prevent loss of suction.
- Affix the sucker to the surgical drape using drape attachment rings.
- Adjust roller pump speed for appropriate suction.
- Place the sucker into the desired cavity for fluidremoval.
- Suckers that are not hand held at all times must be properly affixed to the surgical drape using drape attachment rings to prevent inadvertent removal during operating procedures.

#### **CAUTION**

- Lumen patency should be verified on each sucker prior to use
- Excessive force placed on the sucker tips may causebreakage.
- Check the security of all luer connections. All accessory lines connected to the device must be connected tightly in order to prevent accidental introduction of air into the device or blood loss.

### F. USE

Various styles of suckers are available to maximize ease and efficiency of fluid removal based on the anatomical and physiological condition of the surgical field.

Note: Due to the complex nature of the surgical procedures in which these cannulae are employed, the selection of a particular cannula and its actual method of use are left to the discretion of the individual surgeon.

Orientation of the cannulae and associated tubing should allow for smooth radius bends to preclude kinking of the perfusion circuit.

#### **WARNING**

- The blood to be treated should contain anticoagulant. Always apply and maintain a correct dose and accurate monitoring of the anticoagulant before, during and after the bypass.

#### **CAUTION**

- The following complication may occur upon usage of the indicated suckers:  
– abnormal hemolysis

### G. DEVICE REPLACEMENT

A spare device must always be available during perfusion. After 6 hours of use with blood or if situations occur, which lead the person responsible for perfusion to determine that the safety of the patient may be compromised, consider device replacement.

#### **CAUTION**

- Use sterile methods during all replacement procedure.

### H. MEDICAL DEVICES FOR USE WITH THE INTRACARDIAC SUCKERS

All tubing used to make the circuit connections must be of a diameter which is compatible with the dimensions of the connectors on the device (1/4").

Currently SORIN GROUP ITALIA is not aware of any contraindications to the use of the device with occlusive or non-occlusive peristaltic pumps or with centrifugal pumps. The use of other types of pump must be agreed with SORIN GROUP ITALIA.

#### **CAUTION**

The user should observe the warnings and cautions and follow instruction for Use accompanying the separate device.

### I. RETURN OF USED PRODUCTS

Should the user be dissatisfied with anything related to the quality of the product, the product distributor or the authorized local SORIN GROUP ITALIA representative should be notified.

All parameters considered critical by the user must be reported with particular care and urgency. The following is the minimum information that should be provided:

- Detailed description of the event and, if pertinent, the conditions of the patient;
- Identification of the product involved;

- Lot number of the product involved;
- Availability of the product involved;
- All the indications the user considers useful in order to understand the origin of the elements of dissatisfaction.

SORIN GROUP ITALIA reserves the right to authorize, if necessary, recall of the product involved in the notification for assessment. If the product to be returned is contaminated, it must be treated, packed and handled in conformity with the provisions of the legislation in force in the country where the product was used.

#### **CAUTION**

It is the responsibility of the health care institution to adequately prepare and identify the product for return shipment. Do not return products that have been exposed to blood borne infectious diseases.

#### **ONLY for US customers**

If for any reason the product must be returned to the manufacturer, a returned good authorisation (RGA) number is required from LivaNova USA, Inc. prior to shipping .

If the product has been in contact with blood or blood fluids, it must be thoroughly cleaned and disinfected before packing. It should be shipped in either the original carton or an equivalent carton to prevent damage during shipment, and it should be properly labelled with an RGA number and an indication of the biohazardous nature of the content in the shipment.

#### **CAUTION**

It is the responsibility of the health care institution to adequately prepare and identify the product for return shipment. Do not return products that have been exposed to blood borne infectious diseases.

The shipping address for returned goods in the US is:

**LivaNova USA, Inc.**  
14401 West 65<sup>th</sup> Way  
Arvada, CO USA 80004

### J. LIMITED WARRANTY

*This Limited Warranty is in addition to any statutory rights of the Purchaser pursuant to applicable law.*

SORIN GROUP ITALIA warrants that all reasonable care has been taken in the manufacture of this medical device, as required by the nature of the device and the use for which the device is intended.

SORIN GROUP ITALIA warrants that the medical device is capable of functioning as indicated in the current instructions for use when used in accordance with them by a qualified user and before any expiry date indicated on the packaging.

However, SORIN GROUP ITALIA cannot guarantee that the user will use the device correctly, nor that the incorrect diagnosis or therapy and/or that the particular physical and biological haematological characteristics of an individual patient, do not affect the performance and effectiveness of the device with damaging consequences for the patient, even though the specified instructions for use have been respected.

SORIN GROUP ITALIA, whilst emphasizing the need to adhere strictly to the instructions for use and to adopt all the precautions necessary for the correct use of the device, cannot assume any responsibility for any loss, damage, expense, incidents or consequences arising directly or indirectly from the improper use of this device.

SORIN GROUP ITALIA undertakes to replace the medical device in the event that it is defective at the time of placing on the market or whilst being shipped by SORIN GROUP ITALIA up to the time of delivery to the final user unless such defect has been caused by mishandling by the purchaser.

The above replaces all other warranties explicit or implicit, written or verbal, including warranties of merchantability and fitness for purpose. No person, including any representative, agent, dealer, distributor or intermediary of SORIN GROUP ITALIA or any other industrial or commercial organization is authorized to make any representation or warranty concerning this medical device except as expressly stated herein. SORIN GROUP ITALIA disclaims any warranty of merchantability and any warranty of fitness for purpose with regard to this product other than what is expressly stated herein. The purchaser undertakes to comply with the terms of this Limited Warranty and in particular agrees, in the event of a dispute or litigation with SORIN GROUP ITALIA, not to make claims based on alleged or proven changes or alterations made to this Limited Warranty by any representative, agent, dealer, distributor or other intermediary.

The existing relations between the parties to the contract (also in the case that it is not drawn up in writing) to whom this Warranty is given as well as every dispute related to it or in any way connected to it as well as anything related to it or any dispute concerning this Warranty, its interpretation and execution, nothing excluded and/or reserved, are regulated exclusively by the Italian law and jurisdiction. The court chosen is the Court of Modena (Italy).

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**Sorin Group Italia S.r.l.**