Valleylab™

Single-Use Pencils

REF E2515**REF** E2515H**REF** E2516**REF** E2516H**en** Instructions for Use**fr** Mode d'emploi**pt** Instruções de uso**da** Brugervejledning**de** Gebrauchsanleitung**es** Instrucciones de uso**fi** Käyttöohjeet**it** Istruzioni per l'uso**nl** Gebruiksaanwijzing**no** Bruksanvisning**sv** Bruksanvisning**ru** Инструкция по применению**zh** 使用说明

Part No. PT00135747

COVIDIEN, COVIDIEN with logo, and Covidien logo and Positive Results for Life are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company, ™ brands are trademarks of their respective owner.

© 2011 Covidien.
 Covidien llc,
15 Hampshire Street, Mansfield, MA 02048
USA.
EC/REP Covidien Ireland Limited,
IDA Business & Technology Park, Tullamore.

Made in Mexico. Printed in Mexico.

www.covidien.com

REV 10/2020

STERILE R

Single use

Rx
ONLY

Do not re-sterilize

Consult
instructions
for use

2797

Valleylab™

REF E2515 **Rocker Switch Pencil
10' (3 m) Cord**

REF E2515H **Rocker Switch Pencil,
Holster
10' (3 m) Cord**

REF E2516 **Button Switch Pencil
10' (3 m) Cord**

REF E2516H **Button Switch Pencil,
Holster
10' (3 m) Cord**

For use with a max peak voltage of 5600 V.

Not made with natural rubber latex

Do not use if package is opened
or damaged

Indications for Use

The Valleylab Handswitching Pencil is a handswitching electrode intended for use in general electrosurgical procedures.

Warning

This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in bio-incompatibility, infection, or product failure risks to the patient.

Fire Hazard Do not place active accessories near or in contact with flammable materials, such as gauze or surgical drapes, flammable gases, or high levels of oxygen. Electrosurgical accessories that are activated or hot from use can cause a fire. When not using active accessories, place them in a holster or in a clean dry, nonconductive, and highly visible area that is not in contact with the patient. Inadvertent contact with the patient may result in burns.

Danger: Explosion Hazard Do not use electrosurgery in the presence of flammable anesthetics.

Fire Hazard Both oxygen (O₂) and nitrous oxide (N₂O) support combustion. Avoid atmospheres enriched with O₂ and N₂O. Enriched atmospheres may result in fires and burns to the patient or surgical personnel.

Warning

Fire/Explosion Hazard The following substances contribute to increased fire and explosion hazards in the operating room:

- Flammable substances, such as alcohol-based skin prepping agents and tinctures
- Naturally occurring flammable gases that may accumulate in body cavities, such as the bowel
- O₂-enriched atmospheres
- Oxidizing agents, such as N₂O atmospheres

The sparking and heating associated with electrosurgery can provide an ignition source. Observe fire precautions at all times. When using electrosurgery in the same room with any of these substances or gases, prevent their accumulation or pooling under surgical drapes, or within the area where electrosurgery is performed.

Fire/Explosion Hazard Verify all oxygen circuit connections are leak free before and during the use of electrosurgery. Verify that endotracheal tubes are leak free, and the cuff is properly sealed to prevent oxygen leaks. Enriched oxygen atmospheres may result in fires and burns to patients or the surgical team.

Position surgical electrode cables to avoid contact with the patient or other leads.

When not in use, place active accessories in a holster or in a clean, dry, nonconductive, and highly visible area not in contact with the patient. Inadvertent contact with the patient may result in burns.

Before installing or removing the electrode, ensure the handset is not connected to the electrosurgical generator, or the generator is off or in standby mode, if available.

Confirm proper electrosurgical generator power settings before proceeding with surgery. Use the lowest power settings to achieve the desired surgical effect.

Warning

Some surgeons may elect to “buzz the hemostat” during surgical procedures. It is not recommended, and the hazards of such a practice probably cannot be eliminated. Burns to the surgeon’s hands are possible. To minimize the risk:

- Do not lean on the patient, the table, or the retractors while buzzing the hemostat.
- Activate cut rather than coag. Cut has a lower voltage than coag.
- Use the lowest power setting possible for the minimum time necessary to achieve hemostasis.
- Activate the generator after the accessory makes contact with the hemostat. Do not arc to the hemostat.
- Firmly grasp as much of the hemostat as possible before activating the generator. This disperses the current over a larger area and minimizes the current concentration at the fingertips.
- “Buzz the hemostat” below hand level (as close as possible to the patient) to reduce the opportunity for current to follow alternate paths through the surgeon’s hands.
- When using a stainless steel blade electrode, place the **flat** surface against the hemostat or other metal instrument.
- When using a coated or nonstick blade electrode, place the **edge** of the electrode against the hemostat or other metal instrument.

Precaution


Ensure the electrode is securely seated in the pencil. An improperly installed electrode may result in injury to the patient or surgical team by arcing at the electrode and pencil connection.

Inspect accessories and cords for breaks, cracks, nicks, or other damage before every use. If damaged, do not use. Failure to observe this caution may result in injury or electrical shock to the patient or surgical team.


Before Surgery

1. Remove the pencil (and holster) from the package using aseptic technique.
2. Change the electrode orientation if necessary. Grasp the electrode by the insulating sleeve, pull the electrode out of the pencil, rotate, and reinsert it into the pencil. The blade electrode (supplied) is designed to lock in place to prevent rotation.

The pencil also accepts standard 3/32” (2.4 mm) diameter non-locking electrodes.


3. Remove the tip protector from the electrode. Insert the pencil connector into the *handswitching* receptacle of a Covidien generator.


An adapter may be necessary to connect the pencil to generators not manufactured by Covidien. Refer to the *Surgical Products Catalog* for the appropriate adapter.

4. Attach the safety holster to the sterile drape with tape or by pulling the drape material through the slot as shown (E2400 Accessory Holster supplied with E2515H/E2516H).


During Surgery

Place the pencil in the holster when not in use. If available, place the generator in standby mode until needed for surgery.

After Surgery

Discard the pencil, electrode, and holster after use. They are not designed to withstand resterilization. **Do not resterilize.**

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