


**REF**

**PS3551C4NS, PS3551C8NS**

**CE**

<b>en</b>	Monopolar instrument cables
<b>es</b>	Cables de instrumentos monopolares
<b>pt</b>	Cabos de instrumentos monopolares
<b>de</b>	Monopolare Instrumentenkabel
<b>fr</b>	Câbles d'instruments monopolaires
<b>it</b>	Cavi strumenti monopolari
<b>cs</b>	Monopolární přístrojové kabely



Purple Surgical International Limited / Purple Surgical UK limited  
2 Chestnut House, Farm Close, Shenley, Hertfordshire,  
WD7 9AD, England  
Tel: +44 (0) 1923 839333  
E-Mail: [mail@puplesurgical.com](mailto:mail@puplesurgical.com)  
Web: [www.puplesurgical.com](http://www.puplesurgical.com)

**EC REP**

Advena Limited  
Tower Business Centre, 2<sup>nd</sup> floor, Tower Street,  
Swatar, BKR 4013, Malta



**en - The following should be thoroughly read before use.**

**IMPORTANT NOTES**

This document is designed to provide instructions on how this product/s is used. It is not a reference to surgical techniques.

**INDICATIONS**

The Purple Surgical 'Monopolar instrument cables' should be used with suitable monopolar devices and an appropriate electrosurgical generator.

**PRODUCT CHARACTERISTICS**

The Instrument/s are designed for Monopolar use only

- Device conforms to Electrical Safety Standard IEC 60601-2-2.

- The maximum rated voltage this device has been tested in conformance to is 3,500V/peak.

**HOW SUPPLIED**

Supplied in Non sterile format to Kit packer, for sterilisation as part of a kit prior to market release to intended user.

**SHELF LIFE**

4 years from date of sterilisation - use by date as denoted on labelling in format YYYY-MM

**INDICATIONS FOR STERILISATION PROCESSING**

The device performances have been qualified as suitable with consideration to the following method, and limitations, of sterilisation:

— Ethylene Oxide Sterilisation - Maximum Temperature of 59°C and Minimum Pressure of 50mBar

Note\* It is the responsibility of the kit packer to validate the sterilisation process employed, with respect to the conditions of their presentation of the packaging system.

**COMPATIBILITY AND CONNECTIVITY**

The Purple Surgical 'Monopolar instrument cables' are used to connect a monopolar device by means of the 4mm female socket on the cable to the 4mm Monopolar Instrument post on the device.

The cables are available with either a 4mm (PS3551CA) or 8mm (PS3551CB) monopolar pin for connection to the electrosurgical generator.

**CONTRAINDICATIONS**

- This device is not to be used for connection to bipolar instrumentation.

- DO NOT USE on patients who have electronic implants such as cardiac pacemakers without first consulting a qualified professional (e.g., cardiologist). A possible hazard exists because interference with the action of the electronic implant may occur, or the implant may be damaged.

**PRECAUTIONS / WARNINGS**

(1) DO NOT use if package is received opened or damaged.

(2) Procedures should be performed only by clinicians with adequate training and knowledge. Medical literature should be consulted for techniques, hazards, contraindications and complications prior to the procedure.

(3) Connect adaptors and accessories to the electrosurgical unit only when the unit is off. Failure to do so may result in an injury or electrical shock to the patient or operating room personnel.

(4) Refer to appropriate electrosurgical system user manual for indications and instructions to ensure that all safety precautions are followed.

(5) Laparoscopic instruments may vary in diameter from manufacturer to manufacturer. When laparoscopic instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility prior to initiation of procedure.

(6) The user should review the instructions for use of all devices, being used in conjunction with these product/s.

(7) During monopolar surgery, the PATIENT should not come into contact with metal parts which are earthed or which have an appreciable capacitance to earth (for example operating table supports, etc.).

(8) During monopolar surgery, Skin-to-skin contact (for example between the arms and body of the PATIENT) should be avoided, for example by insertion of dry gauze.

(9) During monopolar surgery, the PATIENT leads should be positioned in such a way that contact with the PATIENT or other leads is avoided.

Temporarily unused ACTIVE ELECTRODES should be stored in a location that is isolated from the PATIENT.

(10) During monopolar surgery, the output power selected should be as low as possible for the intended purpose, and it shall not exceed the maximum rated voltage of the cable and any connected device

(11) During monopolar surgery, the use of flammable anaesthetics or oxidizing gases such as nitrous oxide (N<sub>2</sub>O) and oxygen should be avoided if a surgical procedure is carried out in the region of the thorax or the head, unless these agents are sucked away.

(12) Non-flammable agents should be used for cleaning and disinfection wherever possible. Flammable agents used for cleaning or disinfecting, or as solvents of adhesives should be allowed to evaporate before the application of HF surgery. There is a RISK of pooling of flammable solutions under the PATIENT or in body depressions such as the umbilicus, and in body cavities such as the vagina. Any fluid pooled in these areas should be mopped up before HF SURGICAL EQUIPMENT is used.

(13) Attention should be called to the danger of ignition of endogenous gases. Some materials, for example cotton and gauze, when saturated with oxygen may be ignited by sparks produced in NORMAL USE of the HF SURGICAL EQUIPMENT.

(14) For PATIENTS with electrically conductive implants, a possible HAZARD exists due to concentration or re-direction of HF currents. In case of doubt, qualified advice should be obtained.

(15) During Monopolar Surgery, it is recommended that the use of smoke-plume extraction is employed.

**ADVERSE REACTIONS**

No known issues.

**DIRECTIONS FOR USE**

(1) Visually inspect the sterile barrier system pouch prior to its opening so as to ensure it is not already open or compromised. If a sterile barrier failure is found, the contained product shall be considered non-sterile and it shall not be used.

(2) Using aseptic technique, remove the product/s from its packaging.

(3) The 4mm female socket on the cable should be connected to a standard 4mm monopolar diathermy post on any monopolar instrument.

(4) The male connector on the cable should be connected to the relevant monopolar power outlet on an appropriate electrosurgical generator, whilst the generator is switched off.

**SINGLE USE PRECAUTIONS**

This product/s are designed and sold for single use only. Re-processing and, or, re-sterilisation is not permitted.

The effects of any unauthorised re-processing or re-sterilisation can result in the following complications:

1. Cross contamination due to ineffective re-processing/re-sterilisation.

2. Mechanical fatigue, and associated failure, due to the effects of the re-processing / re-sterilisation method.

**DISPOSAL**

Discard after single use. This product/s is to be disposed of as controlled medical waste according to national guidelines.

**SERIOUS INCIDENT REPORTING**

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.



en: LABELLING SYMBOLOGY - Standard ISO 15223-1, unless otherwise stated  
 es: SIMBOLOGÍA DEL ETIQUETADO - Norma ISO 15223-1 salvo que se indique lo contrario  
 pt: SÍMBOLOS UTILIZADOS NO RÓTULO - Norma ISO 15223-1, excepto indicação em contrário  
 de: KENNZEICHNUNGS SYMBOLOGIE - Standardmäßig laut ISO 15223-1, falls nicht anders vermerkt  
 fr: SYMBOLES D'ÉTIQUETAGE - Norme ISO 15223-1, sauf indication contraire  
 It: SIMBOLI SULL'ETICHETTA - Standard ISO 15223-1, se non diversamente specificato  
 cs: SYMBOLY OZNAČENÍ - norma ISO 15223-1 pokud není uvedeno jinak

	<p><b>EC REP</b></p>	<p><b>MD</b></p>	<p><b>CE</b></p>
<p>en: Manufacturer                  es: Fabricante                  pt: Fabricante                  de: Hersteller                  fr: Fabricant                  It: Produttore                  cs: Výrobce</p>	<p>en: Authorized representative                  es: Representante autorizado                  pt: Representante autorizado                  de: Bevollmächtigten                  fr: Mandataire                  It: Rappresentante autorizzato                  cs: Zplnomocněný zástupce</p>	<p>en: Medical device                  es: Producto sanitario                  pt: Dispositivo médico                  de: Medizinprodukt                  fr: Dispositif médical                  It: Dispositivo medico                  cs: Zdravotnický prostředek</p>	<p>en: CE Mark                  es: Marca CE                  pt: Marca CE                  de: CE-Kennzeichnung                  fr: CE Marquage                  It: Marchio CE                  cs: Značka CE</p>
<p><b>REF</b></p>	<p><b>LOT</b></p>		
<p>en: Catalogue number                  es: Número de catálogo                  pt: Número de catálogo                  de: Katalognummer                  fr: Numéro de catalogue                  It: Numero di catalogo                  cs: Katalogové číslo</p>	<p>en: Batch code                  es: Código de lote                  pt: Código de lote                  de: Chargen-Code                  fr: Code de lot                  It: Codice lotto                  cs: Šíko farže</p>	<p>en: Consult instructions for use                  es: Consultar las instrucciones de uso                  pt: Consultar instruções de utilização                  de: Gebrauchsanleitung beachten                  fr: Consulter le mode d'emploi                  It: Consultare le Istruzioni per l'uso                  cs: Viz pokyny pro použití</p>	<p>en: Do not reuse                  es: No reutilizar                  pt: Não reutilizar                  de: Nicht wiederverwenden                  fr: Ne pas réutiliser                  It: Non riutilizzare                  cs: Nepoužívejte opakovaně</p>
<p>en: Keep dry                  es: Mantener en un lugar seco                  pt: Manter seco                  de: Trocken aufbewahren.                  fr: À conserver dans un endroit sec                  It: Conservare all'asciutto                  cs: Uchovávejte v suchu</p>	<p>en: Keep away from sunlight                  es: Mantener alejado de la luz del sol                  pt: Manter ao abrigo da luz solar                  de: Fern von Sonnenlicht aufbewahren.                  fr: À conserver à l'abri de la lumière                  It: Tenere lontano dalla luce del sole                  cs: Chraňte před slunečním světlem</p>	<p>en: Do not use if package is damaged                  es: No usar si el envase está dañado                  pt: Não usar caso o embalagem esteja danificada                  de: Nicht verwenden, wenn das Paket beschädigt ist                  fr: Ne pas utiliser si l'emballage est endommagé                  It: Non usare se la confezione è danneggiata                  cs: Nepoužívejte, pokud je obal poškozen</p>	<p>en: Caution                  es: Precaución                  pt: Atenção                  de: Vorsicht                  fr: Avertissement                  It: Attenzione                  cs: Upozornění</p>
<p>en: Type BF applied part                  es: Escribir la pieza aplicada de tipo BF                  pt: Peca aplicada de tipo BF                  de: Anwendungsteil vom Typ BF                  fr: Partie appliquée de type BF                  It: Parte applicata di tipo BF                  es: Aplikované součást typu BF</p>	<p>en: Use-by date                  es: Fecha de uso recomendado                  pt: Data de validade                  de: Verfallsdatum                  fr: Date limite d'utilisation                  It: Usare entro la data                  cs: Spotřebujte do</p>	<p>en: Non Sterile                  es: No estéril                  pt: Não estéril                  de: Nicht steril                  fr: Non stérile                  It: Non Sterile                  cs: Není sterilní</p>	<p>en: Importer                  es: Importador                  pt: Importador                  de: Importeur                  fr: Importateur                  It: Importatore                  es: Dovezice</p>
<p><b>CH REP</b></p>	<p><b>UK CA</b></p>		
<p>en: Authorized representative in Switzerland                  es: Mandatario en Suiza                  pt: Representante autorizado na Suíça                  de: Bevollmächtigten in der Schweiz                  fr: Mandataire en Suisse                  It: Mandatario in Svizzera                  cs: Zplnomocněný zástupce ve Švýcarsku</p>	<p>en: UKCA Mark                  es: Marcado UKCA                  pt: Marcado UKCA                  de: UKCA-Kennzeichnung                  fr: Marquage UKCA                  It: Marchio UKCA                  cs: Značka UKCA se schváleným číslem organizace</p>		