



Introducing Our Full Range of Personal Protective Equipment

PPE



Unisurge

Unisurge have now expanded our PPE portfolio to ensure that we can help all trusts across every department of a hospital. These include Masks, Eyewear, Visors, Headwear, Gloves along with Gowns, Scrub suits, Warm up jackets and overshoes. This offering of products is now being taken by NHSSC and NHS Trusts. These items are built to ensure that hospitals have the high quality required for all patients and staff.

Unisurge work closely with trusts across the UK, we offer a range of high quality products at competitive prices. Medical Professionals Nationwide request our bespoke Procedure Packs, Swabs, Needles, Syringes, Polyware and Single use Instruments along with an extensive portfolio of Drapes and Gowns. Unisurge's offices, warehouse and sterilisation facilities are all based in the UK, this gives us the flexibility to deliver stock on our own vehicles when our customers need them.





PPE



The right products in the right hands at the right price

At Unisurge our philosophy is quite simple:

To provide our customers with the service and expertise they would expect from a supplier of products of such critical importance to the day-to-day running of their healthcare department.

All products listed
are latex free





Masks

A comprehensive range of medical surgical masks that fully comply with the relevant standards.


Requirements for surgical masks according to European standards

	Type IIR	N95	FFP2	FFP3
Bacterial filtration efficiency (BFE, %)	>98	>95	>94	>99
Breathing Resistance (PA/cm ²)	<49			
Breathing Resistance 30L / Min (mbar)			0.7	1
Breathing Resistance 95L / Min (mbar)			2.4	3
Inhalation (mmH2O)		≤ 35		
Exhalation (mmH2O)		≤ 25		
Fluid Resistance (mm Hg)	>120			
Inward Leakage (%)			<8	>2
Fit Test Required	No	Yes	Yes	Yes
Standard	EN 14683	NIOSH	EN 149	EN 149




Mask Type IIR with Ear Loops Non Sterile

Code	Description	Bacterial filtration efficiency	Breathing Resistance (PA/cm ²)	Fluid resistance (mm Hg)	EN14683 type	Case Qty
F931607	Face mask, blue	>98	<49	>120	IIR	20x50




Mask Type IIR with 4 Ties Non Sterile

Code	Description	Bacterial filtration efficiency	Breathing Resistance (PA/cm ²)	Fluid resistance (mm Hg)	EN14683 type	Case Qty
F961606	Face mask, blue	>98	<49	>120	IIR	20x50



Mask Type IIR with Anti fog Visor and Ear Loops Non Sterile

Code	Description	Bacterial filtration efficiency	Breathing Resistance (PA/cm ²)	Fluid resistance (mm Hg)	EN14683 type	Case Qty
F931605	Face mask with visor	>98	<49	>120	IIR	20x20



Mask Type IIR with Anti fog Visor and Ties Non Sterile

Code	Description	Bacterial filtration efficiency	Breathing Resistance (PA/cm ²)	Fluid resistance (mm Hg)	EN14683 type	Case Qty
F931604	Face mask with visor	>98	<49	>120	IIR	20x20




Mask FFP2 Valved Ear Loops Non Sterile

Code	Description	Bacterial filtration efficiency	Breathing Resistance 30L / Min (mbar)	Breathing Resistance 95L / Min (mbar)	EN149 type	Case Qty
F931620	Face mask	>94	0.7	2.4	FFP2	20x20




Mask FFP3 Valved Ear Loops Non Sterile

Code	Description	Bacterial filtration efficiency	Breathing Resistance 30L / Min (mbar)	Breathing Resistance 95L / Min (mbar)	EN149 type	Case Qty
F931622	Face mask	>99	1	3	FFP3	20x5



Mask FFP2 Unvalved Ear Loops Non Sterile

Code	Description	Bacterial filtration efficiency	Breathing Resistance 30L / Min (mbar)	Breathing Resistance 95L / Min (mbar)	EN149 type	Case Qty
F931621	Face mask	>94	0.7	2.4	FFP2	20x20




Mask N95 Shaped Ear Loops Non Sterile

Code	Description	Bacterial filtration efficiency	Inhalation	Exhalation	NIOSH type	Case Qty
F931623	Face mask	>95	≤ 35	≤ 25	N95	20x20

Eyewear

Our range of eye wear and visors all have anti fog coating to ensure clear vision and comply to European Standard EN166: 2001



Goggles Protective with Headband Non Sterile

Code	Description	Size	Colour	Case Qty
F931615	Protective goggles with optional valve	Universal	Clear	20x10



Protective Plastic Glasses with Frame Non Sterile

Code	Description	Size	Colour	Case Qty
F931616	Protective glasses	Universal	Various	5x10

Visors

Featuring full face coverage, anti-fog plastic, adjustable headbands and complies to European standard (EU) 2016/425



P.E.T Visor with Polyurethane Non Latex Elastic Head Cushioning Non Sterile

Code	Description	Size	Colour	Case Qty
F931610	Protective faceshield with headband	33 x 22cm	Purple/Clear	200

Headwear

Our portfolio is designed with comfort and protection in mind. They fit well, provide full hair coverage and the temperature regulation properties offer user comfort and meet the PPE Standard (EU) 2016/425.

HEADWEAR



Surgeons Hood Non Sterile

Code	Description	Colour	Case Qty
F931501	Cap theatre hood	Blue	10x100



Nurses Cap Non Sterile

Code	Description	Colour	Case Qty
F931502	Cap theatre elasticated back	Blue	20x100
F931503	Cap theatre elasticated back	Green	20x100
F931504	Cap theatre elasticated back	Pink	20x100
F931505	Cap theatre elasticated back	White	20x100
F931506	Cap theatre elasticated back	Yellow	20x100
F931507	Cap theatre elasticated back	Red	20x100



Bouffant Non Sterile

Code	Description	Colour	Case Qty
F931508	Cap theatre bouffant	Blue	10x100
F931510	Cap theatre bouffant	Green	10x100



Patient Cap Non Sterile

Code	Description	Colour	Case Qty
F931509	Cap theatre bouffant	White	10x100



Surgeons Cap Tie Back High Front Non Sterile

Code	Description	Colour	Case Qty
F931516	Cap theatre tie back high front	Red	6x100
F931517	Cap theatre tie back high front	Yellow	6x100
F931518	Cap theatre tie back high front	Blue	6x100
F931519	Cap theatre tie back high front	Green	6x100
F931520	Cap theatre tie back high front	Pink	6x100
F931521	Cap theatre tie back high front	White	6x100

Gloves

Medical gloves provide maximum protection against cross contamination during invasive procedures, and fully comply with the requirements of the European Standard EN455 parts 1, 2 and 3. Available both sterile and non-sterile.



Gloves Sterile


Code	Description	Size	Case Qty
F840001	Gloves examination nitrile powder free cuffed and walled sterile pair	Small	8x50 pairs
F840002	Gloves examination nitrile powder free cuffed and walled sterile pair	Medium	8x50 pairs
F840003	Gloves examination nitrile powder free cuffed and walled sterile pair	Large	8x50 pairs
F840004	Gloves examination nitrile powder free cuffed and walled sterile pair	Extra large	8x50 pairs

Gloves Non Sterile

Code	Description	Size	Case Qty
F840005	Gloves examination nitrile powder free non-sterile	Small	10x100
F840006	Gloves examination nitrile powder free non-sterile	Medium	10x100
F840007	Gloves examination nitrile powder free non-sterile	Large	10x100
F840008	Gloves examination nitrile powder free non-sterile	Extra large	10x100


Gowns

All Uni-Form surgical gowns fully comply with the requirements of the European Standard EN13795.




Gown Isolation PP/PE Non Woven 38GSM Non Sterile

Code	Description	Size	Colour	Case Qty
F931029	Gown isolation / cover yellow with knitted cuffs non sterile	Universal	Yellow	100




Gown Isolation/Cover PP Non Woven 30GSM Non Sterile

Code	Description	Size	Colour	Case Qty
F931028	Gown isolation / cover blue non sterile	Universal	Blue	100




Protective Apron/Gown with Thumb Loops Non Sterile

Code	Description	Size	Colour	Case Qty
F935034	Protective apron/gown with long sleeve & thumb loops non-sterile	Universal	Blue	5x20




Gown Surgical Standard Sterile

Code	Description	Size	Case Qty
F934001	Gown surgical standard SMS	Petite / Small	32
F934002	Gown surgical standard SMS	Medium	32
F934003	Gown surgical standard SMS	Large	32
F934004	Gown surgical standard SMS	Large, X Long	30
F934005	Gown surgical standard SMS	X Large	30
F934006	Gown surgical standard SMS	X Large, X Long	30
F934007	Gown surgical standard SMS	XX Large	22
F934008	Gown surgical standard SMS	XX Large, X Long	22




Gown Surgical Reinforced Sterile

Code	Description	Size	Case Qty
F934009	Gown surgical reinforced SMS	Petite / Small	30
F934010	Gown surgical reinforced SMS	Medium	30
F934011	Gown surgical reinforced SMS	Large	30
F934012	Gown surgical reinforced SMS	Large, X Long	28
F934013	Gown surgical reinforced SMS	X Large	28
F934014	Gown surgical reinforced SMS	X Large, X Long	28
F934015	Gown surgical reinforced SMS	XX Large	20
F934016	Gown surgical reinforced SMS	XX Large, X Long	20



Gown Surgical Impervious Sterile

Code	Description	Size	Case Qty
F934017	Gown surgical impervious SMS	Petite / Small	30
F934018	Gown surgical impervious SMS	Medium	30
F934019	Gown surgical impervious SMS	Large	30
F934020	Gown surgical impervious SMS	Large, X Long	28
F934021	Gown surgical impervious SMS	X Large	28
F934022	Gown surgical impervious SMS	X Large, X Long	28
F934023	Gown surgical impervious SMS	XX Large	20
F934024	Gown surgical impervious SMS	XX Large, X Long	20




Gown Chemotherapy Knitted Cuffs Full Back Tie Neck Polycoated Spun Bond, Virtually Lint Free Non Sterile

Code	Description	Size	Case Qty
F935148	Gown chemotherapy knitted cuffs full back tie neck closure non sterile	Large	60
F935149	Gown chemotherapy knitted cuffs full back tie neck closure non sterile	X-Large	60


Test Chemotherapy Drug	Minimum Breakthrough Detection Time (Minutes)
Carmustine	No breakthrough up to 240 minutes
Cisplatin	No breakthrough up to 240 minutes
Cyclophosphamide (Cytoxan)	No breakthrough up to 240 minutes
Dacarbazine (DTIC)	No breakthrough up to 240 minutes
Doxorubicin Hydrochloride	No breakthrough up to 240 minutes
Etoposide (Toposar)	No breakthrough up to 240 minutes
Fluorouracil	No breakthrough up to 240 minutes
Mitomycin C	No breakthrough up to 240 minutes
Paclitaxel (Taxol)	No breakthrough up to 240 minutes
Thiotepa	No breakthrough up to 240 minutes
Vincristine Sulfate	No breakthrough up to 240 minutes





Scrub Top with Pockets Short Sleeve Non Sterile

Code	Description	Size	Case Qty
F931039	Scrub top with pockets short sleeve, super soft	Medium	60
F931040	Scrub top with pockets short sleeve, super soft	Large	60
F931041	Scrub top with pockets short sleeve, super soft	XL	60



Scrub Pant Elastic Waist with Pockets Non Sterile

Code	Description	Size	Case Qty
F931042	Scrub pant elastic waist with pockets, super soft	Medium	60
F931043	Scrub pant elastic waist with pockets, super soft	Large	60
F931044	Scrub pant elastic waist with pockets, super soft	XL	60




Warm Up Jacket Knitted Cuffs Collar Non Sterile

Code	Description	Size	Case Qty
F931035	Warm up jacket knitted cuffs collar 3 pockets, neck loop	Large	24
F931036	Warm up jacket knitted cuffs collar 3 pockets, neck loop	X-Large	24
F931037	Warm up jacket knitted cuffs collar 3 pockets, neck loop, half sleeves	Large	48
F931038	Warm up jacket knitted cuffs collar 3 pockets, neck loop, half sleeves	X-Large	48



Overshoes

CPE overshoes with latex free elasticated edges, manufactured in accordance with PPE standard (EU) 2016/425



Overshoe Non Sterile

Code	Description	Size	Colour	Case Qty
F8999-300	Overshoe CPE	Universal	Blue	20x100







For more details on our comprehensive product range,
unbeatable prices call us now on

084 55 21 31 11

UNISURGE International Limited

Farboud Innovation Park,
Formula Drive, Newmarket,
CB8 0BF United Kingdom

Telephone: +44 (0) 84 55 21 31 11

Fax: +44 (0) 84 55 21 32 10

Email: info@unisurge.com

Web: www.unisurge.com