



Fine Surgical Instruments







the right products
in the right hands
at the right price

To provide our customers with the service and expertise they would expect from a manufacturer of such critically important products used in the day-to-day running of their clinical department.





what makes Unisurge different?

Established in 1988, Unisurge is a trusted manufacturer and supplier of practitioner-specified procedure packs, dressing packs, medical disposable and theatre products to health care professionals. Our range of bespoke single use products is both cost effective and patient specific. Surgical teams can be confident that all the necessary components for a procedure will be supplied, to their specification, in a sterile pack, ready for immediate use.

Our promise to you:

Experienced sales team ready to take your order and advise you on delivery times and other Unisurge products – we aim to surpass all the standards of your current service.

Innovative products constantly being introduced to increase your purchasing efficiency.

A comprehensive range of products always kept in stock ready for fast-track delivery to your door.

Regular deliveries throughout the UK.

All at the most cost effective prices.



our commitment to you...

From our Newmarket headquarters, the development team keeps abreast of the latest advances in surgical practice and takes full account of the consistent rise in patient expectations.

Unisurge offers a comprehensive range of theatre quality surgical instruments. They are available as re-usable as well as single-use products. These include individual supplementary instruments, instrument sets and instruments that can be supplied as part of a sterile procedure pack which can include all the components for any surgical procedure regardless of the speciality. The contents of each of our packs can be practitioner specified to meet the requirement of all modern surgical procedures including drapes and gowns.

Linked to a significant investment in state of the art information technology, we provide a proven and trusted distribution network. Using our own fleet of vehicles, Unisurge is able to ensure that sterile packs arrive at their destination on time and in optimum condition.


**Our aim is to tailor our service to suit your specific needs
and quite simply surpass the expectations you place upon us.**

Signature Instruments...

A comprehensive range of instruments that have been manufactured and hand finished by craftsmen. There might be slight variations in the dimensions of an instrument from those given in the descriptions however, the instruments are inspected to ensure that any variations are within agreed tolerances. Unisurge operates active Quality Control / Assurance to fully comply with the regulatory requirements.

All instruments are manufactured from stainless steel and have been finished to the highest standards.

Unisurge is committed to achieving full compliance relating to quality and labour standards.



needle holder features and benefits

Tips are uniformly smooth and carefully rounded to avoid catching

Jaw serrations are precision cut. The diamond mesh enables a firm grip

Box joint is thoroughly cleaned and polished for a smooth closure

Satin finish for reduced glare

Shanks are slim and light weight

Each instrument is marked with a code number for more convenient re-ordering

Finger rings are generously proportioned and subtly shaped for maximum comfort and balance

Black handles remind you that these are the very best quality single-use instruments



scissors



features and benefits

Tips are either pointed or rounded

Edges are precision ground to ensure an even cut every time

Screw joint surfaces are polished for excellent cutting action

Satin finish for reduced glare

Shanks are slim and lightweight

Each instrument is marked with a code number for more convenient re-ordering

Finger rings are generously proportioned and subtly shaped for maximum comfort and balance

Black handles remind you that these are the very best quality single-use instruments



forceps

features and benefits

Tips are uniformly smooth and carefully rounded to avoid catching

Box joint is thoroughly cleaned and polished for a smooth closure

Satin finish for reduced glare

Shanks are slim and lightweight

Each instrument is marked with a code number for more convenient re-ordering


Finger rings are generously proportioned and subtly shaped for maximum comfort and balance

Black handles remind you that these are the very best quality single-use instruments



retractor

features and benefits



All retractors have smooth rounded edges (except teeth where sharp) optimally sized for surgical access

Ratchet joint for fine finger control allowing acute placement of jaws with a locking mechanism

Satin finish for reduced glare

Each instrument is marked with a code number for more convenient re-ordering

Finger rings are generously proportioned and subtly shaped for maximum comfort and balance

Black handles remind you that these are the very best quality single-use instruments



laparoscopic scissors

features and benefits

Tips are rounded

Edges are precision ground to ensure an even cut every time


Fully insulated to the working tip of the blade

Infinite rotational adjustment for optimum tip control

Ergonomic handles


Range of working lengths – 17, 33 & 44cm

Universal bayonet connector



laparoscopic graspers & dissectors

features and benefits



Comprehensive range of jaw configurations

Available with and without ratchet

Graduated intelligent ratchet action



Fully insulated

Ergonomic handles

Infinite rotational adjustment for optimum tip control

Range of working lengths –
17, 33 & 44cm

Universal bayonet connector



Signature™



other information



Single use. Do not re-use



Lot/Batch number



Date of manufacture

REF

Replaces the use of the words such as code, reference, article number



Sterile. Sterilised by Ethylene Oxide.



Read instructions for use and/or warnings included



Shows "use-by" date for the product as YYYY-MM



Latex-free

contents

ANAESTHETIC

15

CARDIO THORACIC

21

GENERAL

31

Gynaecological, Genito Urinary, Laparoscopy, Obstetrics

NEUROSURGERY

79

ENT, Maxillo Facial, Dental

OPHTHALMIC

97

ORTHOPAEDIC

107

Podiatry

PLASTIC

121

Introducing our Anaesthetic Range

 **Signature**
ANAESTHETIC

Fine Surgical Instruments

Type	Ref No	Description	Case Quantity
Dilators			
	Tracheal		
	F894025	Bowlby Child 10.5cm (4.5")	2 X 25
	F894024	Bowlby Adult 14cm (5.5")	2 X 25
Forceps			
	Intubation		
	F891123	Magill Adult 25cm (10")	50
	F891122	Magill Child 20cm (8")	50
	Dressing		
	F893093	Lister Sinus 15cm (6")	12 X 20
Laryngoscope Blades			
	Fibre Optic		
	F891013	Macintosh No 1	2 X 25
	F891014	Macintosh No 2	2 X 25
	F891015	Macintosh No 3	2 X 25
	F891016	Macintosh No 4	2 X 25
	F891017	Miller No 0	2 X 25
	F891018	Miller No 1	2 X 25
	F891019	Miller No 2	2 X 25
	F891020	Miller No 3	2 X 25
	LED		
	F891001	Macintosh No 1	2 X 25
	F891002	Macintosh No 2	2 X 25
	F891003	Macintosh No 3	2 X 25
	F891004	Macintosh No 4	2 X 25
	F891005	Miller No 0	2 X 25
	F891006	Miller No 1	2 X 25
	F891007	Miller No 2	2 X 25
	F891008	Miller No 3	2 X 25



Type	Ref No	Description	Case Quantity
Laryngoscope Handles			
	Fibre Optic		
	F891030	Standard	2 X 20
	F891032	Pen Light	2 X 20
	LED		
	F891029	Standard	2 X 20
	F891031	Pen Light	2 X 20
Suckers			
	F891040	Yankauer Rose Tip 28cm (11.2")	100
	F894133	Magill Small	8 X 10

Anaesthetic Instrument Sets

Anaesthetic Instrument Sets

I.V Cutdown Set

Ref No	F881001
UOM	10
Contains:	
Forceps Artery Halstead Mosquito Straight 12.5cm (5")	3
Forceps Artery Halstead Mosquito Curved 12.5cm (5")	3
Forceps Adson Dissecting Fine Toothed 12.5cm (5")	1
Forceps Adson Dissecting Fine Non Toothed 12.5cm (5")	1

Scissors Iris Straight 11cm (4.4")	1
Hook Skin Single (Blunt) 15cm (6")	1
Retractor Hook Skin Gillies Large 18cm (7.2")	1
Needle Aneurysm Small 15cm (6")	1
Needle Holder Kilner 13cm (5.2")	1
Handle BP No 3	1
Tray	1
Sterile Field	1

Tracheostomy Set

Ref No	F884001
UOM	10
Contains:	
Forceps Sponge Holding Rampley 24cm (9.5")	2
Clip Towel Cross Action	5
Pin Instrument Mayo 15cm (6")	2
Forceps Dissecting McIndoe Non Toothed 15cm (6")	1
Forceps Dissecting Gillies Toothed 15cm (6")	1
Forceps Artery Halstead Mosquito Straight 12.5cm (5")	5
Forceps Artery Dunhill Curved 12.8cm (5")	5
Forceps Tissue Allis 2:3 Teeth 15cm (6")	1
Handle BP No 3	1

Handle BP No 4	1
Needle Holder Masson 27cm (10.5")	1
Scissors Dressing S/S Straight 12.5cm (5")	1
Scissors Strabismus Straight Blunt 11.5cm (5.6")	1
Scissors McIndoe C.O.F 18cm (6")	1
Hook Skin Single (Blunt) 15cm (6")	1
Gillies Skin Hook 18cm (7")	1
Hook Skin Double (Blunt) 15cm (6")	1
Retractor West Blunt 3:4 Teeth 14cm (5.5")	1
Dilator Tracheal Bowlby Adult	1
Tray	1
Sterile Field	1

Introducing our
Cardio Thoracic Range



Signature
CARDIO THORACIC
Fine Surgical Instruments

Type	Ref No	Description	Case Quantity
Cannulas			
	F893016	Heparin Small Plain End 2mm	12 X 20
	F893017	Heparin Medium Plain End 2.5mm	12 X 20
Cutters			
	Wire		
	F896007	Universal Wire 12cm (4.8")	12 X 10
Forceps			
	Clip		
	F893007	DeBakey Bulldog Small Straight	12 X 20
	F893005	DeBakey Bulldog Medium Straight	12 X 20
	F893004	DeBakey Bulldog Large Straight	12 X 20
	F893006	DeBakey Bulldog Small Curved	12 X 20
	F893008	Dieffenbach Bulldog	12 X 20
	Clamp		
	F893009	Tubing Straight 12.5cm (5")	12 X 20
	Artery		
	F893032	Heiss Straight 20cm (8")	8 X 10
	F893035	O'Shaughnessy 20cm (8")	8 X 10
	F893033	Roberts Curved 23cm (9")	8 X 10
	F892005	Pean Straight 23cm (9")	8 X 10
	F892006	Pean Curved 23cm (9")	8 X 10
	F892007	Pean Curved 25.5cm (10")	8 X 10
	F892008	Pean Curved 30cm (12")	8 X 10
	Ligature		
	F893091	Lahey 20cm (8")	8 X 10
	F893092	Lahey 23cm (9")	8 X 10
	F893090	O'Shaughnessy 22cm (8.8")	50
	Tissue		
	F893102	Duval 23cm (9")	8 X 10

Type	Ref No	Description	Case Quantity
Forceps (Continued)			
	Dissecting		
	F893063	DeBakey 18cm (7.2")	8 X 10
	F893085	DeBakey 25cm (10")	50
	F893086	DeBakey 30cm (12")	50
Needle Holders			
	F893118	Converse 12.5cm (5")	12 X 20
	F893117	Crile Wood 15cm (6")	12 X 20
	F893119	Crile Wood 20cm (8")	8 X 10
	F893120	DeBakey 18cm (7")	8 X 10
	F893121	Derf 12cm (4.75")	12 X 20
	F893122	Foster with R/H Scissors 12.5cm (5")	12 X 20
	F893124	Halsey 12.5cm (5")	12 X 20
	F893125	Halsey Tungsten Carbide 12.5cm (5")	12 X 20
	F893137	Nievert 13cm (5.2")	8 X 10
Retractors			
	Chest		
	F892001	Finochietto Small	5
	F892002	Finochietto Medium	5
	F892003	Finochietto Large	5
	F892004	Tudor Edwards Rib Spreader 23cm (9")	5
Scissors			
	Straight		
	F893205	Nelson Straight 22.5cm (9")	8 X 10
	F893324	Nelson Curved 22.5cm (9")	8 X 10
	Angled		
	F893206	Potts Angled on Flat 45 Degree 19cm (7.5")	8 X 10
	F893207	Potts Angled on Side 25 Degree 19cm (7.5")	8 X 10
	F893208	Potts Angled on Side 45 Degree 19cm (7.5")	8 X 10

Cardio Thoracic Instrument Sets



Cardio Thoracic Instrument Sets

Thoracotomy Basic Set

Ref No	F882001	Handle BP No 5 Graduated	1
UOM	5	Needle Holder Thompson Walker 20cm (8")	2
Contains:		Scissors Metzenbaum C.O.F 23cm (9")	2
Forceps Sponge Holding Rampley 24cm (9.5")	3	Scissors Nelson Curved 22.5cm (9")	1
Pin Instrument Mayo 15cm (6")	2	Allison Lung Retractor 32.5cm (13")	1
Forceps Dissecting Bonney Toothed 18cm (7")	1	Tudor Edwards Rib Spreader 23cm (9")	1
Moynihan Cholecystectomy 22.5cm (9")	2	Tray	1
Forceps Artery Roberts Curved 22cm (8.8")	10	Sterile Field	1
Forceps Tissue Duval 23cm (9.2")	4		

Cardio Thoracic Instrument Sets (Continued)

Cardiovascular Set

Ref No	F882002		
UOM	5		
Contains:			
Forceps Sponge Holding Rampley 24cm (9.5")	8	Scissors Potts Angled On Side 25 Degree	1
Clip Towel Cross Action	12	Scissors Potts Angled On Side 45 Degree	1
Handle BP No 3	3	Scissors Mayo Straight 17.5cm (7")	2
Handle BP No 4	1	Scissors Mayo Curved 16.5cm (6.5")	2
Cutter Universal Wire 12cm (4.8")	2	Scissors Metzenbaum C.O.F 19cm (7.5")	2
Forceps Dissecting Adson Toothed 12.5cm (5")	1	Forceps Tissue Allis 2:3 Teeth 15cm (6")	3
Forceps Dissecting Adson Non Toothed 12.5cm (5")	1	Retractor Travers Self Retaining 20cm (8")	2
Forceps Dissecting Adson Browne Multi Toothed 12cm (4.7")	2	Forceps Artery Kelly Curved 14cm (5.5")	6
Forceps Dissecting DeBakey 18cm (7.2")	4	Retractor Deaver 21.5cm Blade 2.5cm (1")	2
Forceps Dissecting DeBakey Atraumatic 25.5cm (10")	2	Retractor Deaver 25cm Blade 2.5cm (1")	1
Forceps Tissue Allis 3:4 Teeth 15cm (6")	6	Retractor Deaver 33cm Blade 3.8cm (1.5")	2
Forceps Tissue Babcock 16cm (6.7")	4	Retractor Finochietto Large	1
Forceps Artery Halstead Mosquito Straight 12.5cm (5")	1	Retractor Self Retaining Weitlaner 3/4 Sharp Teeth 14cm (5.5")	2
Forceps Artery Halstead Mosquito Curved 12.5cm (5")	6	Suction Tube Yankauer Rose Tip 28cm (11.2")	2
Forceps Artery Crile Straight 14cm (5.5")	2	Suction Tube Frazier 12fg	1
Forceps Artery Crile Curved 14cm (5.5")	12	Suction Tube Pooles Straight 24cm (9.5")	1
Scissors Metzenbaum C.O.F 23cm (9")	1	Needle Holder Mayo Hegar 14cm (5.6")	3
		Needle Holder Mayo 20cm (8")	3
		Needle Holder Mayo Hegar 18cm (7")	5
		Hook Adson Dissecting Sharp 19.5cm (7.7")	2
		Hook Adson Dissecting Blunt 19.5cm (7.7")	1
		Tray	1
		Sterile Field	1

Cardio Thoracic Procedure Packs

Cardio Thoracic Procedure Packs

Abrams Pleural Biopsy Pack

Ref No	F200186
UOM	10
Contains:	
Gallipot Polypropylene 60ml	1
Swab Non Woven 7.5 x 7.5cm 4ply	5
Scalpel Blade No 10	1
Suture Mersilk W775	1
Needle Hypodermic 21g x 1.5" Green	2
Needle Hypodermic 23G x 1.25" Blue	1
Syringe 10ml Luer Slip	1
Syringe 5ml Luer Slip	1
Bag Yellow Poly 50 x 38cm with Adhesive Strip on reverse	1
Gloves Exam Nitrile, Powder Free, Cuffed, Walleded, M/L, Pair	1
Towel Dressing 43 x 38cm	1
Paper Crepe White 90 x 90cm	1
Drape Fenestrated 60 x 75cm, with Off Centre Adhesive Fenestration 10cm	1



Drape Gown Extra Large	1
Needle Abrams Pleural Biopsy	1
Forceps Artery Spencer Wells Curved 12.5cm (5")	1
Forceps Dissecting T.O.E. Serrated 12.5cm (5")	2
Needle Holder Kilner 13cm (5.2")	1
Handle BP No 3 Graduated	1
Tray	1
Sterile Field	1



Cardio Thoracic Procedure Packs (Continued)



Chest Drain Pack

Ref No	F200187
UOM	10
Contains:	
Gallipot Polypropylene 60ml	1
Swab Gauze 10 x 10cm 8ply	5
Scalpel Blade No 11	1
Suture Mersilk W775	1
Needle Hypodermic 21g x 1.5" Green	2
Needle Hypodermic 23g x 1.25" Blue	1
Syringe 10ml Luer Slip Eccentric	1
Syringe 5ml Luer Slip	1
Bag Yellow Poly 45 x 29cm Adhesive Strip on reverse	1
Gloves Exam Nitrile, Powder Free, Cuffed, Walled, M/L Pair	1
Towel Dressing 43 x 38cm	1
Drape Adhesive 60 x 75cm with Off Centre Adhesive Fenestration 10cm	1



Drape Gown Extra Large	1
Forceps Artery Spencer Wells Straight 23cm (9.2")	1
Forceps Artery Roberts Curved 22cm (8.8")	1
Needle Holder Mayo Hegar 14cm (5.6")	1
Forceps Artery Halstead Mosquito Curved 12.5cm (5")	1
Scissors Mayo Straight 17.5cm (7")	1
Tray	1
Sterile Field	1

Cardio Thoracic Procedure Packs

Cardio Thoracic Procedure Packs (Continued)

Pacemaker Pack

Ref No	F200194
UOM	4
Contains:	
Bowl Polypropylene 1000ml	1
Dish Kidney Polypropylene 10" (800ml)	2
Gallipot Polypropylene 60ml	2
Bowl Polypropylene 250ml	1
Swab Gauze XRD 10 x 7.5cm 32ply	10
Needle Hypodermic 19g x 1.5" White	1
Needle Hypodermic 21g x 1.5" Green	1
Needle Hypodermic 25g x 1" Orange	1
Needle Hypodermic 23g x 1" Blue	1
Suture Mersilk W667	1
Suture Vicryl W9719	1
Suture Vicryl W9120	1
Syringe 2ml Luer Slip	1
Syringe 20ml Luer Slip Eccentric	1
Syringe 10ml Luer Slip Eccentric	1
Labels Drug Lidocaine	1
Labels Drug Gentamicin	1
Foam Cubes White 5 x 3 x 2.5cm	1
Scissors Iris Micro Eye Straight 11.25cm (4.5")	1
Forceps Artery Birketts Curved 18cm (7.2")	1
Forceps Artery Halstead Mosquito Curved 12.5cm (5")	3
Scissors Dressing S/S Straight 12.5cm (5")	1
Scissors McIndoe Curved 18cm (7")	1
Retractor Travers Self Retaining 20cm (8")	1



Clamp Towel	6
Dressing Mepore 9 x 15cm	1
Sponge Stick Orange 20cm (8")	2
Basin Liner Single, Folded	1
Bag Plain Top L 250 x 100 x 380mm	1
Dressing Pad 20 x 20cm	1
Towel Scrimex 48 x 40cm	4
Drape Fluoroscope 91 x 51cm	1
Drape Medium Utility N/A 90 x 90cm	3
Drape "U" 193 x 304cm (10 x 101cm Split) & Control Panel	1
Gown Large No Reinforcement	2
Scalpel BP Handle With No 15 Blade	1
Forceps Dissecting Adson Non Toothed 12.5cm (5")	1
Needle Holder Kilner 13cm (5.2")	1
Forceps Dissecting Adson Toothed 12.5cm (5")	1
Forceps Artery Halstead Mosquito Straight 12.5cm (5")	3
Tray	1
Sterile Field	1

Introducing our General Range incorporating:-

- Gynaecological
- Genito Urinary
- Laparoscopy
- Obstetrics



Fine Surgical Instruments

Type	Ref No	Description	Case Quantity
BP Handles			
	F893002	Mini Blade	12 X 20
	F893110	Scalpel Barron	12 X 20
	F893107	Scalpel BP No 3	12 X 20
	F893108	Scalpel BP No 4	12 X 20
	F893111	Scalpel BP No 5	12 X 20
	F893112	Scalpel BP No 7	12 X 20
Dilators			
	F893018	Umbilical	12 X 20
Directors			
	F893019	Fistula With Guard 13cm (5.2")	12 X 20
	F893020	Lockhart Mummery Set of 4	2 X 25
Dissectors			
	F893021	MacDonald Double Ended 19cm (7.5")	8 X 10
	F896008	Watson Cheyne Elevator 19cm (7.5")	8 X 10
Forceps			
	F893087	Epilation 8.75cm (3.5")	12 X 20
	F893088	Epilation 10cm (4")	12 X 20
	F893093	Lister Sinus 15cm (6")	12 X 20
	F893009	Tubing Straight 12.5cm (5")	12 X 20
	F893106	Vasectomy Straight 15cm (6")	8 X 20
	F893109	Vasectomy Straight Blunt Modified 12.5cm (5")	12 X 20
	Clip		
	F893011	Michel Clip Applying & Removing 12cm (4.7")	2 X 50
	F893296	Michel Clip Removing 12cm (4.7")	2 X 100
	Artery		
	F893024	Adson Straight 18.5cm (7.4")	8 X 20
	F893305	Cairn 14cm (5.5")	12 X 20

Type	Ref No	Description	Case Quantity
Forceps (Continued)			
	Artery (Continued)		
	F893025	Crile Curved 14cm (5.5")	12 X 20
	F893026	Crile Straight 14cm (5.5")	12 X 20
	F893027	Dunhill Curved 12.5cm (5")	12 X 20
	F893028	Dunhill Straight 12.5cm (5")	12 X 20
	F893044	Grey Turner Curved 14.5cm (5.75")	8 X 10
	F893045	Grey Turner Curved 18cm (7")	8 X 10
	F893029	Halstead Mosquito Curved 12.5cm (5")	12 X 20
	F893030	Halstead Mosquito Straight 12.5cm (5")	12 X 20
	F893032	Heiss Straight 20cm (8")	8 X 10
	F893104	Kocher Toothed Straight 14cm (5.5")	8 X 10
	F893031	Kocher Toothed Curved 14cm (5.5")	8 X 10
	F896009	Kocher Toothed Straight 20cm (8")	8 X 10
	F893291	Kocher Toothed Curved 20cm (8")	8 X 10
	F893091	Lahey Ligature 20cm (8")	8 X 10
	F893092	Lahey Ligature 23cm (9")	8 X 10
	F893307	Moynihan Cholecystectomy 22.5cm (9")	8 X 10
	F893033	Roberts Curved 22.5cm (9")	8 X 10
	F893038	Spencer Wells Curved 12.5cm (5")	12 X 20
	F893046	Spencer Wells Curved 15cm (6")	8 X 10
	F893039	Spencer Wells Curved 18cm (7")	8 X 10
	F893047	Spencer Wells Curved 20cm (8")	8 X 10
	F893040	Spencer Wells Straight 12.5cm (5")	12 X 20
	F893041	Spencer Wells Straight 15cm (6")	8 X 10
	F893042	Spencer Wells Straight 18cm (7.2")	8 X 10
	F893043	Spencer Wells Straight 19cm (7.6")	8 X 10
	F893036	Spencer Wells Straight 20cm (8")	8 X 10
	F893037	Spencer Wells Straight 23cm (9.2")	8 X 10
	Sponge		
	F893094	Rampléy Sponge Holding 18cm (7")	8 X 10
	F893095	Rampléy Sponge Holding 24cm (9.5")	8 X 10

Type	Ref No	Description	Case Quantity
Forceps (Continued)			
	Tissue		
	F893096	Allis 4:5 Teeth 15cm (6")	8 X 20
	F893099	Allis American Pattern 14.5cm (5.75")	12 X 20
	F893097	Allis Toothed 2:3 15cm (6")	8 X 20
	F893098	Allis Toothed 3:4 15cm (6")	8 X 20
	F893298	Allis Toothed 4:5 26cm (10")	50
	F893290	Allis Toothed 5:6 20cm (8")	8 X 10
	F893100	Babcock 16.5cm (6.5")	8 X 20
	F893105	Babcock 20cm (8")	8 X 10
	F893102	Duval 23cm (9")	8 X 10
	F893103	Judd-Allis 3:4 Teeth 19.5cm (7.75")	8 X 10
	F893101	Lanes 17.5cm (7")	8 X 10
	F893253	Littlewood Toothed 2:3 18cm (7")	8 X 10
	F893089	Patterson Grasping 50cm (20")	5
	Towel		
	F893012	Backhaus 11cm (4.4")	12 X 20
	F893013	Ball and Socket	12 X 20
	F893014	Cross Action	12 X 10
	Dissecting		
	F893287	Adson Fine Non Toothed 12.5cm (5")	12 X 20
	F893286	Adson Fine Toothed 12.5cm (5")	12 X 20
	F893059	Adson Non Toothed 12.5cm (5")	12 X 20
	F893060	Adson Toothed 1:2 Bayonet 11cm (4.5")	12 X 10
	F893058	Adson Toothed 1:2 Bayonet 19cm (7.5")	8 X 10
	F893057	Adson Toothed 12.5cm (5")	12 X 20
	F893064	Gillies Toothed 15cm (6")	8 X 20
	F893084	Jean Non Toothed 13.75cm (5.5")	12 X 20
	F893067	Lane Non Toothed 18cm (7")	8 X 10
	F893069	Lane Toothed 1:2 Teeth 12.5cm (5.2")	12 X 10
	F893068	Lane Toothed 1:2 Teeth 15cm (6")	12 X 10
	F893070	Lane Toothed 1:2 Teeth 18cm (7")	8 X 10

Type	Ref No	Description	Case Quantity
Forceps (Continued)			
	Dissecting (Continued)		
	F893071	McIndoe Non Toothed 15cm (6")	12 X 20
	F893302	McIndoe Toothed 1:2 15cm (6")	8 X 20
	F893074	Ramsey Toothed 1:2 Teeth 18cm (7")	8 X 20
	F893077	T.O.E 12.5cm (5")	12 X 20
	F893073	Treves Non Toothed 12.5cm (5")	8 X 20
	F893081	Treves Toothed Straight 12.5cm (5")	12 X 20
	F893080	Treves Toothed 15cm (6")	12 X 20
	F893082	Waugh Non Toothed 18cm (7")	8 X 10
	F893083	Waugh Toothed 1:2 20cm (8")	8 X 10
	Biopsy		
	F893061	Lloyd Davies Biopsy 35cm (14")	5
	F893284	Lloyd Davies Biopsy 50cm (20")	5
	F893062	Patterson Biopsy Shaft 35cm (14")	5
	F893052	Patterson Biopsy 50cm (20")	5
	Diathermy		
	F893279	Bipolar 11.25cm (4.5") and 3m Lead	50
	F893280	Bipolar 15cm (6") and 3m Lead	50
	F893281	Bipolar 20cm (8") and 3m Lead	50
	Splinter		
	F893048	Arthurs	12 X 20
	F893049	Hunter	12 X 20
	F893050	Martins	12 X 20
	Stone		
	F893318	Desjardins Grasping Small 23cm (9")	25
	F893317	Desjardins Grasping Large 23cm (9")	25
Mirrors			
	Dental		
	F894042	No 4	2 X 50
	F894043	No 5	2 X 50

Type	Ref No	Description	Case Quantity
Needles			
	F893114	Aneurysm Small 15cm (6")	8 X 10
	F893115	Abrams Pleural Biopsy	8 X 10
	F893116	Post Mortem Half (4.4")	12 X 20
Needle Holders			
	F893118	Converse 12.5cm (5")	12 X 20
	F893117	Crile Wood 15cm (6")	12 X 20
	F893119	Crile Wood 20cm (8")	8 X 10
	F893122	Foster with R/H Scissors 12.5cm (5")	12 X 20
	F893123	Gillies 16cm (6.25")	8 X 10
	F893124	Halsey 12.5cm (5")	12 X 20
	F893125	Halsey Tungsten Carbide 15cm (6")	12 X 20
	F893126	Kilner 13.5cm (5.25")	12 X 20
	F893127	Lawrence 15cm (6")	8 X 20
	F893320	Lloyd Davis 25cm (10")	8 X 10
	F893308	Masson 27cm (10.5")	8 X 10
	F893128	Mayo 18cm (7")	8 X 10
	F893129	Mayo 20cm (8")	8 X 10
	F893130	Mayo 22.5cm (9")	8 X 10
	F893131	Mayo Hegar 13cm (5.2")	12 X 20
	F893132	Mayo Hegar 14cm (5.6")	12 X 20
	F893133	Mayo Hegar 16cm (6.4")	8 X 10
	F893134	Mayo Hegar 18cm (7.2")	8 X 10
	F893136	Mayo Hegar Fine 18cm (7.2")	8 X 10
	F893137	Nievert 13cm (5.2")	8 X 10
	F893139	Thomson Walker 20cm (8")	8 X 10
	F893140	Wright 14cm (5.6")	12 X 20
Pins			
	Mayo		
	F893142	Mayo Instrument Pin 15cm (6")	8 X 10
Probes			
	F893143	Nickel Silver With Eye 12.5cm (5")	12 X 20
	F893146	Nickel Silver With Eye 15cm (6")	12 X 20
	F893144	Wards Double Ended 15cm (6")	8 X 10

Type	Ref No	Description	Case Quantity
Proctoscopes			
	F893147	Gabriel Adult 6.4cm (2.5")	20
	F893148	Kelly 14cm (5.5")	20
	F893149	Naunton Morgan Rectal 16cm (6.5")	12 X 20
Retractors			
	Hook		
	F893156	Gillies Skin 18cm (7")	8 X 20
	F893157	McIndoe Skin 19cm (7.6")	8 X 10
	F893155	Skin Single (Blunt) 15cm (6")	8 X 10
	F893153	Skin Double (Blunt) 15cm (6")	8 X 10
	Hand Held		
	F893309	Allison Lung Retractor 32.5cm (13")	5
	F893151	Czerny 17.5cm (7")	25
	F893163	Deaver 21.5cm (8.5") Blade 2.5cm (1")	20
	F893172	Deaver 25cm (10") Blade 2.5cm (1")	20
	F893152	Deaver 30cm (12") Blade 3.8cm (1.5")	20
	F893162	Kelly 21.5cm (8.5") Blade 3.8cm (1.5")	20
	F893161	Kilner Catspaw Double Ended 15cm (6")	8 X 10
	F893164	Langenbeck Small 21.5cm (8.5") Blade 25 x 6mm	8 X 10
	F893165	Langenbeck Medium 21.5cm (8.5") Blade 35 x 11mm	50
	F893321	Langenbeck Medium Narrow 21.5cm (8.5") Blade 35 x 8mm	50
	F893322	Langenbeck Medium Wide 21.5cm (8.5") Blade 35 x 15mm	50
	F893166	Langenbeck Large 21.5cm (8.5") Blade 41 x 11mm	50
	F893323	Langenbeck X-Large Long 21.5cm (8.5") Blade 55 x 11mm	50
	F893313	Lloyd Davis 29cm (11.5")	5
	F893167	Morris 24.5cm (9.75") Blade Small 70 x 50mm	50
	F893325	Morris 24.5cm (9.75") Blade Medium 70 x 65mm	50
	F894143	Rake 17cm (6.7") long 5.5cm head	8 X 5
	Self Retaining		
	F893310	Balfour Abdominal	5
	F893158	Jackson Burrows 15.6cm (6.5")	20
	F893312	Joll Thyroid 15cm (6")	50
	F893170	Travers 4:5 Blunt 20cm (8")	25

Type	Ref No	Description	Case Quantity
Retractors (Continued)			
	Self Retaining (Continued)		
	F893168	Weitlaner 3:4 Blunt Teeth 14cm (5.5")	2 X 50
	F893169	Weitlaner 3:4 Sharp Teeth 14cm (5.5")	2 X 50
	F893173	West Blunt 3:4 Teeth 14cm (5.5")	50
Scissors			
	Dressing		
	F893178	S/B Straight 12.5cm (5")	12 X 20
	F893179	S/S Straight 12.5cm (5")	12 X 20
	F893176	B/B Straight 12.5cm (5")	12 X 20
	F893180	S/S Straight 20cm (8")	8 X 10
	F893177	B/B Straight 20cm (8")	8 X 10
	F893185	Lister Bandage 14cm (5.5")	12 X 10
	F893187	Lister Bandage 15cm (6")	8 X 10
	F893186	Lister Bandage 18.75cm (7.5")	8 X 10
	Dissecting		
	F893189	Lloyd Davies Rectal Curved 25cm (10")	5
	F893190	Mayo C.O.F 16.5cm (6.5")	8 X 10
	F893191	Mayo C.O.F 19cm (7.5")	8 X 10
	F893192	Mayo C.O.F 23cm (9")	8 X 10
	F893194	Mayo Straight 14cm (5.5")	8 X 20
	F893195	Mayo Straight 15cm (6")	8 X 20
	F893196	Mayo Straight 16.5cm (6.5")	8 X 10
	F893213	Mayo Straight 19cm (7.5")	8 X 10
	F893197	Mayo Straight 23cm (9")	8 X 10
	F893198	McIndoe C.O.F 12.5cm (5")	12 X 20
	F893199	McIndoe C.O.F 18cm (6")	8 X 10
	F893200	McIndoe C.O.F 19cm (7.5")	8 X 10
	F893301	McIndoe C.O.F Tungsten Carbide 23cm (7")	8 X 10
	F893211	Metzenbaum C.O.F 18cm (7")	8 X 10
	F893202	Metzenbaum C.O.F 19cm (7.5")	8 X 10
	F893203	Metzenbaum C.O.F 23cm (9")	8 X 10
	F893300	Metzenbaum C.O.F Tungsten Carbide 18cm (7")	8 X 10
	F893215	Metzenbaum Fine Curved 15cm (6")	8 X 10

Type	Ref No	Description	Case Quantity
Scissors (Continued)			
	Dissecting (Continued)		
	F893201	Metzenbaum Straight 14cm (5.5")	8 X 10
	F893292	Metzenbaum Straight 23cm (9")	8 X 10
	F895048	Strabismus Curved Blunt 11.5cm (4.6")	12 X 20
	F895049	Strabismus Straight Blunt 11.5cm (4.6")	12 X 20
	Cutter		
	F896007	Universal Wire 12cm (4.8")	12 X 10
Suckers			
	F893218	Pooles 24cm (9.5")	100
Syringes			
	Dental		
	F894061	Dental Cartridge	8 X 10
Trocars & Cannulas			
	F893219	Nelson 5.5mm	8 X 10

General Instrument Sets



General Instrument Sets

General Set

Ref No	F883019	Handle BP No 3	2
UOM	5	Handle BP No 4	1
Contains:		Needle Holder Mayo Hegar 14cm (5.6")	1
Forceps Sponge Holding Rampley 24cm (9.5")	5	Needle Holder Mayo Hegar 18cm (7")	2
Clip Towel Cross Action	5	Scissors Mayo Straight 14cm (5.5")	2
Pin Instrument Mayo 15cm (6")	3	Scissors Mayo C.O.F. 19cm (7.6")	1
Forceps Dissecting T.O.E. Serrated 12.5cm (5")	2	Scissors McIndoe Curved 18cm (7")	1
Forceps Dissecting Lane Toothed 1:2 Teeth 12.5cm (5")	1	Scissors Dressing S/S Straight 12.5cm (5")	1
Forceps Dissecting Lane Toothed 1:2 Teeth 18cm (7")	1	Needle Aneurysm Small 15cm (6")	1
Forceps Dissecting Gillies Toothed 15cm (6")	1	Probe Nickel Silver 15cm (6")	1
Forceps Dissecting McIndoe Non Toothed 15cm (6")	1	Curette Volkmann Small D/E 22cm Size A 22cm (8.6")	1
Forceps Artery Dunhill Curved 12.8cm (5")	10	Suction Tube Poole Straight 24cm (9.5")	1
Forceps Artery Spencer Wells Straight 18cm (7.2")	10	Forceps Artery Roberts Curved 22cm (8.8")	1
Forceps Artery Spencer Wells Curved 18cm (7.2")	2	Langenbeck Small 21.5cm (8.5") Blade 25 x 6mm	2
Forceps Moynihan Cholecystectomy 22.5cm (9")	2	Langenbeck Large 21.5cm (8.5") Blade 41 x 11mm	2
Forceps Tissue Allis 2:3 Teeth 15cm (6")	2	Retractor Morris 24.5cm (9.75") Blade Small 70 x 50mm	2
Forceps Tissue Babcock 16cm (6.7")	2	Retractor Travers Self Retaining 20cm (8")	1
Forceps Tissue Littlewood 2:3 Teeth 19cm (7.5")	2	Retractor Balfour Abdominal	1
		Tray	1
		Sterile Field	1



General Instrument Sets (Continued)

Large General Set

Ref No	F883020		
UOM	5		
Contains:			
Forceps Sponge Holding Rampley 24cm (9.5")	5	Forceps Dissecting McIndoe Non Toothed 15cm (6")	1
Clip Towel Cross Action	10	Scissors Dressing S/S Straight 12.5cm (5")	1
Pin Instrument Mayo 15cm (6")	5	Scissors McIndoe Curved 18cm (7")	1
Forceps Dissecting Lane Toothed 1:2 Teeth 18cm (7")	2	Scissors Nelson Curved 22.5cm (9")	1
Forceps Dissecting Waugh Non Toothed 20cm (8")	1	Needle Aneurysm Small 15cm (6")	1
Forceps Dissecting Waugh Toothed 20cm (8")	1	Probe Nickel Silver 15cm (6")	1
Forceps Dissecting Gillies Toothed 15cm (6")	1	Dissector Watson Cheyne 19cm (7.5")	1
Forceps Dissecting McIndoe Non Toothed 15cm (6")	1	Curette Volkmann Small D/E Size A 22cm (8.6")	1
Forceps Artery Halstead Mosquito Curved 12.5cm (5")	5	Dissector McDonald D/E 19cm (7.5")	1
Forceps Artery Spencer Wells Curved 18cm (7.2")	10	Forceps Dressing Lister Sinus 12.5cm (5")	1
Forceps Moynihan Cholecystectomy 22.5cm (9")	5	Suction Tube Poole Straight 24cm (9.5")	1
Forceps Tissue Allis 2:3 Teeth 15cm (6")	2	Langenbeck Large 21.5cm (8.5")	2
Forceps Tissue Babcock 16cm (6.7")	2	Retractor Morris Small 2 x 1.5cm Blade, 9.75cm (3.8")	1
Forceps Tissue Duval 23cm (9.2")	2	Retractor Kelly Blade 38mm 21.5cm (1.5")	2
Forceps Tissue Lane 17.5cm (7")	2	Retractor Balfour Abdominal	1
Handle BP No 4	2	Retractor Travers Self Retaining 20cm (8")	1
Handle BP No 3	1	Forceps Artery Crile Curved 14cm (5.5")	5
Needle Holder Mayo 20cm (8")	1	Forceps Artery Cairn Curved 14cm (5.5")	5
Scissors Mayo Straight 17.5cm (7")	2	Needle Holder Mayo Hegar 15cm (6")	1
Scissors Mayo Curved 16.5cm (6.5")	1	Scissors Metzenbaum C.O.F 19cm (7")	1
		Diathermy Bipolar Forceps Imm Tip & 3m Lead 15cm (6")	1
		Diathermy Bipolar Forceps Imm Tip & 3m Lead 20cm (8")	1
		Tray	1
		Sterile Field	1

General Instrument Sets



General Instrument Sets (Continued)

Minor Set

Ref No	F883021	Handle BP No 3	1
UOM	5	Handle BP No 4	1
Contains:		Needle Holder Mayo Hegar 18cm (7")	1
Forceps Sponge Holding Rampley 24cm (9.5")	2	Scissors Dressing S/S Straight 12.5cm (5")	1
Clip Towel Backhaus 11cm (4.5")	5	Scissors Mayo Straight 17.5cm (7")	1
Pin Instrument Mayo 15cm (6")	2	Scissors McIndoe Curved 18cm (7")	1
Forceps Adson Dissecting Fine Non Toothed 12.5cm (5")	1	Forceps Dressing Lister Sinus 12.5cm (5")	1
Forceps Adson Dissecting Fine Toothed 12.5cm (5")	1	Needle Aneurysm Small 15cm (6")	1
Forceps Dissecting McIndoe Non Toothed 15cm (6")	2	Scissors Iris Straight 11cm (4.4")	1
Forceps Dissecting Gillies Toothed 15cm (6")	2	Curette Volkmann Small D/E Size A 22cm (8.6")	1
Forceps Artery Crile Curved 14cm (5.5")	5	Probe Nickel Silver 15cm (6")	1
Forceps Artery Halstead Mosquito Curved 12.5cm (5")	5	Retractor West Blunt 3:4 Teeth 14cm (5.5")	1
Forceps Tissue Allis 2:3 Teeth 15cm (6")	2	Retractor Hook Single Blunt	2
Forceps Tissue Babcock 16cm (6.3")	2	Retractor Kilner Catspaw D/E 15cm (6")	2
Forceps Tissue Lane 17.5cm (7")	2	Retractor Langenbeck Small 21.5cm (8.5")	2
		Blade 25 x 6mm	
		Tray	1
		Sterile Field	1

Intestinal Set

Ref No	F883022	Forceps Moynihan Cholecystectomy 22.5cm (9")	2
UOM	5	Forceps Ligature Lahey 20cm (8")	2
Contains:		Forceps Tissue Judd-Allis 3:4 Teeth 19cm (7.6")	4
Forceps Dissecting Waugh Non Toothed 20cm (8")	1	Needle Holder Lloyd Davis	1
Forceps Dissecting Waugh Toothed 20cm (8")	1	Scissors Nelson Curved 22.5cm (9")	2
Forceps Dissecting DeBakey 25cm (10")	1	Retractor Blade Deaver 33cm, Blade 3.8cm (1.5")	1
Forceps Artery Roberts Curved 22.5cm (9")	5	Tray	1
Forceps Artery Spencer Wells Straight 18cm (7.2")	10	Sterile Field	1



General Instrument Sets (Continued)

Suture Set

Ref No	F883026	Forceps Dissecting Non Toothed Treves 12.5cm (5")	
UOM	10	Scissors Dressing S/S Straight 12.5cm (5")	
Contains:		Needle Holder Kilner 13cm (5.2")	
Forceps Artery Halstead Mosquito Straight 12.5cm (5")		Tray	
Forceps Artery Halstead Mosquito Curved 12.5cm (5")		Sterile Field	
Forceps Dissecting T.O.E. Serrated 12.5cm (5")			

Skin Biopsy Set

Ref No	F887006	Forceps Dissecting T.O.E. Serrated 12.5cm (5")	
UOM	10	Hook Skin Single (Blunt) 15cm (6")	
Contains:		Scissor/Needle Holder Fosters R/H 12.5cm (5")	
Forceps Artery Halstead Mosquito Curved 12.5cm (5")	2	Scissors Iris Straight 11cm (4.4")	
Forceps Dissecting McIndoe Non Toothed 15cm (6")		Handle BP No 3	
Forceps Dissecting Gillies Toothed 15cm (6")		Tray	
		Sterile Field	

Gall Bladder Set

Ref No	F883024	Scissors Nelson Curved 22.5cm (9")	
UOM	10	Retractor Self Retaining Millin 80 x 25mm	
Contains:		Desjardins Stone Grasping Small 23cm (9")	
Forceps Dissecting Waugh Toothed 20cm (8")		Tray	
Forceps Dissecting Waugh Non Toothed 20cm (8")		Sterile Field	
Forceps Moynihan Cholecystectomy 22.5cm (9")	2		

General Procedure Packs

Abrams Pleural Biopsy Pack

Ref No	F200186
UOM	10
Contains:	
Gallipot Polypropylene 60ml	1
Swab Non Woven 7.5 x 7.5cm 4ply	5
Scalpel Blade No 10	1
Suture Mersilk W775	1
Needle Hypodermic 21g x 1.5" Green	2
Needle Hypodermic 23G x 1.25" Blue	1
Syringe 10ml Luer Slip	1
Syringe 5ml Luer Slip	1
Bag Yellow Poly 50 x 38cm with Adhesive Strip on reverse	1
Gloves Exam Nitrile, Powder Free, Cuffed, Walled, M/L, Pair	1
Towel Dressing 43 x 38cm	1
Paper Crepe White 90 x 90cm	1



Drape Fenestrated 60 x 75cm, with Off Centre Adhesive Fenestration 10cm	1
Drape Gown Extra Large	1
Needle Abrams Pleural Biopsy	1
Forceps Artery Spencer Wells Curved 12.5cm (5")	1
Forceps Dissecting T.O.E. Serrated 12.5cm (5")	2
Needle Holder Kilner 13cm (5.2")	1
Handle BP No 3 Graduated	1
Tray	1
Sterile Field	1

Biopsy Pack

Ref No	F883004
UOM	30
Contains:	
Forceps Artery Halstead Mosquito Straight 12.5cm (5")	1
Forceps Dissecting Adson Non Toothed 12.5cm (5")	1
Forceps Dissecting Adson Toothed 12.5cm (5")	1
Needle Holder Derf 12cm (4.8")	1
Scissors Dressing S/B Straight 12.5cm (5")	1
Scissors Kilner Curved 11cm (4.4")	1
Towel Dressing 43 x 38cm	1



Swab Non Woven 7.5 x 7.5cm 4ply	5
Paper Crepe White (Folded) 60 x 60cm	1
Scalpel BP Blue Handle With No 15 Blade	1
Tray	1
Sterile Field	1

General Procedure Packs (Continued)

Basic Minor Op Pack

Ref No	F883005
UOM	30
Contains:	
Forceps Artery Halstead Mosquito Straight 12.5cm (5")	
Forceps Dissecting Gillies Toothed 15cm (6")	
Forceps Dissecting McIndoe Non Toothed 15cm (6")	
Handle BP No 3 Graduated	
Needle Holder Kilner 13cm (5.2")	
Scissors Iris Curved 11.25cm (4.5")	
Scissors Iris Straight 11cm (4.4")	
Blade Carbon Steel No 15	



Towel Dressing 43 x 38cm	2
Swab Non Woven 7.5 x 7.5cm 4ply	5
Drape Fenestrated 60 x 75cm Off Centre S/A Fen 10cm	
Tray	
Sterile Field	

Standard Minor Op Pack

Ref No	F883006
UOM	18
Contains:	
Forceps Artery Halstead Mosquito Curved 12.5cm (5")	
Forceps Artery Halstead Mosquito Straight 12.5cm (5")	
Forceps Dissecting Adson Non Toothed 12.5cm (5")	
Forceps Dissecting Adson Toothed 12.5cm (5")	
Handle BP No 3 Graduated	
Needle Holder Lawrence 15cm (6")	
Retractor Skin Hook Gillies Fine 18cm (7.2")	
Scissors Dressing B/B Straight 12.5cm (5")	
Scissors Iris Straight 11cm (4.4")	



Scissors Strabismus Curved Blunt 11.5cm (4.6")	
Blade Carbon Steel No 15	
Towel Scrimex 48 x 40cm	2
Swab Gauze 10 x 10cm 12ply	5
Balls Non Woven Medium	3
Drape Small S/A Utility 76 x 76cm	2
Tray	
Sterile Field	

General Procedure Packs (Continued)

Suture Standard Pack

Ref No	F883008
UOM	30
Contains:	
Forceps Dissecting T.O.E. Serrated 12.5cm (5")	1
Forceps Dissecting Toothed Treves 15cm (6")	1
Needle Holder Kilner 13cm (5.2")	1
Scissors Iris Straight 11cm (4.4")	1
Towel Dressing 43 x 38cm	2
Swab Non Woven 7.5 x 7.5cm 4 ply	5
Tray	1
Sterile Field	1



Suture Standard Fine Pack

Ref No	F883009
UOM	30
Contains:	
Forceps Dissecting Iris Non Toothed 10cm (4")	1
Forceps Dissecting Iris Toothed 10cm (4")	1
Needle Holder Kilner 13cm (5.2")	1
Scissors Iris Straight 11cm (4.4")	1
Towel Dressing 43 x 38cm	2
Swab Non Woven 7.5 x 7.5cm	5
Tray	1
Sterile Field	1



General Procedure Packs (Continued)

Suture Basic Pack

Ref No	F883010
UOM	50
Contains:	
Forceps Dissecting T.O.E. Serrated 12.5cm (5")	1
Forceps Dissecting Toothed Treves 12.5cm (5")	1
Needle Holder Kilner 13cm (5.2")	1
Scissors Iris Straight 11cm (4.4")	1
Towel Dressing 43 x 38cm	1
Tray	1
Sterile Field	1



Suture Basic Fine Pack

Ref No	F883011
UOM	50
Contains:	
Forceps Iris Non Toothed 10cm (4")	1
Forceps Iris Toothed 10cm (4")	1
Needle Holder Kilner 13cm (5.2")	1
Scissors Iris Straight 11cm (4.4")	1
Towel Dressing 43 x 38cm	1
Tray	1
Sterile Field	1



General Procedure Packs

General Procedure Packs (Continued)

Suture Gold Standard Pack

Ref No	F883015
UOM	30
Contains:	
Forceps Artery Halstead Mosquito Curved 12.5cm (5")	1
Forceps Dissecting Iris Non Toothed 10cm (4")	1
Forceps Dissecting Iris Toothed 10cm (4")	1
Needle Holder Kilner 13cm (5.2")	1
Scissors Dressing S/S Straight 12.5cm (5")	1
Towel Dressing 43 x 38cm	1
Paper Blue Tissue/PE 90 x 90cm	1
Swab Non Woven 10 x 10cm 4Ply	5
Gloves Exam Nitrile, Powder Free, Cuffed, Walleted, M/L Pair	1



Bag Yellow Folded 45 x 29cm with Adhesive Strip	1
Tray	1
Sterile Field	1

Suture Gold Fine Pack

Ref No	F883016
UOM	30
Contains:	
Forceps Artery Halstead Mosquito Curved 12.5cm (5")	2
Forceps Dissecting Iris Non Toothed 10cm (4")	1
Forceps Dissecting Iris Toothed 10cm (4")	1
Needle Holder Kilner 13cm (5.2")	1
Scissors Iris Straight 11cm (4.4")	1
Towel Dressing 43 x 38cm	1
Swab Non Woven 10 x 10cm 4Ply	5
Gloves Exam Nitrile, Powder Free, Cuffed, Walleted, M/L, Pair	1
Bag Yellow Folded 45 x 29cm with Adhesive Strip	1



Paper Blue Tissue/PE 50 x 45cm	1
Tray	1
Sterile Field	1

General Procedure Packs (Continued)

Minor Op Pack

Ref No	F883014
UOM	25
Contains:	
Forceps Artery Halstead Mosquito Curved 12.5cm (5")	2
Forceps Dissecting Adson Non Toothed 12.5cm (5")	1
Handle BP No 3 Graduated	1
Needle Holder Kilner 13cm (5.2")	1
Scissors Iris Straight 11cm (4.4")	1
Blade Carbon Steel No 11	1
Swab Non Woven 5 x 5cm 4ply	5
Drape Fenestrated Deluxe N/A 70 x 90cm with 10cm Fen	1



Tray	1
Sterile Field	1

Minor Op Pack

Ref No	F883013
UOM	30
Contains:	
Forceps Dissecting Adson Toothed Ex Fine Micro 12.5cm (5")	1
Handle BP No 3 Graduated	1
Needle Holder Halsey Serrated 12.5cm (5")	1
Scissors Iris Straight 11cm (4.4")	1
Scissors Strabismus Curved Blunt 11.5cm (4.6")	1
Towel Dressing 43 x 38cm	1
Swab Non Woven 7.5 x 7.5cm	10



Tray	1
Sterile Field	1

General - Gynaecological

Type	Ref No	Description	Case Quantity
------	--------	-------------	---------------

Gynaecological

Cannulas

	Uterine		
	F893220	Currie	100
	F893221	Leech Wilkinson Small	8 X 20
	F893222	Leech Wilkinson Medium	8 X 20
	F893223	Leech Wilkinson Large	8 X 20
	F893224	Spackman 36cm (14.4")	30

Curette Sims


	F893227	Small Sharp 15cm (6")	8 X 10
	F893228	Medium Sharp 16cm (6.4")	8 X 10
	F893229	Large Sharp 19cm (7.4")	8 X 10
	F893326	Small Double Ended Sharp/Blunt 15cm (6")	8 X 10
	F893327	Medium Double Ended Sharp/Blunt 16cm (6.4")	8 X 10
	F893328	Large Double Ended Sharp/Blunt 19cm (7.4")	8 X 10

Catheters

	F893225	Silver Female 14fg	8 X 10
---	---------	--------------------	--------

Dilators

	Cervical		
	F893230	Hegar Double Ended 1/2mm 22cm (8.8")	8 X 20
	F893231	Hegar Double Ended 3/4mm 22cm (8.8")	8 X 20
	F893232	Hegar Double Ended 5/6mm 22cm (8.8")	8 X 20
	F893233	Hegar Double Ended 7/8mm 22cm (8.8")	8 X 20
	F893234	Hegar Double Ended 9/10mm 22cm (8.8")	8 X 20
	F893235	Hegar Double Ended 11/12mm 22cm (8.8")	8 X 20
	F893236	Hegar Double Ended 13/14mm 22cm (8.8")	8 X 10
	F893237	Hegar Double Ended 15/16mm 22cm (8.8")	8 X 10
	F893238	Hegar Double Ended 17/18mm 22cm (8.8")	8 X 10
	F883003	Hegar Double Ended Set of 9, 1/2...17/18mm 22cm (8.8")	20

Type	Ref No	Description	Case Quantity
Forceps			
	Biopsy		
	F893240	Eppendorfer Biopsy	5
	F893239	Tischler Morgan Biopsy	20
	Clamp		
	F893247	Bonney Polyp and Ovum 24cm (9.5")	8 X 10
	F893248	Kogan with Ratchet 24cm (9.5")	5
	F893094	Rampley Sponge Holding 18cm (7")	8 X 10
	F893095	Rampley Sponge Holding 24cm (9.5")	8 X 10
	Tissue		
	F893258	Bozemann 25cm (10")	50
	F893252	Green Armytage 20cm (8")	8 X 10
	F893242	Vulsellum Lower Toothed 1:1 24cm (9.5")	100
	F893255	Vulsellum Teale Uterine Curved Toothed 3:4 18cm (7")	8 X 10
	F893256	Vulsellum Teale Uterine Curved Toothed 3:4 23cm (9")	8 X 10
	F893257	Vulsellum Teale Uterine Straight Toothed 3:4 23cm (9")	8 X 10
	F893314	Maingot Hysterectomy 20cm (8")	50
	Dissecting		
	F893241	Bonney Non Toothed 18cm (7.2")	8 X 10
	F893254	Bonney Toothed 18cm (7")	8 X 10
Hooks			
	F893251	Saunders IUCD Removal 32cm (12.5")	8 X 25
	F893278	Emmett Cervical Stay Hook	8 X 20
Retractors			
	Hand Held		
	F893154	Doyen Curved 25cm (12") Blade 8cm (3.5")	20
	F893262	Doyen Curved 25cm (12") Blade 10cm (4")	20
	Self Retaining		
	F893311	Gelpi Vulva 16.5cm (6.5")	50

General - Gynaecological

Type	Ref No	Description	Case Quantity
Scissors			
	Dissecting		
	F893190	Mayo C.O.F 16.5cm (6.5")	8 X 10
	F893191	Mayo C.O.F 19cm (7.5")	8 X 10
	F893192	Mayo C.O.F 22cm (8.5")	8 X 10
	F893194	Mayo Straight 14cm (5.5")	8 X 20
	F893196	Mayo Straight 16.5cm (6.5")	8 X 10
	F893213	Mayo Straight 19cm (7.5")	8 X 10
	F893197	Mayo Straight 22cm (8.5")	8 X 10
	F893266	Sims Uterine Straight 20cm (8")	8 X 10
	F893267	Sims Uterine Curved 20cm (8")	8 X 10
	Episiotomy		
	F893263	Bergman 23cm (9")	50
	F893264	Braun Stadler 14cm (5.5")	12 X 10
	F893265	Mayo 24cm (9.5")	100
	F893268	Vant 24cm (9.5")	8 X 10
Screws			
	F893315	Myomectomy Screw 19cm (7")	8 X 10
Sounds			
	Umbilical		
	F893018	Umbilical 12.5cm (5")	12 X 20
	Sound		
	F893269	Galabin Malleable 29cm (11.5")	100
	F893270	Horrocks Long 29cm (11.5")	100
	F893271	Simpson 29cm (11.5")	100
Specula			
	F810004	Cusco Extra Small	40
	F810005	Cusco Small	40
	F810006	Cusco Medium	40
	F810007	Cusco Medium Long	40

Type	Ref No	Description	Case Quantity
Specula (Continued)			
	F810008	Cusco Large	40
	F810009	Winterton	20
	F810002	Cusco Medium Laser with Smoke Tube Blackened	20
	F810003	Cusco Large Laser with Smoke Tube Blackened	20
	F893259	Sims Small Single Ended	20
	F893272	Sims Small Double Ended	2 X 20
	F893260	Sims Medium Single Ended	20
	F893273	Sims Medium Double Ended	20
	F893261	Sims Large Single Ended	20
	F893274	Sims Large Double Ended	20
	F893316	Auward	10
Suckers			
	F893218	Poole 24cm (9.5")	100
	F893217	Simpson Smith 31fg	8 X 10
Trocars			
	F893219/HRT	HRT Trocar & Cannula Hormone Implant Set	8 X 10

Gynaecological Instrument Sets



Gynaecological Instrument Sets

Perineal Set

Ref No	F883023	Forceps Dissecting T.O.E. Serrated 12.5cm (5")	1
UOM	10	Scissors Mayo Straight 17.5cm (7")	1
Contains:		Scissors Mayo Straight 19cm (7.6")	1
Forceps Sponge Holding Rampley 24cm (9.5")	3	Needle Holder Mayo Hegar 18cm (7")	1
Forceps Artery Dunhill Curved 12.8cm (5")	12	Handle BP No 3	2
Forceps Artery Spencer Wells Curved 18cm (7.2")	12	Clip Towel Cross Action	4
Forceps Dissecting Toothed Treves 12.5cm (5")	1	Tray	1
		Sterile Field	1

Laparoscopy Set

Ref No	F883035	Forceps Dissecting Toothed Treves 12.5cm (5")	1
UOM	10	Handle BP No 3	1
Contains:		Scissors Mayo Straight 17.5cm (7")	1
Forceps Sponge Holding Rampley 24cm (9.5")	2	Needle Verres 120mm	1
Clip Towel Backhaus 11cm (4.5")	4	Tray	1
Forceps Dissecting T.O.E. Serrated 12.5cm (5")	1	Sterile Field	1

Hysterosalpingogram & Hydrotubation Set

Ref No	F883036	Speculum Sims D/E Medium	1
UOM	25	Forceps Artery Roberts Curved 22cm (8.8")	3
Contains:		Tray	1
Speculum Sims D/E Small	1	Sterile Field	1

Gynaecological Instrument Sets (Continued)

Abdominal Hysterectomy Set

Ref No	F883030		
UOM	10		
Contains:			
Forceps Sponge Holding Rampley 24cm (9.5")	5	Forceps Tissue Littlewood 2:3 Teeth 19cm (7.5")	2
Clip Towel Backhaus 11cm (4.5")	4	Handle BP No 4	2
Pin Instrument Mayo 15cm (6")	4	Needle Holder Mayo Hegar 18cm (7")	1
Forceps Dissecting Bonney Non Toothed 18cm (7")	1	Needle Holder Mayo 20cm (8")	2
Forceps Dissecting Bonney Toothed 18cm (7")	1	Scissors Dressing S/S Straight 12.5cm (5")	1
Forceps Artery Spencer Wells Curved 12.5cm (5")	10	Scissors Mayo Straight 19cm (7.6")	1
Forceps Artery Spencer Wells Straight 20cm (8")	5	Scissors Mayo C.O.F. 19cm (7.6")	1
Forceps Artery Kocher Toothed Straight 14cm (5.5")	6	Forceps Vulsellum Teale Uterine Curved 3:4 Teeth 23cm (9")	2
Forceps Artery Kocher Toothed Curved 20cm (8")	6	Forceps Tissue Maingot Hysterectomy 20cm (8")	2
Forceps Moynihan Cholecystectomy 22.5cm (9")	4	Myomectomy Screw 19cm (7")	1
Forceps Tissue Allis 2:3 Teeth 15cm (6")	4	Forceps Dressing Lister Sinus 12.5cm (5")	1
Forceps Tissue Judd-Allis 3:4 Teeth 19cm (7.6")	2	Suction Tube Poole Straight 24cm (9.5")	1
		Retractor Gelpi Vulva 16.5cm (6.5")	1
		Forceps Ligature Lahey 20cm (8")	2
		Diathermy Bipolar Forceps Imm Tip & 3m Lead 15cm (6")	1
		Tray	1
		Sterile Field	1

Gynaecological Instrument Sets



Gynaecological Instrument Sets (Continued)

Gynaecological Basic Set

Ref No	F883028	Forceps Tissue Allis 2:3 Teeth 15cm (6")	3
UOM	10	Forceps Tissue Littlewood 2:3 Teeth 19cm (7.5")	3
Contains:		Handle BP No 4	2
Forceps Sponge Holding Rampley 24cm (9.5")	6	Needle Holder Mayo Hegar 18cm (7")	2
Clip Towel Backhaus 11cm (4.5")	6	Scissors Dressing B/B Straight 12.5cm (5")	1
Pin Instrument Mayo 15cm (6")	4	Scissors Mayo Curved 16.5cm (6.5")	1
Forceps Dissecting Bonney Non Toothed 18cm (7")	1	Scissors Mayo Straight 19cm (7.6")	1
Forceps Dissecting Bonney Toothed 18cm (7")	2	Scissors Mayo C.O.F. 19cm (7.6")	1
Forceps Dissecting T.O.E. Serrated 12.5cm (5")	3	Forceps Lower Vulsellum 1:1 Teeth 25cm (10")	1
Forceps Dissecting Lane Toothed 1:2 Teeth 12.5cm (5")	1	Suction Tube Poole Straight 24cm (9.5")	1
Forceps Artery Dunhill Curved 12.8cm (5")	12	Retractor Morris 24.5cm (9.75") Blade Small 70 x 50mm	2
Forceps Artery Spencer Wells Curved 18cm (7.2")	12	Retractor Doyen Gynae Curved 10cm (4")	1
Forceps Artery Kocher Toothed Straight 14cm (5.5")	10	Tray	1
		Sterile Field	1

Vaginal Suture Set

Ref No	F883029	Handle BP No 3	1
UOM	10	Handle BP No 4	1
Contains:		Needle Holder Mayo Hegar 18cm (7")	1
Forceps Sponge Holding Rampley 24cm (9.5")	2	Scissors Dressing S/S Straight 12.5cm (5")	1
Clip Towel Cross Action	4	Scissors Mayo Straight 17.5cm (7")	1
Forceps Dissecting T.O.E. Serrated 12.5cm (5")	1	Scissors Mayo Curved 16.5cm (6.5")	1
Forceps Dissecting Toothed Treves 12.5cm (5")	1	Forceps Vulsellum Teale Uterine Curved 3:4 Teeth 23cm (9")	1
Forceps Artery Curved Kocher 1:2 Teeth 21cm (8")	2	Speculum Sims D/E Medium	1
Forceps Artery Crile Straight 14cm (5.5")	2	Auvard Speculum	1
Forceps Tissue Allis 2:3 Teeth 15cm (6")	2	Tray	1
		Sterile Field	1

Gynaecological Instrument Sets (Continued)

Vaginal Repair & Hysterectomy Set

Ref No	F883032		
UOM	10		
Contains:			
Forceps Dissecting Toothed Treves 12.5cm (5")	1	Forceps Artery Spencer Wells Straight 20cm (8")	6
Forceps Dissecting T.O.E. Serrated 12.5cm (5")	1	Forceps Moynihan Cholecystectomy 22.5cm (9")	4
Forceps Dissecting Bonney Non Toothed 18cm (7")	1	Forceps Tissue Lane 17.5cm (7")	4
Forceps Dissecting Bonney Toothed 18cm (7")	1	Handle BP No 4	2
Forceps Artery Kocher Toothed Straight 14cm (5.5")	20	Needle Holder Thomson Walker 20cm (8")	2
Forceps Artery Cairn 14cm (5.5")	10	Clip Towel Backhaus 11cm (4.5")	4
Forceps Artery Curved Kocher 1:2 Teeth 21cm (8")	6	Scissors Dressing B/B Straight 12.5cm (5")	1
		Scissors Mayo Straight 14cm (5.5")	1
		Scissors Mayo Curved 16.5cm (6.5")	1
		Retractor Doyen Gynae Curved 10cm (4")	2
		Tray	1
		Sterile Field	1

Tubal Ligation Set

Ref No	F883034		
UOM	10		
Contains:			
Forceps Sponge Holding Rampley 24cm (9.5")	5	Needle Holder Mayo Hegar 18cm (7")	2
Clip Towel Cross Action	5	Scissors Mayo Straight 14cm (5.5")	1
Forceps Dissecting T.O.E. Serrated 12.5cm (5")	2	Scissors Mayo Curved 16.5cm (6.5")	1
Forceps Dissecting Lane Toothed 1:2 Teeth 18cm (7")	2	Scissors Dressing S/S Straight 12.5cm (5")	1
Forceps Artery Crile Straight 14cm (5.5")	5	Needle Aneurysm Small 15cm (6")	2
Forceps Artery Spencer Wells Curved 18cm (7.2")	10	Retractor Langenbeck Small 21.5cm (8.5") Blade 25 x 6mm	2
Forceps Dissecting Lane Toothed 1:2 Teeth 12.5cm (5")	2	Retractor Langenbeck Medium 21.5cm (8.5") Blade 35 x 11mm	2
Forcep Clip Michel Applier 12cm (4.7")	1	Retractor Morris 24.5cm (9.75") Blade Small 70 x 50mm	2
Handle BP No 4	1	Retractor Balfour Abdominal	1
		Tray	1
		Sterile Field	1

Gynaecological Instrument Sets



Gynaecological Instrument Sets (Continued)

D & C Set

Ref No	F883038	Sound Uterine Galabin Malleable 29cm (11.5")	1
UOM	10	Speculum Sims Medium D/E	1
Contains:		Curette Sims Large D/E 19cm (7.6")	1
Forceps Sponge Holding Rampley 24cm (9.5")	3	Curette Sims Small D/E 15cm (6")	1
Clip Towel Cross Action	4	Auvard Speculum	1
Forceps Dissecting Bonney Non Toothed 18cm (7")	1	Catheter Silver Female	1
Forceps Dissecting Bonney Toothed 18cm (7")	1	Forceps Polyp and Ovum Bonneys 24cm (9.6")	1
Forceps Tissue Vulsellum Teale Curved Toothed 3:4 23cm (9")	2	Tray	1
Dilators Hegar D/E Set of 9 1/2...17/18 22cm (8.8")	1	Sterile Field	1

Prolapse Set

Ref No	F883033	Forceps Tissue Littlewood 2:3 Teeth 19cm (7.5")	10
UOM	10	Handle BP No 4	2
Contains:		Needle Holder Mayo Hegar 18cm (7")	2
Clip Towel Backhaus 11cm (4.5")	8	Scissors Mayo Straight 16.5cm (6.5")	2
Forceps Dissecting Bonney Non Toothed 18cm (7")	1	Scissors Mayo Curved 16.5cm (6.5")	2
Forceps Dissecting Bonney Toothed 18cm (7")	1	Forceps Lower Vulsellum 1:1 Teeth 25cm (10")	1
Forceps Artery Crile Straight 14cm (5.5")	10	Tray	1
Forceps Artery Spencer Wells Curved 18cm (7.2")	5	Sterile Field	1

Gynaecological Instrument Sets (Continued)

Prostaglandin Set

Ref No	F883031	Speculum Cusco Medium	1
UOM	10	Sound Uterine Galabin Malleable 29cm (11.5")	1
Contains:		Speculum Sims D/E Medium	1
Forceps Sponge Holding Rampley 24cm (9.5")	1	Tray	1
Forceps Lower Vulsellum 1:1 Teeth 25cm (10")	2	Sterile Field	1

Gynaecological Procedure Packs

Gynaecological Procedure Packs

IUCD Pack

Ref No	F883001
UOM	14
Contains:	
Forceps Sponge Holding Rampley 24cm (9.5")	1
Forceps Vulsellum Teale Uterine Curved 3:4 Teeth 23cm (9")	1
Scissors Mayo C.O.F. 19cm (7.6")	1
Speculum Cusco Medium	1
Balls Non Woven Medium	5
Towel Dressing 43 x 38cm	1



Paper Blue Tissue/PE 50 x 50cm	1
Tray	1
Sterile Field	1

IUCD Pack (With Sound)

Ref No	F883002
UOM	8
Contains:	
Forceps Sponge Holding Rampley 24cm (9.5")	1
Forceps Vulsellum Teale Uterine Curved 3:4 Teeth 23cm (9")	1
Scissors Mayo C.O.F. 19cm (7.6")	1
Sound Uterine Galabin Malleable 29cm (11.5")	1
Speculum Cusco Medium	1
Gallipot Polypropylene 120ml	1
Balls Non Woven Medium	5
Dish Kidney Polypropylene 10"	1
Paper Blue Tissue/PE 50 x 50cm	1



Towel Dressing 43 x 38cm	1
Tray	1
Sterile Field	1

Gynaecological Procedure Packs (Continued)

Essential Instrument Delivery Pack

Ref No	F810131
UOM	6
Contains:	
Foam Swab White 38 x 38 x 25mm	4
Bowl Polypropylene 500ml	1
Gallipot Polypropylene 60ml	1
Clamp Umbilical Cord	1
Towel Dressing 43 x 38cm	4
Baby Wrap White 100 x 100cm	1
Dressing Pad 10 x 20cm	5
Forceps Artery Spencer Wells Straight 18cm (7.2")	2



Scissors Mayo Straight 23cm (9")	1
Scissors U.S.A. Pattern Umbilical Straight 10cm (4")	1
Tray	1
Sterile Field	1

Essential Episiotomy Pack

Ref No	F810130
UOM	18
Contains:	
Towel Dressing 43 x 38cm	1
Gallipot Polypropylene 60ml	2
Swab XRD Tampon and Tape	1
Drape Fenestrated 122 x 152cm With Central Adhesive Fenestration 12 x 15cm	1
Needle Holder Mayo 20cm (8")	1
Forceps Dissecting Bonney Non Toothed 18cm (7")	1
Forceps Artery Spencer Wells Straight 12.5cm (5")	1



Scissors Mayo Straight 16.5cm (6.5")	1
Tray	1
Sterile Field	1

Gynaecological Procedure Packs

Gynaecological Procedure Packs (Continued)

IUCD Instrument Pack

Ref No	F883018
UOM	14
Contains:	
Forceps Sponge Holding Rampley 24cm (9.5")	1
Forceps Vulsellum Teale Uterine Curved 3:4 Teeth 23cm (9")	1
Scissors Mayo C.O.F. 19cm (7.6")	1
Cotton Wool Balls Large	5
Dish Kidney Polypropylene 10"	1
Towel Dressing 43 x 38cm	1
Gallipot Polypropylene 120ml	1



Tray	1
Sterile Field	1

Delivery Pack

Ref No	F810122
UOM	5
Contains:	
Bowl Placenta 3000cc	1
Pulp Dish Kidney White	2
Swab Non Woven 10 x 10cm 4Ply	5
Clamp Umbilical Cord	1
Dressing Pad 10 x 10cm	2
Baby Wrap White 100 x 100cm	1
Towel Scrimex 48 x 40cm	2
Sterile Field	1



Gynaecological Procedure Packs (Continued)

Abdominal Pack

Ref No F300945

UOM 5

Contains:

Swab Gauze XRD 30 x 30cm 12ply Taped	10
Swab Gauze XRD 10 x 7.5cm 32ply	10
Needle Counter 15 Count Magnet & Foam Strips	1
Blade Carbon Steel No 15	2
Blade Carbon Steel No 11	1
Needle Hypodermic 19g x 1.5" White	1
Needle Hypodermic 21g x 1.5" Green	1
Syringe 20ml Luer Slip Eccentric	1
Tip Cleaner	1
Skin Closure Steri-Strip 12mm x 10cm (1 Sheet of 3)	2
Dressing Primapore 8.3 x 6cm	2
Towel Scrimex 48 x 40cm	4
Cover Light Handle	2
Blade Carbon Steel No 10	2
Gown Large No Reinforcement	2
Drape Top S/A 150 x 240cm Reinforced with Tube Holders	1
Drape Bottom S/A 175 x 180cm Reinforced with Tube Holders	1
Drape Side S/A 90 x 90cm Reinforced	2
Chloraprep 10.5ml	2
Skin Marker Pen Regular & Ruler	1
Scissors Metzenbaum C.O.F 19cm (7.5")	1
Forcep Sponge Holding 19cm (7.5")	2




Diathermy Quiver	1
Forceps Artery Kocher Toothed Straight 14cm (5.5")	2
Retractor Kelly 21.5cm (1.5") Blade 38mm	1
Handle BP No 3	1
Monopolar Diathermy Forcep 17.8cm (7") 8mm (Bovie) With Lead 3m	1
Needle Holder Mayo 20cm (8")	1
Clip Towel Ball & Socket	2
Forceps Artery Halstead Mosquito Curved 12.5cm (5")	2
Forceps Adson Dissecting Fine Toothed 12.5cm (5")	2
Handle Scalpel Barron	2
Needle Holder Crile Wood 15cm (6")	2
Retractor Skin Hook Gillies Fine 18cm (7.2")	2
Retractor Langenbeck Medium 21.5cm (8.5") Blade 35 x 11mm	2
Scissors Dressing S/B Straight 12.5cm (5")	1
Tray	1
Sterile Field	2

General - Genito Urinary

Type	Ref No	Description	Case Quantity
------	--------	-------------	---------------

Genito Urinary

Forceps

	Clamp		
	F893001	Denis Browne (Bead) 14cm (5.5")	12 X 10
	F893276	Thompson Walker Penile 12cm (4.75")	2 X 25
	F893109	Vasectomy Modified Straight Blunt 12.5cm (5")	12 X 20
	F893106	Vasectomy Straight 15cm (6")	8 X 20
	Tissue		
	F893051	Soonawalla 15cm (6")	8 X 10
	Stone		
	F893318	Desjardins Grasping Small 23cm (9")	25
	F893317	Desjardins Grasping Large 23cm (9")	25

Introducers

	Catheter		
	F893275	De Pezzer Curved	50

Retractors

	Self Retaining		
	F893319	Millin 80 x 25mm	5

Genito Urinary Instrument Sets

Vasectomy Set

Ref No	F883027	Forceps Artery Halstead Mosquito Curved 12.5cm (5")	5
UOM	10	Needle Holder Kilner 13cm (5.2")	1
Contains:		Scissors Mayo Straight 14cm (5.5")	1
Forceps Dissecting T.O.E. Serrated 12.5cm (5")	1	Forceps Sponge Holding Rampley 24cm (9.5")	1
Forceps Dissecting Toothed Treves 12.5cm (5")	1	Tray	1
Forceps Vasectomy Straight 15.5cm (6")	2	Sterile Field	1
Forceps Artery Halstead Mosquito Straight 12.5cm (5")	5		

Genito Urinary Procedure Packs

Genito Urinary Procedure Packs

Vasectomy Pack

Ref No	F883012
UOM	8
Contains:	
Forceps Artery Halstead Mosquito Curved 12.5cm (5")	2
Forceps Artery Halstead Mosquito Straight 12.5cm (5")	2
Forceps Dissecting Adson Toothed 12.5cm (5")	1
Forceps Dissecting Gillies Toothed 15cm (6")	1
Forceps Sponge Holding Rampley 24cm (9.5")	1
Forceps Tissue Allis 3:4 15cm (6")	2
Forceps Vasectomy Straight 15.5cm (6")	2
Handle BP No 3 Graduated	1
Needle Holder Crile Wood 15cm (6")	1
Scissors Dressing S/B Straight 12.5cm (5")	1
Scissors Iris Straight 11cm (4.4")	1



Scissors Strabismus Straight Blunt 11.5cm (5.6")	1
Gallipot Polypropylene 60ml	1
Dish Kidney Polypropylene 10"	1
Drape Fenestrated Deluxe N/A 70 x 90cm (10cm Fen)	1
Towel Scrimex 48 x 40cm	2
Gown Large No Reinforcement	1
Tray	1
Sterile Field	1

Type	Ref No	Description	Case Quantity
------	--------	-------------	---------------

Laparoscopy

Instruments



Dissector

F898023	Maryland Dissector 5 x 440mm (Long)	12
F898002	Maryland Dissector 5 x 330mm (Standard)	12
F898024	Maryland Dissector 5 x 170mm (Short)	12



Grasper With Ratchet

F898025	Atraumatic Fenestrated Duckbill With Ratchet 5 x 440mm (Long)	12
F898003	Atraumatic Fenestrated Duckbill With Ratchet 5 x 330mm (Standard)	12
F898026	Atraumatic Fenestrated Duckbill With Ratchet 5 x 170mm (Short)	12
F898027	Babcock With Ratchet 5 x 440mm (Long)	12
F898018	Babcock With Ratchet 5 x 330mm (Standard)	12
F898028	Babcock With Ratchet 5 x 170mm (Short)	12
F898029	Clinch Fenestrated With Ratchet 5 x 440mm (Long)	12
F898017	Clinch Fenestrated With Ratchet 5 x 330mm (Standard)	12
F898030	Clinch Fenestrated With Ratchet 5 x 170mm (Short)	12
F898031	Atraumatic Tapered With Ratchet 5 x 440mm (Long)	12
F898019	Atraumatic Tapered With Ratchet 5 x 330mm (Standard)	12
F898032	Atraumatic Tapered With Ratchet 5 x 170mm (Short)	12



Grasper Without Ratchet

F898033	Atraumatic Fenestrated Duckbill 5 x 440mm (Long)	12
F898020	Atraumatic Fenestrated Duckbill 5 x 330mm (Standard)	12
F898034	Atraumatic Fenestrated Duckbill 5 x 170mm (Short)	12
F898035	Clinch Fenestrated 5 x 440mm (Long)	12
F898021	Clinch Fenestrated 5 x 330mm (Standard)	12
F898036	Clinch Fenestrated 5 x 170mm (Short)	12
F898037	Atraumatic Tapered 5 x 440mm (Long)	12
F898022	Atraumatic Tapered 5 x 330mm (Standard)	12
F898038	Atraumatic Tapered 5 x 170mm (Short)	12



Hook

F898011	L Hook Blunt High Frequency Electrode 5 x 330mm	12
F898012	L Hook Cutting High Frequency Electrode 5 x 330mm	12
F898039	J Hook Blunt High Frequency Electrode 5 x 330mm	12

General - Laparoscopy

Type	Ref No	Description	Case Quantity
Reposables			
	Reposable		
	F898100	Reposable Handle Without Ratchet	5
	F898101	Metzenbaum Curved High Insulation 5 x 330mm Shaft for Repposable Handle	12
	Scissors		
	F898041	Metzenbaum Curved High Insulation 5 x 440mm (Long)	12
	F898016	Metzenbaum Curved High Insulation 5 x 330mm (Standard)	12
	F898042	Metzenbaum Curved High Insulation 5 x 170mm (Short)	12
	F898101	Metzenbaum Curved High Insulation 5 x 330mm Shaft for Repposable Handle	12
Sundries			
	Cable		
	F898043	Monopolar Cable 8mm x 3m	12
	Needle		
	F898004	Verres Needle 120mm	12
	F898040	Verres Needle 150mm	12
	Retrieval		
	F898044	Endo Retrieval Bag 250ml 10mm	12
	F898009	Endo Retrieval Bag 400ml 10mm	12
	F898045	Endo Retrieval Bag 800ml 10mm	12
	F898048	Endo Retrieval Pouch 75 x 150mm	12
	Solution		
	F898046	Anti-fog Solution	20
	Tubing		
	F898005	Suction/Irrigation With Trumpet Valve 10mm	12
	F898008	Suction/Irrigation With Trumpet Valve 5mm	12
	F898047	Tubing Insufflation & Filter Luer Lock Ends 3m	12

Laparoscopy Procedure Packs

Laparoscopic Cholecystectomy Pack

Ref No	F300952
UOM	2
Contains:	
Swab Gauze XRD 10 x 7.5cm 32ply	10
Needle Blunt Fill Red Safety 18g x 1.5"	1
Needle Hypodermic 23g x 1.25" Blue	1
Syringe 10ml Luer Lock	1
Syringe 20ml Luer Lock	1
Dressing Mepore 6 x 7cm	4
Cover Light Handle Green	2
Towel Scrimex 48 x 40cm	2
Swab Gauze XRD 15 x 2.5cm 16ply	5
Scalpel BP Handle With No 15 Blade	1
Tip Cleaner	1
Sharps Pad	1
Tube Laparoscopic Trumpet Suction Irrigation Set	1
Anti-fog Solution	1
Tubing Insufflation & Filter Luer Lock Ends 3m	1
Drape Camera Video 15 x 220cm With Perforated Tip	1
Drape Instrument Pouch With Three Compartments 30 x 50cm	1
Drape Side S/A 90 x 90cm Reinforced	2
Drape Bottom S/A 175 x 180cm Reinforced With Tube Holders	1



Drape Top S/A 150 x 240cm Reinforced With Tube Holders	1
Drape Trolley Cover 140 x 190cm Fan Folded	1
Drape Mayo Cover 75 x 144cm Reinforced	1
Suture Polysorb CL-673G	2
Suture Polysorb SL-5678	1
Laparoscopic Dissector Maryland 5 x 330mm	1
Laparoscopic Grasper Clinch Fenestrated 5 x 330mm	1
Laparoscopic Grasper Clinch Fenestrated 5 x 330mm With Ratchet	1
Laparoscopic Scissor Curved Metzenbaum 5 x 330mm	1
Pouch Tissue Retrieval 75 x 150mm	1
Tray	2
Sterile Field	1

Laparoscopy Procedure Packs

Laparoscopy Procedure Packs (Continued)

Laparoscopic General Appendix Pack

Ref No	F300953
UOM	2
Contains:	
Swab Gauze XRD 15 x 2.5cm 16ply	5
Swab Gauze XRD 10 x 7.5cm 32ply	5
Swab Gauze XRD 22.5 x 22.5cm 12ply Untaped	5
Sharps Box 30 Count With Foam Strips	1
Cover Light Handle Green	2
Syringe 20ml Luer Lock	1
Dressing Mepore 6 x 7cm	3
Needle Blunt Fill Red Safety 18g x 1.5"	1
Needle Hypodermic 21g x 1.5" Green	1
Scalpel BP Handle With No 15 Blade	1
Dish Kidney Polypropylene 10"	2
Towel Scrimex 48 x 40cm	2
Bowl Polypropylene 250ml	2
Tip Cleaner	1
Tubing Insufflation & Filter Luer Lock Ends 3m	1
Anti-fog Solution	1
Drape Side S/A 90 x 90cm Reinforced	2
Drape Bottom S/A 175 x 180cm Reinforced With Tube Holders	1
Drape Top S/A 150 x 240cm Reinforced With Tube Holders	1
Drape Trolley Cover 140 x 190cm Fan Folded	1



Drape Camera Video 15 x 220cm With Perforated Tip	1
Drape Instrument Pouch With Three Compartments 30 x 50cm	1
Suture Polysorb CL-673G	2
Suture Polysorb SL-5678	1
Laparoscopic Dissector Maryland 5 x 330mm	1
Laparoscopic Grasper Fenestrated Duckbill Atraumatic 5 x 330mm	1
Laparoscopic Monopolar Electrode L Hook Blunt 5 x 330mm	1
Laparoscopic Grasper Tapered Atraumatic 5 x 330mm	1
Laparoscopic Grasper Babcock 5 x 330mm With Ratchet	1
Laparoscopic Scissor Curved Metzenbaum 5 x 330mm	1
Tray	2
Sterile Field	1

Laparoscopy Procedure Packs (Continued)

Laparoscopic General Colorectal Pack

Ref No F300954

UOM 2

Contains:

Swab Gauze XRD 10 x 7.5cm 32ply	10
Swab Gauze XRD 30 x 30cm 12ply Taped	5
Swab Gauze XRD 15 x 2.5cm 16ply	5
Scalpel BP Handle With No 10 Blade	1
Scalpel BP Handle With No 15 Blade	1
Needle Blunt Fill Red Safety 18g x 1.5"	1
Needle Hypodermic 23g x 1.25" Blue	1
Syringe 20ml Luer Lock	1
Syringe 10ml Luer Lock	1
Sharps Pad	1
Cover Light Handle Green	2
Dressing Mepore 9 x 15cm	1
Dressing Mepore 6 x 7cm	3
Towel Scrimex 48 x 40cm	2
Needle Hypodermic 21g x 1.5" Green	1
Tip Cleaner	1
Suction Tube Yankauer Midi	1
Tubing Insufflation & Filter Luer Lock Ends 3m	1
Anti-fog Solution	1
Tube Laparoscopic Trumpet Suction Irrigation Set	1
Drape Camera Video 15 x 220cm With Perforated Tip	1
Drape Instrument Pouch With Three Compartments 30 x 50cm	2



Drape Side S/A 90 x 90cm Reinforced	3
Drape Bottom S/A 175 x 180cm Reinforced With Tube Holders	2
Drape Lithotomy Leggings 78 x 121cm (Each)	2
Drape Mayo Cover 75 x 144cm Reinforced	1
Drape Underbuttock 95 x 135cm	1
Suture Polysorb CL-673G	1
Suture Polysorb SL-5678	2
Connecting Tube 7mm x 3m FFM	1
Laparoscopic Dissector Maryland 5 x 330mm	1
Laparoscopic Grasper Fenestrated Duckbill Atraumatic 5 x 330mm	1
Laparoscopic Scissor Curved Metzenbaum 5 x 330mm	1
Laparoscopic Monopolar Electrode L Hook Blunt 5 x 330mm	1
Tray	2
Sterile Field	1

Laparoscopy Procedure Packs

Laparoscopy Procedure Packs (Continued)

Laparoscopic General Pack

Ref No	F300955
UOM	5
Contains:	
Scissors Dressing B/B Straight 12.5cm (5")	1
Retractor Langenbeck Large 22cm (8.5") Blade 4l x 11mm	2
Needle Holder Mayo Long 23cm (9")	2
Forceps Sponge Holding Rampley 24cm (9.5")	2
Forceps Artery Spencer Wells Straight 23cm (9.2")	2
Forceps Dissecting Gillies Toothed 15cm (6")	1
Forceps Tissue Littlewoods 2:3 Teeth 19cm (7.5")	2
Scissors Mayo C.O.F 19cm (7.6")	1
Scissors McIndoe Curved 18cm (7")	1
Tubing Insufflation & Filter Luer Lock Ends 3m	1
Anti-fog Solution	1
Drape Camera Video 15 x 220cm With Perforated Tip	1



Laparoscopic Scissors Curved Metzendaum 5 x 330mm	1
Laparoscopic Dissector Maryland 5 x 330mm	1
Laparoscopic Grasper Tapered Atraumatic 5 x 330mm	2
Tray	1
Sterile Field	1

Laparoscopy Procedure Packs (Continued)

Laparoscopic Gynae Diagnostic Pack

Ref No	F300956
UOM	4
Contains:	
Swab Gauze XRD 10 x 7.5cm 32ply	5
Cover Light Handle Green	1
Sharps Pad	1
Cannula Spackman Intra Uterine 36cm (14.4")	1
Blade Carbon Steel No 11	1
Needle Hypodermic 21g x 1.5" Green	1
Syringe 20ml Luer Lock	1
Dressing Mepore 6 x 7cm	1
Anti-fog Solution	1
Tubing Insufflation & Filter Luer Lock Ends 3m	1
Needle Verres 120mm	1
Drape Camera Video 15 x 220cm With Perforated Tip	1
Drape Tape Strip 10 x 50cm	1
Drape Instrument Pouch With Three Compartments 30 x 50cm	1



Drape Top S/A 150 x 240cm Reinforced With Tube Holders	1
Drape Side S/A 90 x 90cm Reinforced	3
Drape Lithotomy Leggings 78 x 121cm (Each)	2
Drape Underbuttock 95 x 135cm	1
Tray	1
Sterile Field	1

Type	Ref No	Description	Case Quantity
Obstetrics			
Dilators			
	Umbilical		
	F893018	Umbilical 12.5cm (5")	12 X 20
Forceps			
	Delivery		
	F893243	Anderson	5
	F893244	Keilland	5
	F893245	Simpson	5
	F893246	Wrigley	5
Scissors			
	Umbilical		
	F893212	U.S.A Pattern Straight 10cm (4")	12 X 20
Suckers			
	F893218	Poole 24cm (9.5")	100



Obstetrics Instrument Sets

Caesarian Set

Ref No	F883037
UOM	10
Contains:	
Forceps Sponge Holding Rampley 24cm (9.5")	5
Forceps Dissecting Bonney Non Toothed 18cm (7")	2
Forceps Dissecting Bonney Toothed 18cm (7")	2
Forceps Artery Crile Curved 14cm (5.5")	10
Forceps Artery Spencer Wells Curved 18cm (7.2")	10
Forceps Tissue Allis 2:3 Teeth 15cm (6")	2
Forceps Tissue Littlewood 2:3 Teeth 19cm (7.5")	3
Handle BP No 4	2
Needle Holder Mayo Hegar 18cm (7")	2

Scissors Dressing S/S Straight 12.5cm (5")	1
Scissors Mayo Curved 18cm (7.6")	1
Scissors Mayo Straight 19cm (7.6")	1
Forceps Obstetric Wrigley	1
Forceps Uterine Green Armytage 20cm (8")	6
Retractor Doyen Gynae Curved 10cm (4")	2
Retractor Morris 24.5cm (9.75") Blade Small 70 x 50mm	1
Retractor Langenbeck Medium 21.5cm (8.5") Blade 35 x 11mm	1
Forcep Clip Michel Applier 12cm (4.7")	1
Clip Towel Cross Action	5
Tray	1
Sterile Field	1

Obstetrics Procedure Packs

Caesarian Pack

Ref No	F300178
UOM	3
Contains:	
Bag Yellow Poly 50 x 38cm Adhesive Strip on reverse	1
Balls Non Woven Medium	10
Gallipot Polypropylene 120ml	1
Maternity Pad Unlooped Size 2	1
Swab Gauze XRD 45 x 45cm 12ply	5
Swab Gauze XRD 10 x 7.5cm 32ply	5
Towel Scrimex 48 x 40cm	1
Baby Blanket 80 x 80cm	1
Catheter Foley Aquamatic 14fg (10ml Balloon)	1
Drape Mayo Cover 58 x 137cm	1
Suction Tube Yankauer Wide Bore Maxi Without VC	1
Blade Carbon Steel No 22	1
Blade Scalpel Remover	1



Dressing Mepore 9 x 25cm	1
Drape Utility N/A 75 x 90cm	2
Drape C Section 254 x 305cm Full Incise, Fluid Collection Pouch	1
Connecting Tube 7mm x 3m FFM	1
Sharps Pad	1
Handle BP No 4	1
Bag Urine Drainage 2L and Tap	1
Paper Crepe White 137 x 137cm	1
Tray	1
Sterile Field	1

Introducing our Neurosurgery Range incorporating:-

- ENT
- Maxillo Facial
- Dental



Fine Surgical Instruments

Type	Ref No	Description	Case Quantity
Braces			
	F894067	Hudson with 14mm Perforator and 16mm Burr	5
Clips			
	F895006	Serrafine Bulldog Curved 38mm	12 X 20
Elevators			
	F894069	Adson Dural Straight 7mm	8 X 10
	F894070	Pennybacker 10mm	8 X 10
Forceps			
	Artery		
	F893025	Crile Curved 14cm (5.5")	12 X 20
	F893026	Crile Straight 14cm (5.5")	12 X 20
	F893032	Heiss Straight 20cm (8")	8 X 10
	Clip		
	F894159	Raney Clip Applicator and Remover 15cm (6")	8 X 10
	Dissecting		
	F893059	Adson Non Toothed 12.5cm (5")	12 X 20
	F893057	Adson Toothed 12.5cm (5")	12 X 20
	F893060	Adson Toothed 1:2 Bayonet 11cm (4.5")	8 X 10
	F893058	Adson Toothed 1:2 Bayonet 19cm (7.5")	8 X 10
	F894158	Adson Toothed Bayonet 19cm (7.5")	8 X 10
	F893064	Gillies Toothed 15cm (6")	8 X 20
	F893083	Waugh Toothed 1:2 20cm (8")	8 X 10
	F893082	Waugh Non Toothed 18cm (7")	8 X 10
	F893287	Adson Fine Non Toothed 12.5cm (5")	12 X 20
	F893286	Adson Fine Toothed 12.5cm (5")	12 X 20
	F893071	McIndoe Non Toothed 15cm (6")	12 X 20
	F893302	McIndoe Toothed 1:2 15cm (6")	8 X 20

Type	Ref No	Description	Case Quantity
Needle Holders			
	F893126	Kilner 13.5cm (5.25")	12 X 20
Punches			
	Rongeur		
	F894068	Kerrison 18cm Shaft 45° 0.3mm Lite	50
Retractors			
	Self Retaining		
	F894160	Cone Hemi-Laminectomy 30cm (12")	25
	F893171	Travers 4:5 Blunt 14cm (5.5")	50
	F893170	Travers 4:5 Blunt 20cm (8")	25
	F893150	Weitlaner Small Curved Scalp 14cm (5.5")	2 X 50
	F893168	Weitlaner 3:4 Blunt Teeth 14cm (5.5")	2 X 50
	F893169	Weitlaner 3:4 Sharp Teeth 14cm (5.5")	2 X 50
	Hand Held		
	F893164	Langenbeck Small 21.5cm (8.5") Blade 25 x 6mm	8 X 10
	F893165	Langenbeck Medium 21.5cm (8.5") Blade 35 x 11mm	50
	F893321	Langenbeck Medium Narrow 21.5cm (8.5") Blade 35 x 8mm	50
	F893322	Langenbeck Medium Wide 21.5cm (8.5") Blade 35 x 15mm	50
	F893166	Langenbeck Large 21.5cm (8.5") Blade 41 x 11mm	50
	F893323	Langenbeck X-Large Long 21.5cm (8.5") Blade 55 x 11mm	50
	F893151	Czerny 17.5cm (7")	25
	F894163	Waltons Malleable Small 20cm (8") x 25mm (1")	8 X 20
	F894164	Waltons Malleable Medium 20cm (8") x 38mm (1.5")	8 X 20
	F894165	Waltons Malleable Large 20cm (8") x 50mm (2")	8 X 20
	Hook		
	F894166	Adson Dissecting Sharp 19.5cm (7.7")	8 X 20
	F894167	Adson Dissecting Blunt 19.5cm (7.7")	8 X 20
	Dissector		
	F893021	MacDonald Double Ended 19cm (7.5")	8 X 10
	F894162	Olivecrona Double Ended 23cm (9")	8 X 10

Type	Ref No	Description	Case Quantity
Saws			
	F894065	Gigli Complete with 2 Handles	2 X 25
Scissors			
	F893204	Micro Curved 9cm (3.6")	8 X 10
	F894161	McKissock Curved 14.5cm (5.5")	12 X 20

Type	Ref No	Description	Case Quantity
ENT			
Bipods			
	F894002	Draffin (Pair) Paediatric 40cm (16")	25
	F894001	Draffin (Pair) Adult 45cm (18")	25
	F894107	Magauran Plate for Draffin Bipod	25
Blades			
Myringotome			
	F894003	Politzer	8 X 10
Curettes			
	F894010	Angel James Non Guarded 12mm	8 X 10
	F894011	Beckmann Adenoid Non Guarded Size 10mm	8 X 10
	F894012	Beckmann Adenoid Non Guarded Size 14 mm	8 X 10
	F894013	Beckmann Adenoid Non Guarded Size 16mm	8 X 10
	F894014	Beckmann Adenoid Non Guarded Size 18mm	8 X 10
	F894016	St Clair Thomson Adenoid Guarded Size 8mm	8 X 10
	F894017	St Clair Thomson Adenoid Guarded Size 10mm	8 X 10
	F894018	St Clair Thomson Adenoid Guarded Size 12mm	8 X 10
	F894019	St Clair Thomson Adenoid Guarded Size 14mm	8 X 10
	F894020	St Clair Thomson Adenoid Guarded Size 16mm	8 X 10
	F894021	St Clair Thomson Adenoid Guarded Size 18mm	8 X 10
Probes			
	F894015	Jobson Horne Ear 17.5cm (7")	8 X 10
	F893015	Macewan 18cm (7")	8 X 10
Cutters			
Wire			
	F896007	Universal Wire 12cm (4.8")	12 X 10

Type	Ref No	Description	Case Quantity
Dilators			
	Tracheal		
	F894024	Bowlby Adult 14cm (5.5")	2 X 25
	F894025	Bowlby Child 10.5cm (4.5")	2 X 25
Dissectors			
	Tonsil		
	F894026	Gwynne Evans 6mm	8 X 10
	F894114	Howarth Sharp/Blunt 21.5cm (8.5")	8 X 10
	F894062	Mitchell Trimmer	8 X 10
Elevators			
	F893022	Howarth Periosteal 23cm (9")	8 X 10
	F894034	Ficklins	8 X 20
	F894035	Freer	8 X 20
	F894036	Hills	8 X 20
Forceps			
	Artery		
	F894079	Birkett Curved 18cm (7")	8 X 10
	F894078	Birkett Straight 18cm (7")	8 X 10
	F894076	Negus Small Curved 19cm (7.5")	8 X 10
	F894077	Negus Large Curved 19cm (7.5")	8 X 10
	F894081	Wilson 19cm (7.5")	50
	Dissecting		
	F893023	Adson Browne Multi Toothed 12cm (4.8")	12 X 20
	F893287	Adson Fine Non Toothed 12.5cm (5")	12 X 20
	F893286	Adson Fine Toothed 12.5cm (5")	12 X 20
	F893064	Gillies Toothed 15cm (6")	8 X 20
	F893071	McIndoe Non Toothed 15cm (6")	12 X 20
	F893302	McIndoe Toothed 1:2 15cm (6")	8 X 20

Type	Ref No	Description	Case Quantity
Forceps (Continued)			
	Dissecting (Continued)		
	F894149	Struempel Nasal Biopsy 14cm (5.5")	12 X 20
	F894112	Tilley Nasal 14cm (5.5")	8 X 10
	F893082	Waugh Non Toothed 18cm (7")	8 X 10
	F893083	Waugh Toothed 1:2 20cm (8")	8 X 10
	F894152	Ormerod Cup Fine 9cm (3.5")	12 X 10
	Dressing		
	F894073	Tilley Aural 14cm (5.5")	8 X 20
	F894080	Lucs Grasping Oval Jaw 14cm (5.5")	2 X 25
	F894118	Hartmann Cup 8.5cm (3.3")	12 X 10
	F894072	Hartmann Crocodile 14cm (5.5")	12 X 10
	F894071	Hartmann Micro Serrated Black 9cm (3.5")	12 X 10
	F893093	Lister Sinus 15cm (6")	12 X 20
	F894115	Cawthorne Crocodile Serrated 10cm (4")	8 X 10
	F894074	Denis Browne Tonsil Holding 20cm (8")	50
Gags & Blades			
	F894085	Boyle Davis Adult and Doughty Blade 8.75cm (3.5")	25
	F894086	Boyle Davis Adult and Doughty Blade 10cm (4")	25
	F894087	Boyle Davis Adult and Doughty Blade 11.25cm (4.5")	25
	F894082	Boyle Davis Child and Doughty Blade 6.25cm (2.5")	25
	F894083	Boyle Davis Child and Doughty Blade 7.5cm (3")	25
	F894084	Boyle Davis Child and Doughty Blade 8.75cm (3.5")	25
	F894088	Doyens Adult	50
	F894089	Doyens Child	50
Handles			
	F893002	Mini Blade	12 X 20

Type	Ref No	Description	Case Quantity
Hooks			
	F894109	Cawthorne 14cm (5.5")	8 X 20
	F894090	St Barts Wax 14cm (5.5")	8 X 20
	F894091	Tracheal Double Blunt 18cm (7")	8 X 20
	F894092	Tracheal Single Blunt 18cm (7")	8 X 20
	F894066	Formby D/E Scoop & Hook	8 X 10
Mirrors			
	F894094	Laryngeal & Handle 8mm	2 X 50
	F894095	Laryngeal & Handle 10mm	2 X 50
	F894096	Laryngeal & Handle 12mm	2 X 50
	F894097	Laryngeal & Handle 14mm	2 X 50
	F894098	Laryngeal & Handle 16mm	2 X 50
	F894099	Laryngeal & Handle 18mm	2 X 50
	F894100	Laryngeal & Handle 20mm	2 X 50
	F894101	Laryngeal & Handle 22mm	2 X 50
	F894102	Laryngeal & Handle 24mm	2 X 50
	F894103	Laryngeal & Handle 26mm	2 X 50
	F894104	Post Nasal & Handle 6mm	5 X 20
	F894105	Post Nasal & Handle 8mm	5 X 20
	F894106	Post Nasal & Handle 10mm	5 X 20
Needle Holders			
	F893126	Kilner 13.5cm (5.25")	12 X 20
Probes			
	F894015	Jobson Horne Ear 17.5cm (7")	8 X 10
	F894108	Jobson Horne 13cm (5")	12 X 20
Punches			
	Rongeur		
	F894068	Kerrison 18cm Shaft 45° 0.3mm	8 X 10

Type	Ref No	Description	Case Quantity
------	--------	-------------	---------------

Pushers



Ligature

F894075	Negus Knot Tier 20cm (8")	8 X 10
---------	---------------------------	--------

Retractors



Hand Held

F894121	Aufricht 16cm (6.25")	8 X 20
F893156	Gillies Skin Hook 18cm (7")	8 X 20
F894119	Hurds Pillar 21.5cm (8.5")	8 X 10
F894120	Mollison Pillar 21.5cm (8.5")	8 X 10
F893160	Kilner Alae 8cm (3.25")	12 X 10
F893161	Kilner Catspaw Double Ended 15cm (6")	8 X 10
F893156	Gillies Skin Hook 18cm (7")	8 X 20
F893157	McIndoes Skin Hook 19cm (7.6")	8 X 10
F893155	Hook Single (Blunt) 15cm (6")	8 X 10
F893153	Hook Skin Double (Blunt) 15cm (6")	8 X 10
F893164	Langenbeck Small 21.5cm (8.5") Blade 25 x 6mm	8 X 10
F893165	Langenbeck Medium 21.5cm (8.5") Blade 35 x 11mm	50
F893321	Langenbeck Medium Narrow 21.5cm (8.5") Blade 35 x 8mm	50
F893322	Langenbeck Medium Wide 21.5cm (8.5") Blade 35 x 15mm	50
F893166	Langenbeck Large 21.5cm (8.5") Blade 41 x 11mm	50
F893323	Langenbeck X-Large Long 21.5cm (8.5") Blade 55 x 11mm	50
F894023	Lacks Tongue Adult	2 X 25
F894022	Lacks Tongue Child	2 X 25



Self Retaining

F893168	Weitlaner 3:4 Blunt Teeth 14cm (5.5")	2 X 50
F893169	Weitlaner 3:4 Sharp Teeth 14cm (5.5")	2 X 50
F893173	West 3:4 Blunt Teeth 14cm (5.5")	50

Scissors



Dissecting

F894123	Aufricht C.O.F 15cm (6")	12 X 10
F893184	Kilner S/S Curved 11cm (4.5")	12 X 20
F893182	Kilner S/S Curved 14cm (5.5")	12 X 20

Type	Ref No	Description	Case Quantity
Scissors (Continued)			
	Dissecting(Continued)		
	F893183	Kilner S/S Straight 11cm (4.5")	12 X 20
	F893181	Kilner S/S Straight Fine 11cm (4.5")	12 X 20
	F896014	Stevens Tenotomy Curved Fine 11cm (4.4")	12 X 20
	F894122	Wilson Tonsil 20cm (8")	8 X 10
	F893214	Woods Tonsil 19cm (7.5")	8 X 10
	Stitch		
	F893188	Littauer Stitch 13.5cm (5.25")	12 X 20
Snares			
	Tonsil		
	F894124	Eve With Wire	25
Specula			
	Nasal		
	F894154	Thudichum Size 0	12 X 10
	F894126	Thudichum Size 1	12 X 10
	F894127	Thudichum Size 2	12 X 10
	F894128	Thudichum Size 3	12 X 10
	F894129	Thudichum Size 4	12 X 10
	F894125	Thudichum Size 5	12 X 10
	F894130	Thudichum Size 6	12 X 10
	F894131	Thudichum Size 7	12 X 10
	Aural		
	F894132	Tumarkin 4.5cm (1.8")	12 X 10
	F894151	Unispec Oval Set of 5 3.5/4.5/5.5/6.5/7.5mm	12 X 10
	F894150	Unispec Round Set of 5 3/4/5/6/7mm	12 X 10
Suckers			
	F894134	Frazier 9fg with VC	8 X 10
	F894153	Frazier 12fg with VC	8 X 10
	F894148	Lempert with VC & Obturator	8 X 10

Type	Ref No	Description	Case Quantity
Suckers (Continued)			
	F894133	Magill Small	8 X 10
	F891040	Yankauer Rose Tip 28cm (11.2")	100
	F894135	Zoellner 2mm Without VC 16cm (6")	2 X 50
	F894136	Zoellner 2mm With VC 16cm (6")	2 X 50
	F894137	Zoellner 2mm Without VC 17.5cm (7")	2 X 50
	F894138	Zoellner 2mm With VC 17.5cm (7")	2 X 50
Syringes			
	F894061	Dental Cartridge	8 X 10
Trocars			
	F894139	Tilley Lichwitz Trocar 8Fg 19.5cm (7.8")	8 X 10

Neurosurgery - Maxillo Facial

Type	Ref No	Description	Case Quantity
Maxillo Facial			
Forceps			
	Dissecting		
	F893064	Gillies Toothed 15cm (6")	8 X 20
Needle Holders			
	F893126	Kilner 13.5cm (5.25")	12 X 20
Retractors			
	Hand Held		
	F893160	Kilner Alae 8cm (3.25")	12 X 10
	F893161	Kilner Catspaw Double Ended 15cm (6")	8 X 10

Type	Ref No	Description	Case Quantity
Dental			
Burnishers			
	F894141	Ball Ended	8 X 20
Carvers			
	F894027	Wards	8 X 20
Chisels			
	F894005	Cottle 10mm x 17.5cm (7")	8 X 10
	F894006	Cottle 12mm x 17.5cm (7")	8 X 10
	F894007	Coupland Size 1	8 X 10
	F894008	Coupland Size 2	8 X 10
	F894009	Coupland Size 3	8 X 10
Condensers			
	F894028	Condenser Amalgam	8 X 20
Elevators			
	F894029	Cryer Left Size 30	8 X 20
	F894030	Cryer Right Size 30	8 X 20
	F894032	Warwick James Left	8 X 20
	F894031	Warwick James Right	8 X 20
	F894033	Warwick James Straight	8 X 20
Forceps			
	Dissecting		
	F894063	College 15cm (6")	8 X 20
	F894156	Oral Fickling 18cm (7.2")	8 X 10
	F894157	Oral Fickling 1x2 Teeth 18cm (7.2")	8 X 10
	Extraction		
	F894038	Dental 123	25
	F894039	Dental 76N	25

Neurosurgery - Dental

Type	Ref No	Description	Case Quantity
Mirrors			
	Dental		
	F894042	No 4	2 X 50
	F894043	No 5	2 X 50
Probes			
	F894059	Dental C.I.T.P.N	8 X 20
	F894058	Dental I4W	8 X 20
	F894044	Dental I54	8 X 20
	F894045	Dental 73	8 X 20
	F894047	Dental Caulx	8 X 20
	F894048	Dental Exec I33	8 X 20
	F894049	Dental H6	8 X 20
	F894050	Dental L9	8 X 20
	F894053	Dental No I52	8 X 20
	F894051	Dental No 5	8 X 20
	F894052	Dental No 6	8 X 20
	F894055	Dental PF179	8 X 20
	F894054	Dental PF1	8 X 20
	F894046	Dental Round 8mm	8 X 20
	F894056	Dental Standard	8 X 20
	F894057	Dental Flat Plastic	8 X 20
	F896033	Dental Swan Neck	12 X 20
Retractors			
	F894144	Minnesota I4cm (5.5") Long 2cm Wide	12 X 5
Scoops			
	F894066	Formby D/E Scoop & Hook	8 X 20
	F894060	Evacuator	8 X 20

Neurosurgery Instrument Sets

Dental Examination Set

Ref No	F887001
UOM	2 X 50
Contains:	
Probe Dental Standard	1
Dental Mirror No 4	1



Dental Filling Set

Ref No	F887002
UOM	20
Contains:	
Mirror Dental No 4	1
Probe Dental Standard	1
Probe Flat Plastic	1
Forceps College 15cm (6")	1
Spoon Excavator	1
Carver Wards	1
Condenser Amalgam 16.5cm (6.5")	1
Burnisher Ball Ended	1
Towel Dressing 43 x 38cm	1
Dental Rolls No 2	6



Tray	1
Sterile Field	1

Dental Temporary Filling Set

Ref No	F887003
UOM	50
Contains:	
Mirror Dental No 4	1
Probe Dental Standard	1
Probe Flat Plastic	1
Forceps College 15cm (6")	1



Neurosurgery Instrument Sets

Neurosurgery Instrument Sets (Continued)

Tracheostomy Set

Ref No	F884001
UOM	10
Contains:	
Forceps Sponge Holding Rampley 24cm (9.5")	2
Clip Towel Cross Action	5
Pin Instrument Mayo 15cm (6")	2
Forceps Dissecting McIndoe Non Toothed 15cm (6")	1
Forceps Dissecting Gillies Toothed 15cm (6")	1
Forceps Artery Halstead Mosquito Straight 12.5cm (5")	5
Forceps Artery Dunhill Curved 12.8cm (5")	5
Forceps Tissue Allis 2:3 Teeth 15cm (6")	1
Handle BP No 3	1

Handle BP No 4	1
Needle Holder Masson 27cm (10.5")	1
Scissors Dressing S/S Straight 12.5cm (5")	1
Scissors Strabismus Straight Blunt 11.5cm (5.6")	1
Scissors McIndoe C.O.F 18cm (7")	1
Hook Skin Single (Blunt) 15cm (6")	1
Gillies Skin Hook 18cm (7")	1
Hook Skin Double (Blunt) 15cm (6")	1
Retractor West Blunt 3:4 Teeth 14cm (5.5")	1
Dilator Tracheal Bowlby Adult	1
Tray	1
Sterile Field	1

Thyroidectomy Set

Ref No	F883025
UOM	10
Contains:	
Forceps Sponge Holding Rampley 24cm (9.5")	3
Clip Towel Cross Action	6
Pin Instrument Mayo 15cm (6")	2
Forceps Dissecting Gillies Toothed 15cm (6")	2
Forceps Dissecting Lane Toothed 1:2 Teeth 12.5cm (5")	2
Forceps Dissecting Bonney Non Toothed 18cm (7")	2
Forceps Artery Halstead Mosquito Straight 12.5cm (5")	10
Forceps Artery Halstead Mosquito Curved 12.5cm (5")	10

Forceps Tissue Allis 2:3 Teeth 15cm (6")	5
Scissors Mayo Straight 17.5cm (7")	2
Scissors Mayo Curved 16.5cm (6.5")	1
Forceps Moynihan Cholecystectomy 22.5cm (9")	2
Forceps Ligature Lahey 20cm (8")	2
Needle Aneurysm Small 15cm (6")	1
Hook Skin Single (Blunt) 15cm (6")	2
Forceps Artery Curved Kocher 1:2 Teeth 21cm (8")	1
Joll Thyroid 15cm (6")	1
Tray	1
Sterile Field	1

Introducing our Ophthalmic Range



Fine Surgical Instruments

Type	Ref No	Description	Case Quantity
Clamps			
	Chalazion		
	F895009	Lambert 10mm Round Plate	12 X 20
	F895020	Lambert 15mm Round Plate	12 X 20
	F895013	Desmarres 13mm Oval Plate	12 X 20
	F895010	Desmarres 16mm Oval Plate	12 X 20
	F895011	Desmarres 18mm Oval Plate	12 X 20
	F895012	Desmarres 20mm Oval Plate	12 X 20
	F895067	Desmarres 31mm Oval Plate	12 X 20
	F895030	Greene Tarsal Cyst 13mm Round	12 X 20
Curettes			
	Spoon		
	F895093	Lambert Size A Round 1.5mm	12 X 20
	F895068	Lambert Size B Round 2mm	12 X 20
	F895002	Lambert Size C Oval 3 x 2.2mm	12 X 20
	F895003	Lambert Size D Oval 3.8 x 2.7mm	12 X 20
	F895004	Lambert Size E Oval 4.6 x 3mm	12 X 20
	F895005	Lambert Size F Oval 5.7 x 3.7mm	12 X 20
	F895066	Meibomian Cyst Double Ended 1mm & 2mm	12 X 20
Dilators			
	Lacrimal		
	F895007	Foster	12 X 20
	F895008	Nettleship	12 X 20
Handles			
	F893002	Mini Blade	12 X 20
Hooks			
	Strabismus		
	F895031	Moorfields 14cm (5.6")	12 X 20
	F895077	Kilner (Gillies) 14cm (5.6")	12 X 20
	F895073	Graefe Small 14cm (5.6")	12 X 20
	F895072	Graefe Medium 14cm (5.6")	12 X 20
	F895071	Graefe Large 14cm (5.6")	12 X 20

Type	Ref No	Description	Case Quantity
Markers			
	F895062	Scleral Double Ended 3.5mm & 4mm	12 X 20
Needle Holders			
Combination			
	F893122	Foster with R/H Scissors 12.5cm (5")	12 X 20
Micro			
	F895037	Barraquer Curved With Lock 10cm (4")	12 X 20
	F895036	Barraquer Curved Without Lock 10cm (4")	12 X 20
	F895034	Castroviejo Curved With Lock 14cm (5.5")	12 X 20
	F895032	Castroviejo Curved Without Lock 14cm (5.5")	12 X 20
	F895035	Castroviejo Straight With Lock 14cm (5.5")	12 X 20
	F895033	Castroviejo Straight Without Lock 14cm (5.5")	12 X 20
	F895038	Silcock With Lock 14.5cm (5.8")	12 X 20
	F895039	Troutman Curved With Lock 14cm (5.5")	12 X 20
Probes			
Lacrimonal			
	F895083	Double Ended 0000/000 Liebreichs	12 X 20
	F895069	Double Ended 00/0 Liebreichs	12 X 20
	F895040	Double Ended 1/2 Liebreichs	12 X 20
	F895085	Double Ended 3/4 Liebreichs	12 X 20
	F895064	Foster Set of Three 0.6mm, 0.8mm, 1.0mm	12 X 20
Forceps			
Eye			
	F895028	Graefe Iris Toothed 1:2 8.75cm (3.5")	12 X 20
	F893065	Iris Non Toothed 10cm (4")	12 X 20
	F893066	Iris Toothed 10cm (4")	12 X 20
	F895019	Iris Toothed Fine 10cm (4")	12 X 20
	F895015	Moorfields Non Toothed 11.5cm (4.6")	12 X 20
	F895014	St Martins 1:2 Teeth and Platform 0.3mm 9cm (3.6")	12 X 20
	F895024	St Martins Non Toothed 9cm (3.6")	12 X 20

Type	Ref No	Description	Case Quantity
Forceps (Continued)			
	Fixation		
	F895070	Capsulorhexis Forceps Long	12 X 20
	F895016	Curved Non Toothed 11cm (4.4")	12 X 20
	F895029	Jayle 1:2 Teeth 11.25cm (4.5")	12 X 20
	F895017	Serrated Non Toothed Straight 12cm (4.8")	12 X 20
	F895018	Straight 1:2 Teeth 11cm (4.4")	12 X 20
	Tying		
	F895026	Birks (Tennant) Micro Curved 10cm (4")	12 X 20
	F895027	Birks (Tennant) Micro Straight 10cm (4")	12 X 20
	Corneal		
	F895022	St Martins 1:2 Teeth and Platform 0.15mm 9cm (3.6")	12 X 20
	F895014	St Martins 1:2 Teeth and Platform 0.3mm 9cm (3.6")	12 X 20
	F895024	St Martins Non Toothed 9cm (3.6")	12 X 20
	F895092	Colibri 0.3mm Teeth, No Platform, 7.5cm (5")	12 X 20
	Suture		
	F895021	Moorfields 10cm (4")	12 X 20
	F895023	Moorfields 12cm (4.8")	12 X 20
	Jeweller		
	F897001	Angled Smooth No 7 11cm (4.5")	12 X 20
	F895074	Notched Forceps 11cm (4.5")	12 X 20
	F895080	Straight Very Fine Smooth No 5 11cm (4.5")	12 X 20
	F895079	Straight Fine Smooth No 4 11cm (4.5")	12 X 20
	F897002	Straight Medium Smooth No 3 11cm (4.5")	12 X 20
	Diathermy		
	F895025	Bipolar Forceps 5" and Lead 3m	50
	Epilation		
	F893087	Epilation 8.75cm (3.5")	12 X 20
	F893088	Epilation 10cm (4")	12 X 20

Type	Ref No	Description	Case Quantity
Retractors			
	Hand Held		
	F893156	Gillies Skin Hook 18cm (7")	8 X 20
	F895042	Watson Orbital Malleable Set of 4	50
	Lid		
	F895094	Desmarres Eye (Small) 12mm	8 X 10
	F895095	Desmarres Eye (Medium) 14mm	8 X 10
	F895041	Desmarres Eye (Large) 16mm	8 X 10
	F895057	Stallard Small	12 X 20
	F895056	Stallard Large	12 X 20
	F895082	Jaffe (Pair)	12 X 20
	Chalazion		
	F895013	Desmarres 13mm Oval Plate	12 X 20
	F895010	Desmarres 16mm Oval Plate	12 X 20
	F895011	Desmarres 18mm Oval Plate	12 X 20
	F895012	Desmarres 20mm Oval Plate	12 X 20
	F895067	Desmarres 31mm Oval Plate	12 X 20
Scissors			
	F895045	Iris C.O.F 11.25cm (4.5")	12 X 20
	F895046	Iris Straight 11cm (4.5")	12 X 20
	F895047	Iris Straight Micro 11cm (4.5")	12 X 20
	F895049	Strabismus Straight Blunt 11.5cm (4.6")	12 X 20
	F895048	Strabismus Curved Blunt 11.5cm (4.6")	12 X 20
	F896014	Stevens Tenotomy Curved Fine 11cm (4.4")	12 X 20
	F895078	Stevens Tenotomy Straight Fine 11cm (4.4")	12 X 20
	F895052	Westcott Tenotomy Straight Sharp 11cm (4.4")	12 X 20
	F895081	Westcott Tenotomy Straight Blunt 11cm (4.4")	12 X 20
	F895044	Enucleation Curved 11cm (4.4")	12 X 20

Type	Ref No	Description	Case Quantity
------	--------	-------------	---------------

Scissors (Continued)

	Micro		
	F895087	Castroviejo Corneal Micro 7mm Blade 10cm (4")	12 X 20
	F895088	Castroviejo Corneal Micro 9mm Blade 10cm (4")	12 X 20
	F895043	Castroviejo Corneal Micro 12mm Blade 10cm (4")	12 X 20
	F895050	Vannas Corneal 9mm Blade Curved 8cm (3.2")	12 X 20
	F895091	Vannas Corneal 12mm Blade Curved 8cm (3.2")	12 X 20
	F895051	Vannas Corneal 9mm Blade Straight 8cm (3.2")	12 X 20
	F895090	Vannas Corneal 12mm Blade Straight 8cm (3.2")	12 X 20

Sizers

	Orbital		
	F895084	Sizer 14mm	12 X 20
	F895086	Sizer 16mm	12 X 20
	F895053	Sizer 18mm	12 X 20
	F895054	Sizer 20mm	12 X 20
	F895055	Sizer 22mm	12 X 20

Spatulas

	F895089	Elliot	12 X 20
---	---------	--------	---------

Specula

	Self Retaining		
	F895061	Barraquer Open 40 x 26.5mm	12 X 20
	F895058	Barraquer Open 44 x 24mm	12 X 20
	F895075	Kratz Barraquer Open 40 x 33mm Angled	12 X 20
	F895076	Kratz Barraquer Open 40 x 33mm Straight	12 X 20

Ophthalmic Instrument Sets

Sub Tenon's Set

Ref No	F885002	Forceps Dissecting St. Martins Non Toothed 9cm (3.6")	
UOM	20	Tray	
Contains:		Sterile Field	
Scissors Westcott 11cm (4.4")			
Needle Holder Barraquer Micro, Fine, Curved, No Lock 10cm (4")			

Chalazion Set

Ref No	F885003	Scalpel BP Handle With No 11 Blade	
UOM	20	Tray	
Contains:		Sterile Field	
Retractor Desmarres Eyelid 16mm			
Curette D/E 2mm and 3mm			

Enucleation Strabismus Set

Ref No	F885004	Scissors Strabismus Curved Blunt 11.5cm (4.6")	
UOM	20	Elliot's Spatula	
Contains:		Forceps Artery Halstead Mosquito Straight 12.5cm (5")	
Forceps Moorfields Eye 10cm (4")		Scissors Westcott 11cm (4.4")	
Forceps Sponge Holding Rampley 18cm (7.2")		Needle Holder Barraquer Micro, Fine, Curved, No Lock 10cm (4")	
Hook Strabismus 14cm (5.6")		Gallipot 60ml	
Needle Holder Castroviejo Curved with Lock 14cm (5.5")		Tray	
Clip Bulldog Silver 60mm		Sterile Field	
Forceps Dissecting St. Martins Non Toothed 9cm (3.6")			

Ophthalmic Procedure Packs

Ophthalmic Procedure Packs

Phaco Pack

Ref No	F500071
UOM	20
Contains:	
Bowl Polypropylene 250ml	1
Gallipot Polypropylene 60ml	1
Swab Non Woven 10 x 10cm 4Ply	10
Swab Visispears Lint Free	5
Knife Slit 2.75mm Angled Bevel	1
Cannula Hydrodissector Nucleus 27g x 7/8" Sauter	1
Cannula Lacrimal 21g x 1 3/8"	1
Cannula Anterior Chamber Rycroft 30g x 22mm	1
Syringe 3ml Luer Lock	1
Eye Shield	1
Drape Arm Sleeve 55cm	2
Drape Universal Ophthalmic with Full Incise, Pouch and Wick 140 x 140cm	1
Set of Drug Labels	1



Sharps Pad	1
Speculum Barraquer	1
Forceps Dissecting Moorfields 12cm (4.6")	1
Scissors Strabismums Straight Blunt 11.5cm (5.6")	1
Scissors Westcott 11cm (4.4")	1
Forceps Sponge Holding Rampley 18cm (7.2")	1
Forceps Colibri 0.3mm, Teeth, No platform, 7.5 cm (5")	1
Tray	1
Sterile Field	1

Sub Tenon's Pack

Ref No	F500072
UOM	20
Contains:	
Swab Gauze 5 x 5cm 8ply	5
Gallipot 60ml	1
Syringe 5ml Luer Lock	1
Cannula Sub Tenon Curved 19g x 25mm	1
Eye Pad 8 x 6cm Non Woven	1
Needle Hypodermic 21g x 1.5" Green	1
Drape 76 x 76cm with Central Adhesive Fenestration 6 x 8cm	1



Forceps Moorfields Eye 10cm (4")	1
Scissors Westcott 11cm (4.4")	1
Speculum Barraquer	1
Tray	1
Sterile Field	1

Ophthalmic Procedure Packs (Continued)

AMD Pack

Ref No F500073

UOM 20

Contains:

Braunstein Fixed Caliper 3.5 & 4mm	1
Speculum Barraquer Adult	1
Forceps Moorfields Eye 10cm (4")	1
Drape Fenestrated 40 x 40cm with Adhesive Fenestration 6cm	1
Scissors Iris Micro Eye Straight 11.25cm (4.5")	1
Swab Non Woven 7.5 x 7.5cm 4ply	5
Gallipot 60ml	2
Forceps Sponge Holding Rampley 18cm (7.2")	1



Swab Absorbent Eye On A Stick	5
Cotton Tipped Applicator Small Head 15cm Long	3
Marker Scleral D/E 3.5mm & 4mm	1
Tray	1
Sterile Field	1

Oculoplastic Pack

Ref No F500074

UOM 20

Contains:

Bowl Polypropylene 250ml	1
Gallipot Polypropylene 60ml	2
Applicator Cotton Tipped Small Head 15cm Long	10
Swab Gauze XRD 10 x 10cm 16ply	10
Forceps Multipurpose Blue	4
Drape Side Adhesive & Reinforced 90 x 90cm	2
Drape Bottom Adhesive & Reinforced 175 x 180cm	1
Sharps Pad	1
Forceps St Martin Tying 1:2 Teeth & Platform	1
Forceps Sponge Holding Rampley 18cm (7.2")	1
Scissors Westcott 11cm (4.4")	1



Forcep Halstead Mosquito Curved Fine Tips 12.5cm (5")	1
Forceps Desmarres Chalazion 13mm Oval	1
Scalpel BP Handle With No 15 Blade	1
Needle Holder Castroviejo Straight 14cm (5.5")	1
Scissors Stevens Tenotomy Curved 11cm (4.5")	1
Speculum Barraquer	1
Tray	1
Sterile Field	1

Introducing our Orthopaedic Range incorporating:-

- Podiatry



Fine Surgical Instruments

Type	Ref No	Description	Case Quantity
Chisels			
	F896047	Duray 10mm	8 X 10
	F896048	Duray 15mm	8 X 10
	F896049	Duray 20mm	8 X 10
	F896050	Duray 25mm	8 X 10
	F896051	Duray 30mm	8 X 10
Curettes			
	F896003	Volkman Small Double Ended Size A 22cm (8.75")	8 X 20
	F896002	Volkman Medium Double Ended Size B 22cm (8.75")	8 X 20
	F896001	Volkman Large Double Ended Size C 22cm (8.75")	8 X 20
	F896017	Volkman Spoon Bone Small Size B 16cm (6.5")	8 X 20
	F896016	Volkman Spoon Bone Medium Size C 16cm (6.5")	8 X 20
Cutters			
	Bone/Nail/Wire		
	F896056	Jansen Double Action 19cm (7.5")	8 X 10
	F896057	Lempert 19cm (7.5")	8 X 10
	F896006	Liston Bone 19cm (7.5")	25
	F896005	Liston Bone 22cm (8.8")	25
	F896059	McIndoe Double Action 19cm (7.5")	8 X 10
	F896004	Stamm Bone and Nail 15cm (6")	25
	F896007	Universal Wire 12cm (4.8")	12 X 10
Dissectors			
	F893021	MacDonald Double Ended 19cm (7.5")	8 X 10
	F896008	Watson Cheyne Elevator 19cm (7.5")	8 X 10
Forceps			
	Dissecting		
	F893059	Adson Non Toothed 12.5cm (5")	12 X 20
	F893060	Adson Toothed 1:2 Bayonet 11cm (4.5")	12 X 10
	F896045	Charnley Toothed 14.5cm (5.75")	8 X 10
	F893064	Gillies Toothed 15cm (6")	8 X 20

Type	Ref No	Description	Case Quantity
Forceps (Continued)			
	Dissecting (Continued)		
	F893067	Lane Non Toothed 18cm (7")	8 X 10
	F893070	Lane Toothed 1:2 Teeth 18cm (7")	8 X 10
	Artery		
	F893029	Halstead Mosquito Curved 12.5cm (5")	12 X 20
	F893030	Halstead Mosquito Straight 12.5cm (5")	12 X 20
	F893291	Kocher Toothed Curved 20cm (8")	8 X 10
	F896009	Kocher Toothed Straight 20cm (8")	8 X 10
	F893039	Spencer Wells Curved 18cm (7")	8 X 10
	F893043	Spencer Wells Straight 19cm (7.6")	8 X 10
	Diathermy		
	F893280	Bipolar 15cm (6") and 3m Lead	50
Mallets			
	Heavy/Fine		
	F896052	Heath 4.5cm (1.8") (793gms)	25
	F896053	Toffee Hammer 19cm (7") (226gms)	2 X 25
Needle Holders			
	F893117	Crile Wood 15cm (6")	12 X 20
	F893119	Crile Wood 20cm (8")	8 X 10
	F893123	Gillies 16cm (6.25")	8 X 10
	F893126	Kilner 13.5cm (5.25")	12 X 20
	F893128	Mayo 18cm (7")	8 X 10
	F893134	Mayo Hegar 18cm (7")	8 X 10
	F893137	Nievert 13cm (5.2")	8 X 10
	F893139	Thomson Walker 20cm (8")	8 X 10
Osteotomes			
	F896013	Nievert Left 20cm (8")	8 X 10
	F896011	Nievert Right 20cm (8")	8 X 10
	F896012	Nievert Straight 20cm (8")	8 X 10

Type	Ref No	Description	Case Quantity
Punches			
	Rongeur		
	F894068	Kerrison 18cm Shaft 45° 0.3mm	8 X 10
Retractors			
	Self Retaining		
	F897003	Alms Blunt 70mm	12 X 20
	F893170	Travers 4:5 Blunt 20cm (8")	25
	F893168	Weitlaner 3:4 Blunt Teeth 14cm (5.5")	2 X 50
	F893169	Weitlaner 3:4 Sharp Teeth 14cm (5.5")	2 X 50
	F893150	Weitlaner Small Curved Scalp 14cm (5.5")	2 X 50
	F893173	West Blunt 3:4 Teeth 14cm (5.5")	50
	Hand Held		
	F893161	Kilner Catspaw Double Ended 15cm (6")	8 X 10
	F893164	Langenbeck Small 21.5cm (8.5") Blade 25 x 6mm	8 X 10
	F893165	Langenbeck Medium 21.5cm (8.5") Blade 35 x 11mm	50
	F893321	Langenbeck Medium Narrow 21.5cm (8.5") Blade 35 x 8mm	50
	F893322	Langenbeck Medium Wide 21.5cm (8.5") Blade 35 x 15mm	50
	F893166	Langenbeck Large 21.5cm (8.5") Blade 41 x 11mm	50
	F893323	Langenbeck X-Large Long 21.5cm (8.5") Blade 55 x 11mm	50
	F894143	Rake 17cm (6.7") long 5.5cm head	8 X 5
Elevators			
	Bone		
	F896055	Lane 27cm (10.5")	50
	F896058	Bristow 22cm (8.75")	50
	F893022	Howarth Periosteal 23cm (9")	8 X 10
	Hook		
	F896054	Harrison Tendon Blunt	8 X 10
	F896015	Trethowan Bone Hook	8 X 10
Saws			
	F894065	Gigli Complete with 2 Handles	2 X 25



Type	Ref No	Description	Case Quantity
Scissors			
	Dressing		
	F893185	Lister Bandage 14cm (5.5")	12 X 10
	F893187	Lister Bandage 15cm (6")	8 X 10
	F893186	Lister Bandage 18.75cm (7.5")	8 X 10
	Dissecting		
	F893191	Mayo C.O.F 19cm (7.5")	8 X 10
	F893213	Mayo Straight 19cm (7.5")	8 X 10
	F893200	McIndoe C.O.F 19cm (7.5")	8 X 10
	F893301	McIndoe C.O.F Tungsten Carbide 23cm (7")	8 X 10
	F893300	Metzenbaum C.O.F Tungsten Carbide 19cm (7.5")	8 X 10
	F893202	Metzenbaum C.O.F 19cm (7.5")	8 X 10
	F893201	Metzenbaum Straight 14cm (5.5")	8 X 10

Type	Ref No	Description	Case Quantity
------	--------	-------------	---------------

Podiatry

Applicators

	F896030	Tubegauze Medium	12 X 10
	F896031	Tubegauze Large	12 X 10

Burrs

	Podiatry		
	F896043	Large	12 X 20
	F896034	Medium	12 X 20
	F896042	Small	12 X 20

Cutters

	Nail		
	F896018	Nail 10cm (4")	12 X 10
	F896019	Nail Barrel Spring 13.75cm (5.5")	8 X 10
	F896032	Nail Fine 12cm (4.7")	8 X 10
	F896040	Nail Curved Serrated Grips with Lock 14cm (5.5")	12 X 10
	F896035	Nail Curved Smooth Handles Leaf Spring 14cm (5.5")	8 X 10
	F896037	Nail Straight Smooth Handles Leaf Spring 12cm (4.7")	8 X 10
	F896020	Nail Curved Smooth Handles Barrel Spring 13.5cm (5.25")	8 X 10
	F896036	Nail Straight Smooth Handles Barrel Spring 14cm (5.5")	8 X 10
	F896041	Nail Straight Corner Nipper 13cm (5")	12 X 10
	F896038	Nail Straight Silicone Grips 14cm (5.5")	8 X 10

Dressers

	F896021	Foot 20cm (8")	8 X 20
---	---------	----------------	--------

Elevators

	F896022	Nail 13.5cm (5.25")	12 X 20
---	---------	---------------------	---------




Type	Ref No	Description	Case Quantity
Files			
	F896024	Blacks Nail Fine 13.5cm (5.25")	12 X 20
	F896023	Blacks Nail Medium 13.5cm (5.25")	12 X 20
	F896026	Corn Rasp 15cm (6")	8 X 10
Handles			
	F893002	Mini Blade	12 X 20
Scissors			
	F893210	Toe Nail Serrated	12 X 20
Splitters			
	F896028	Anvil Nail 13.5cm (5.25")	12 X 20
	F896029	Thwaites Nail Splitter 13cm (5")	12 X 10

Orthopaedic Instrument Sets

Orthopaedic Instrument Sets

Standard Podiatry Set

Ref No	F886005
UOM	50
Contains:	
Handle BP No 3 Graduated	1
Cutter Nail Smooth Handles Leaf Spring Silicone Sleeves	1
Dresser Foot 20cm (8")	1
File Nail Fine (Blacks) 13.5cm (5.25")	1



Podiatry Assistant Set

Ref No	F886015
UOM	2 X 20
Contains:	
Diamond Deb Foot Dresser 20cm (8")	1
File Nail (Blacks) 13.5cm (5.25")	1
Cutter Nail Smooth Silicone Handles Straight Leaf Spring 12cm (4.7")	1



Podiatry Set

Ref No	F886016
UOM	50
Contains:	
Diamond Deb Foot Dresser 20cm (8")	1
Cutter Nail Smooth Silicone Handles Curved Leaf Spring 12cm (4.7")	1





Orthopaedic Instrument Sets (Continued)

Basic Podiatry Set

Ref No	F886010
UOM	30
Contains:	
Handle BP No 3 Graduated	1
Cutter Nail Smooth Silicone Handles	1
Curved Leaf Spring 12cm (4.7")	1
Diamond Deb Foot Dresser 20cm (8")	1
File Nail Fine (Blacks) 13.5cm (5.25")	1
Tray	1



Sterile Field	1
---------------	---

Nail Surgery Set

Ref No	F886007
UOM	45
Contains:	
Nail Splitter Nipper Thwaites Silicone Grips	1
Elevator Nail	1
File Nail Fine (Blacks) 13.5cm (5.25")	1
Scissors Dressing S/B Straight 12.5cm (5")	1
Forceps Artery Halstead Mosquito Straight 12.5cm (5")	1
Wooden Applicator 15cm	5
Tray	1



Sterile Field	1
---------------	---

Orthopaedic Instrument Sets (Continued)

Nail Surgery Set

Ref No	F886014
UOM	25
Contains:	
Handle BP No 3 Graduated	1
Mini Blade Handle	1
Forceps Artery Halstead Mosquito Straight 12.5cm (5")	1
Elevator Nail	1
File Nail Medium (Blacks) 13.5cm (5.25")	1
Nail Splitter Nipper Thwaites	1



Applicator Tubegauze (Large)	1
Tray	1
Sterile Field	1



Orthopaedic Procedure Packs

Podiatry Debridement Pack

Ref No	F886004
UOM	14
Contains:	
Forceps Artery Halstead Mosquito Straight 12.5cm (5")	1
Forceps Dissecting Toothed Treves 12.5cm (5")	1
Handle BP No 3 Graduated	1
Probe Nickel Silver With Eye 12.5cm (5")	1
Scissors Iris Curved 11.25cm (4.5")	1
Blade Carbon Steel No 10	1



Towel Dressing 43 x 38cm	1
Tray	1
Sterile Field	1

Basic Podiatry Pack

Ref No	F886001
UOM	14
Contains:	
Scissors Toe Nail Serrated	1
Cutter Nail Curved Leaf Spring	1
Dresser Foot 20cm (8")	1
File Nail Medium (Blacks) 13.5cm (5.25")	1
Scalpel BP Handle With No 10 Blade	1
Towel Dressing 43 x 38cm	1



Tray	1
Sterile Field	1

Orthopaedic Procedure Packs (Continued)

Corn Removal Pack

Ref No	F886002
UOM	14
Contains:	
Scissors Toe Nail Serrated	1
Cutter Nail Curved Leaf Spring	1
Dresser Foot 20cm (8")	1
File Nail Medium (Blacks) 13.5cm (5.25")	1
Rasp Corn	1
Scalpel BP Handle With No 10 Blade	1
Towel Dressing 43 x 38cm	1



Tray	1
Sterile Field	1

Zadeks Procedure Pack

Ref No	F886003
UOM	14
Contains:	
Blade Mini Handle	1
Forceps Artery Halstead Mosquito Straight 12.5cm (5")	1
Scissors Toe Nail Serrated	1
Spatula Double Ended	1
File Nail Medium (Blacks) 13.5cm (5.25")	1
Nail Splitter Nipper Thwaites	1
Blade Carbon Steel No 10	1
Applicator Tubegauze (Medium)	1
Swab Non-Woven 7.5 x 7.5cm	5



Towel Dressing 43 x 38cm	1
Tray	1
Sterile Field	1

Introducing our Plastic Range



Fine Surgical Instruments

Type	Ref No	Description	Case Quantity
Elevators			
	Periosteal		
	F893022	Howarth 23cm (9")	8 X 10
Forceps			
	Artery		
	F893025	Crile Curved 14cm (5.5")	12 X 20
	F893026	Crile Straight 14cm (5.5")	12 X 20
	Tissue		
	F893099	Allis American Pattern 14.5cm (5.75")	8 X 20
	F893097	Allis Toothed 2:3 15cm (6")	8 X 20
	F893098	Allis Toothed 3:4 15cm (6")	8 X 20
	F893096	Allis Toothed 4:5 15cm (6")	8 X 20
	Dissecting		
	F893059	Adson Non Toothed 12.5cm (5")	12 X 20
	F893057	Adson Toothed 12.5cm (5")	12 X 20
	F893287	Adson Fine Non Toothed 12.5cm (5")	12 X 20
	F893286	Adson Fine Toothed 12.5cm (5")	12 X 20
	F893060	Adson Toothed 1:2 Bayonet 11cm (4.5")	12 X 10
	F893058	Adson Toothed 1:2 Bayonet 19cm (7.5")	8 X 10
	F893064	Gillies Toothed 15cm (6")	8 X 20
	F893082	Waugh Non Toothed 18cm (7")	8 X 10
	F893083	Waugh Toothed 1:2 20cm (8")	8 X 10
	F893071	McIndoe Non Toothed 15cm (6")	12 X 20
	F893302	McIndoe Toothed 1:2 15cm (6")	8 X 20
	Dressing		
	F893093	Lister Sinus 15cm (6")	12 X 20
Handles			
	F893002	Mini Blade	12 X 20
	F893110	Scalpel Barron	12 X 20
	F897004	Skin Graft Boards	25

Type	Ref No	Description	Case Quantity
Needle Holders			
	F893117	Crile Wood 15cm (6")	12 X 20
	F893119	Crile Wood 20cm (8")	8 X 10
	F893120	DeBakey 18cm (7")	8 X 10
	F893121	Derf 12cm (4.75")	12 X 20
	F893118	Converse 12.5cm (5")	12 X 20
	F893122	Foster with R/H Scissors 12.5cm (5")	12 X 20
	F893123	Gillies 16cm (6.25")	8 X 10
	F893124	Halsey 12.5cm (5")	12 X 20
	F893125	Halsey Tungsten Carbide 15cm (6")	12 X 20
	F893126	Kilner 13.5cm (5.25")	12 X 20
	F893137	Nievert 13cm (5.2")	8 X 10
Retractors			
	Hand Held		
	F893161	Kilner Catspaw Double Ended 15cm (6")	8 X 10
	Hook		
	F893156	Gillies Skin 18cm (7")	8 X 20
	F893157	McIndoe Skin 19cm (7.6")	8 X 10
	F893155	Skin Single (Blunt) 15cm (6")	8 X 10
	F893153	Skin Double (Blunt) 15cm (6")	8 X 10
	Self Retaining		
	F893150	Weitlaner Small Curved Scalp 14cm (5.5")	2 X 50
	F893160	Kilner Alae 8cm (3.25")	12 X 10
	F896046	Alms Blunt 70mm	8 X 10
Scissors			
	Dissecting		
	F893182	Kilner S/S Curved 11cm (4.5")	12 X 20
	F893184	Kilner S/S Curved 14cm (5.5")	12 X 20
	F893183	Kilner S/S Straight 11cm (4.5")	12 X 20
	F893181	Kilner S/S Straight Fine 11cm (4.5")	12 X 20
	F893198	McIndoe C.O.F 12.5cm (5")	12 X 20
	F893199	McIndoe C.O.F 18cm (7")	8 X 10
	F893200	McIndoe C.O.F 19cm (7.5")	8 X 10



Type	Ref No	Description	Case Quantity
Syringes			
	F894061	Dental Cartridge	8 X 10



Plastic Instrument Sets

Skin Graft Set

Ref No	F887005
UOM	10
Contains:	
Handle Scalpel Barron	1
Dissector McDonald D/E 19cm (7.5")	1
Forceps Dissecting McIndoe Non Toothed 15cm (6")	2
Forceps Dissecting Gillies Toothed 15cm (6")	2

Scissors Metzenbaum Straight 14cm (5.5")	1
Clip Towel Backhaus 11cm (4.5")	5
Hook Skin Single (Blunt) 15cm (6")	1
Forceps Sponge Holding Rampley 24cm (9.5")	2
Tray	1
Sterile Field	1

Skin Biopsy Set

Ref No	F887006
UOM	10
Contains:	
Forceps Artery Halstead Mosquito Curved 12.5cm (5")	2
Forceps Dissecting McIndoe Non Toothed 15cm (6")	1
Forceps Dissecting Gillies Toothed 15cm (6")	1

Forceps Dissecting T.O.E. Serrated 12.5cm (5")	1
Hook Skin Single (Blunt) 15cm (6")	1
Scissor/Needle Holder Fosters R/H 12.5cm (5")	1
Scissors Iris Straight 11cm (4.4")	1
Handle BP No 3	1
Tray	1
Sterile Field	1

Suture Set

Ref No	F883026
UOM	10
Contains:	
Forceps Artery Halstead Mosquito Straight 12.5cm (5")	1
Forceps Artery Halstead Mosquito Curved 12.5cm (5")	1
Forceps Dissecting T.O.E. Serrated 12.5cm (5")	1

Forceps Dissecting Toothed Treves 12.5cm (5")	1
Scissors Dressing S/S Straight 12.5cm (5")	1
Needle Holder Kilner 13cm (5.2")	1
Tray	1
Sterile Field	1

Plastic Procedure Packs

Maxillo Facial Pack

Ref No	F700083
UOM	7

Contains:

Swab Gauze XRD 10 x 7.5cm 32ply	5
Needle Solo Supra Dental 27g	1
Skin Marker Pen Fine/Bold & Ruler	1
Needle Holder Halsey Tungsten Carbide 12.5cm (5")	1
Diathermy Bipolar Forceps 15cm (6") and 3m Lead	1
Forcep Sponge Holding 19cm (7.5")	1
Drape Ortho Bar S/A 254 x 152cm	1
Drape Tape 10 x 55cm	1
Drape Utility N/A 112 x 152cm	2
Scissors Iris Straight 11cm (4.4")	1
Forceps Dissecting Iris Toothed 10cm (4")	1
Scissors Kilner Curved 11cm (4.4")	1



Blade Carbon Steel No 15	1
Forceps Adson Dissecting Fine Toothed 12.5cm (5")	1
Retractor Hook Skin Gillies Fine 18cm (7.2")	2
Syringe Dental Cartridge	1
Sharps Box 30 Count Foam Strips	1
Handle BP No 3	1
Tray	1
Sterile Field	1

Plastic Procedure Packs (Continued)

Rhinoplasty Pack

Ref No F700084

UOM 5

Contains:

Forceps Dissecting Adson Browne Multi
Toothed 12cm (4.7") 1

Retractor Aufricht 16cm (6.4") 1

Retractor Kilner Catspaw D/E 15cm (6") 1

Forceps Artery Kocher Toothed Straight
14cm (5.5") 1

Toffee Hammer 19cm (7") (226gms) 1

Diathermy Bipolar Forceps 11.25cm (4.5")
and 3m Lead 1

Forceps Dressing Aural Tilley 14cm (5.5") 1

Retractor Kilner Alae 8cm (3.25") 1

Speculum Killian Septum with Set Screw 1

Scissors Cottle Dorsal 15cm (6") 1

Scissors Cottle Dorsal 14cm (5.5") 1

Forcep Sponge Holding 19cm (7.5") 2

Needle Holder Crile Wood Tungsten Carbide
17.7cm (7") 1

Scissors Iris Straight 11cm (4.4") 1



Scissors Stevens Tenotomy Curved Fine
11cm (4.4") 1

Clip Towel Ball & Socket 2

Forceps Artery Halstead Mosquito Curved
12.5cm (5") 2

Forceps Adson Dissecting Fine Toothed
12.5cm (5") 1

Retractor Hook Skin Gillies Fine 18cm (7.2") 2

Scissors Dressing S/B Straight 12.5cm (5") 1

Syringe Dental Cartridge 1

Tray 1

Sterile Field 1

Plastic Procedure Packs (Continued)

Abdominal Pack

Ref No	F300945
UOM	5
Contains:	
Swab Gauze XRD 30 x 30cm 12ply Taped	10
Swab Gauze XRD 10 x 7.5cm 32ply	10
Needle Counter 15 Count Magnet & Foam Strips	1
Blade Carbon Steel No 15	2
Blade Carbon Steel No 11	1
Needle Hypodermic 19g x 1.5" White	1
Needle Hypodermic 21g x 1.5" Green	1
Syringe 20ml Luer Slip Eccentric	1
Tip Cleaner	1
Skin Closure Steri-Strip 12mm x 10cm (1 Sheet of 3)	2
Dressing Primapore 8.3 x 6cm	2
Towel Scrimex 48 x 40cm	4
Cover Light Handle	2
Blade Carbon Steel No 10	2
Gown Large No Reinforcement	2
Drape Top S/A 150 x 240cm Reinforced with Tube Holders	1
Drape Bottom S/A 175 x 180cm Reinforced with Tube Holders	1
Drape Side S/A 90 x 90cm Reinforced	2
Chloraprep 10.5ml	2
Skin Marker Pen Regular & Ruler	1
Scissors Metzenbaum C.O.F 19cm (7.5")	1
Forcep Sponge Holding 19cm (7.5")	2



Diathermy Quiver	1
Forceps Artery Kocher Toothed Straight 14cm (5.5")	2
Retractor Kelly 21.5cm (1.5") Blade 38mm	1
Handle BP No 3	1
Monopolar Diathermy Forcep 17.8cm (7") 8mm (Bovie) With Lead 3m	1
Needle Holder Mayo 20cm (8")	1
Clip Towel Ball & Socket	2
Forceps Artery Halstead Mosquito Curved 12.5cm (5")	2
Forceps Adson Dissecting Fine Toothed 12.5cm (5")	2
Handle Scalpel Barron	2
Needle Holder Crile Wood 15cm (6")	2
Retractor Skin Hook Gillies Fine 18cm (7.2")	2
Retractor Langenbeck Medium 21.5cm (8.5") Blade 35 x 11mm	2
Scissors Dressing S/B Straight 12.5cm (5")	1
Tray	1
Sterile Field	2

Plastic Procedure Packs (Continued)

Basic Minor Op Pack

Ref No	F883005
UOM	30
Contains:	
Forceps Artery Halstead Mosquito Straight 12.5cm (5")	
Forceps Dissecting Gillies Toothed 15cm (6")	
Forceps Dissecting McIndoe Non Toothed 15cm (6")	
Handle BP No 3 Graduated	
Needle Holder Kilner 13cm (5.2")	
Scissors Iris Curved 11.25cm (4.5")	
Scissors Iris Straight 11cm (4.4")	
Blade Carbon Steel No 15	



Towel Dressing 43 x 38cm	2
Swab Non Woven 7.5 x 7.5cm 4ply	5
Drape Fenestrated 60 x 75cm Off Centre S/A Fen 10cm	
Tray	
Sterile Field	

Standard Minor Op Pack

Ref No	F883006
UOM	18
Contains:	
Forceps Artery Halstead Mosquito Curved 12.5cm (5")	
Forceps Artery Halstead Mosquito Straight 12.5cm (5")	
Forceps Dissecting Adson Non Toothed 12.5cm (5")	
Forceps Dissecting Adson Toothed 12.5cm (5")	
Handle BP No 3 Graduated	
Needle Holder Lawrence 15cm (6")	
Retractor Skin Hook Gillies Fine 18cm (7.2")	
Scissors Dressing B/B Straight 12.5cm (5")	
Scissors Iris Straight 11cm (4.4")	



Scissors Strabismus Curved Blunt 11.5cm (4.6")	
Blade Carbon Steel No 15	
Towel Scrimex 48 x 40cm	2
Swab Gauze 10 x 10cm 12ply	5
Balls Non Woven Medium	3
Drape Small S/A Utility 76 x 76cm	2
Tray	
Sterile Field	

Plastic Procedure Packs (Continued)

Minor Op Pack

Ref No	F883013
UOM	30
Contains:	
Forceps Dissecting Adson Toothed Ex Fine Micro 12.5cm (5")	1
Handle BP No 3 Graduated	1
Needle Holder Halsey Serrated 12.5cm (5")	1
Scissors Iris Straight 11cm (4.4")	1
Scissors Strabismus Curved Blunt 11.5cm (4.6")	1
Towel Dressing 43 x 38cm	1
Swab Non Woven 7.5 x 7.5cm	10



Tray	1
Sterile Field	1

Minor Op Pack

Ref No	F883014
UOM	25
Contains:	
Forceps Artery Halstead Mosquito Curved 12.5cm (5")	2
Forceps Dissecting Adson Non Toothed 12.5cm (5")	1
Handle BP No 3 Graduated	1
Needle Holder Kilner 13cm (5.2")	1
Scissors Iris Straight 11cm (4.4")	1
Blade Carbon Steel No 11	1
Swab Non Woven 5 x 5cm 4ply	5
Drape Fenestrated Deluxe N/A 70 x 90cm with 10cm Fen	1



Tray	1
Sterile Field	1

alphabetical full index

Product	Page
Applicator Tubegauze Large	113
Applicator Tubegauze Medium	113
Bipod Draffin (Pair) Adult 45cm (18")	84
Bipod Draffin (Pair) Child / Paediatric 40cm (16")	84
Bipod Magauran Plate for Draffin Bipod	84
Blade Myringotome Politzer	84
BP Handle Mini Blade	33, 86, 99, 114, 123
BP Handle Scalpel BP No 3	33
BP Handle Scalpel BP No 4	33
BP Handle Scalpel BP No 5	33
BP Handle Scalpel BP No 7	33
Brace Hudson with 14mm Perforator and 16mm Burr	81
Burnisher Ball Ended	92
Burr Podiatry Small	113
Burr Podiatry Medium	113
Burr Podiatry Large	113
Cannula Heparin Small Plain End 2mm	23
Cannula Heparin Medium Plain End 2.5mm	23
Cannula Uterine Currie	51
Cannula Uterine Leech Wilkinson Small	51
Cannula Uterine Leech Wilkinson Medium	51
Cannula Uterine Leech Wilkinson Large	51
Cannula Uterine Spackman 36cm (14.4")	51
Carver Wards	92
Catheter Silver Female 14fg	51
Chisel Cottle 10mm x 17.5cm (7")	92
Chisel Cottle 12mm x 17.5cm (7")	92
Chisel Coupland Size 1	92
Chisel Coupland Size 2	92
Chisel Coupland Size 3	92
Chisel Duray 10mm	109
Chisel Duray 15mm	109
Chisel Duray 20mm	109
Chisel Duray 25mm	109
Chisel Duray 30mm	109
Clamp Chalazion Desmarres 13mm Oval Plate	99
Clamp Chalazion Desmarres 16mm Oval Plate	99
Clamp Chalazion Desmarres 18mm Oval Plate	99
Clamp Chalazion Desmarres 20mm Oval Plate	99
Clamp Chalazion Desmarres 31mm Oval Plate	99
Clamp Chalazion Greene Tarsal Cyst 13mm Round	99
Clamp Chalazion Lambert 10mm Round Plate	99
Clamp Chalazion Lambert 15mm Round Plate	99
Clip Bulldog Serrafine Curved 38mm	81

Product	Page
Condenser Amalgam	92
Curette Angel James Non Guarded 12mm	84
Curette Beckmann Adenoid Non Guarded Size 10mm	84
Curette Beckmann Adenoid Non Guarded Size 14mm	84
Curette Beckmann Adenoid Non Guarded Size 16mm	84
Curette Beckmann Adenoid Non Guarded Size 18mm	84
Curette Sims Small Double Ended Sharp/Blunt 15cm (6")	51
Curette Sims Medium Double Ended Sharp/Blunt 16cm (6.4")	51
Curette Sims Large Double Ended Sharp/Blunt 19cm (7.4")	51
Curette Sims Small Sharp 15cm (6")	51
Curette Sims Medium Sharp 16cm (6.4")	51
Curette Sims Large Sharp 19cm (7.4")	51
Curette St Clair Thomson Adenoid Guarded Size 10mm	84
Curette St Clair Thomson Adenoid Guarded Size 12mm	84
Curette St Clair Thomson Adenoid Guarded Size 14mm	84
Curette St Clair Thomson Adenoid Guarded Size 16mm	84
Curette St Clair Thomson Adenoid Guarded Size 18mm	84
Curette St Clair Thomson Adenoid Guarded Size 8mm	84
Curette Volkmann Small Double Ended Size A 22cm (8.75")	109
Curette Volkmann Medium Double Ended Size B 22cm (8.75")	109
Curette Volkmann Large Double Ended Size C 22cm (8.75")	109
Curette Volkmann Spoon Bone Small Size B 16cm (6.5")	109
Curette Volkmann Spoon Bone Medium Size C 16cm (6.5")	109
Curette Probes Jobson Horne Ear 17.5cm (7")	84
Curette Probes Macewan 18cm (7")	84
Curette Spoon Lambert Size A Round 1.5mm	99
Curette Spoon Lambert Size B Round 2mm	99
Curette Spoon Lambert Size C Oval 3 x 2.2mm	99
Curette Spoon Lambert Size D Oval 3.8 x 2.7mm	99
Curette Spoon Lambert Size E Oval 4.6 x 3mm	99
Curette Spoon Lambert Size F Oval 5.7 x 3.7mm	99
Curette Spoon Meibomian Cyst Double Ended 1mm & 2mm	99
Cutter Bone Jansen Double Action 19cm (7.5")	109
Cutter Bone Lempert 19cm (7.5")	109
Cutter Bone Liston Bone 19cm (7.5")	109
Cutter Bone Liston Bone 22cm (8.8")	109
Cutter Bone McIndoe Double Action 19cm (7.5")	109
Cutter Bone and Nail Stamm 15cm (6")	109
Cutter Wire Universal 12cm (4.8")	23, 40, 84, 109
Cutter Nail 10cm (4")	113
Cutter Nail Barrel Spring 13.75cm (5.5")	113
Cutter Nail Curved Serrated Grips with Lock 14cm (5.5")	113
Cutter Nail Curved Smooth Handles Barrel Spring 13.5cm (5.25")	113

alphabetical full index

Product	Page
Cutter Nail Curved Smooth Handles Leaf Spring 14cm (5.5")	113
Cutter Nail Fine 12cm (4.7")	113
Cutter Nail Straight Corner Nipper 13cm (5")	113
Cutter Nail Straight Silicone Grips 14cm (5.5")	113
Cutter Nail Straight Smooth Handles Barrel Spring 14cm (5.5")	113
Cutter Nail Straight Smooth Handles Leaf Spring 12cm (4.7")	113
Dilator Umbilical	33
Dilator Cervical Hegar Double Ended 1/2mm 22cm (8.8")	51
Dilator Cervical Hegar Double Ended 3/4mm 22cm (8.8")	51
Dilator Cervical Hegar Double Ended 5/6mm 22cm (8.8")	51
Dilator Cervical Hegar Double Ended 7/8mm 22cm (8.8")	51
Dilator Cervical Hegar Double Ended 9/10mm 22cm (8.8")	51
Dilator Cervical Hegar Double Ended 11/12mm 22cm (8.8")	51
Dilator Cervical Hegar Double Ended 13/14mm 22cm (8.8")	51
Dilator Cervical Hegar Double Ended 15/16mm 22cm (8.8")	51
Dilator Cervical Hegar Double Ended 17/18mm 22cm (8.8")	51
Dilator Cervical Hegar Double Ended Set of 9, 1/2... 17/18mm 22cm (8.8")	51
Dilator Umbilical 12.5cm (5")	75
Dilator Lacrimal Foster	99
Dilator Lacrimal Nettleship	99
Dilator Tracheal Bowlby Adult 14cm (5.5")	17, 85
Dilator Tracheal Bowlby Child 10.5cm (4.5")	17, 85
Director Fistula With Guard 13cm (5.2")	33
Director Lockhart Mummery Set of 4	33
Dissector MacDonald Double Ended 19cm (7.5")	33, 82, 109
Retractor Dissector Olivecrona Double Ended 23cm (9")	82
Dissector Watson Cheyne Elevator 19cm (7.5")	33, 109
Dissector Tonsil Gwynne Evans 6mm	85
Dissector Tonsil Howarth Sharp/Blunt 21.5cm (8.5")	85
Dissector Tonsil Mitchell Trimmer	85
Dresser Foot 20cm (8")	113
Elevator Adson Dural Straight 7mm	81
Elevator Cryer Left Size 30	92
Elevator Cryer Right Size 30	92
Elevator Ficklings	85
Elevator Freer	85
Elevator Hills	85
Elevator Howarth Periosteal 23cm (9")	85
Elevator Nail 13.5cm (5.25")	113
Elevator Pennybacker 10mm	81
Elevator Warwick James Left	92
Elevator Warwick James Right	92

Product	Page
Elevator Warwick James Straight	92
Elevator Bone Bristow 22cm (8.75")	111
Elevator Bone Howarth Periosteal 23cm (9")	111
Elevator Bone Lane 27cm (10.5")	111
Elevator Hook Harrison Tendon Blunt	111
Elevator Hook Trethowan Bone Hook	111
Elevator Periosteal Howarth 23cm (9")	123
File Blacks Nail Fine 13.5cm (5.25")	114
File Blacks Nail Medium 13.5cm (5.25")	114
File Corn Rasp 15cm (6")	114
Forcep Artery Adson Straight 18.5cm (7.4")	33
Forcep Artery Birkett Curved 18cm (7")	85
Forcep Artery Birkett Straight 18cm (7")	85
Forcep Artery Cairn 14cm (5.5")	33
Forcep Artery Crile Curved 14cm (5.5")	34, 81, 123
Forcep Artery Crile Straight 14cm (5.5")	34, 81, 123
Forcep Artery Dunhill Curved 12.5cm (5")	34
Forcep Artery Dunhill Straight 12.5cm (5")	34
Forcep Artery Grey Turner Curved 14.5cm (5.75")	34
Forcep Artery Grey Turner Curved 18cm (7")	34
Forcep Artery Halstead Mosquito Curved 12.5cm (5")	34, 110
Forcep Artery Halstead Mosquito Straight 12.5cm (5")	34, 110
Forcep Artery Heiss Straight 20cm (8")	23, 34, 81
Forcep Artery Kocher Toothed Curved 14cm (5.5")	34
Forcep Artery Kocher Toothed Curved 20cm (8")	34, 110
Forcep Artery Kocher Toothed Straight 14cm (5.5")	34
Forcep Artery Kocher Toothed Straight 20cm (8")	34, 110
Forcep Artery Lahey Ligature 20cm (8")	34
Forcep Artery Lahey Ligature 23cm (9")	34
Forcep Artery Lane Toothed 1:2 Teeth 18cm (7")	110
Forcep Artery Moynihan Cholecystectomy 22.5cm (9")	34
Forcep Artery Negus Small Curved 19cm (7.5")	85
Forcep Artery Negus Large Curved 19cm (7.5")	85
Forcep Artery O'Shaugnessy 20cm (8")	23
Forcep Artery Pean Curved 23cm (9")	23
Forcep Artery Pean Curved 25.5cm (10")	23
Forcep Artery Pean Curved 30cm (12")	23
Forcep Artery Pean Straight 23cm (9")	23
Forcep Artery Roberts Curved 22.5cm (9")	34
Forcep Artery Roberts Curved 23cm (9")	23
Forcep Artery Spencer Wells Curved 12.5cm (5")	34
Forcep Artery Spencer Wells Curved 15cm (6")	34
Forcep Artery Spencer Wells Curved 18cm (7")	34, 110
Forcep Artery Spencer Wells Curved 20cm (8")	34

alphabetical full index

Product	Page
Forcep Artery Spencer Wells Straight 12.5cm (5")	34
Forcep Artery Spencer Wells Straight 15cm (6")	34
Forcep Artery Spencer Wells Straight 18cm (7.2")	34
Forcep Artery Spencer Wells Straight 19cm (7.6")	34, 110
Forcep Artery Spencer Wells Straight 20cm (8")	34
Forcep Artery Spencer Wells Straight 23cm (9.2")	34
Forcep Artery Wilson 19cm (7.5")	85
Forcep Biopsy Eppendorfer Biopsy	52
Forcep Biopsy Lloyd Davies Biopsy 35cm (14")	36
Forcep Biopsy Lloyd Davies Biopsy 50cm (20")	36
Forcep Biopsy Patterson Biopsy 50cm (20")	36
Forcep Biopsy Patterson Biopsy Shaft 35cm (14")	36
Forcep Biopsy Tischler Morgan Biopsy	52
Forcep Clamp Backhaus Towel 11cm (4.4")	35
Forcep Clamp Ball and Socket Towel	35
Forcep Clamp Bonney Polyp and Ovum 24cm (9.5")	52
Forcep Clamp Cross Action Towel	35
Forcep Clamp Denis Browne (Bead) 14cm (5.5")	65
Forcep Clamp Kogan with Ratchet 24cm (9.5")	52
Forcep Clamp Thompson Walker Penile 12cm (4.75")	65
Forcep Clamp Vasectomy Modified Straight Blunt 12.5cm (5")	33, 65
Forcep Clamp Vasectomy Straight 15cm (6")	33, 65
Forcep Clamp Tubing Straight 12.5cm (5")	23, 33
Forcep Clip DeBakey Bulldog Small Curved	23
Forcep Clip DeBakey Bulldog Small Straight	23
Forcep Clip DeBakey Bulldog Medium Straight	23
Forcep Clip DeBakey Bulldog Large Straight	23
Forcep Clip Dieffenbach Bulldog	23
Forcep Clip Michel Applying & Removing 12cm (4.7")	33
Forcep Clip Michel Removing 12cm (4.7")	33
Forcep Clip Raney Applicator and Remover 15cm (6")	81
Forcep Corneal Colibri 0.3mm Teeth, No Platform, 7.5cm (5")	101
Forcep Corneal St Martins 1:2 Teeth and Platform 0.15mm 9cm (3.6")	101
Forcep Corneal St Martins 1:2 Teeth and Platform 0.3mm 9cm (3.6")	101
Forcep Corneal St Martins Non Toothed 9cm (3.6")	101
Forcep Delivery Anderson	75
Forcep Delivery Keilland	75
Forcep Delivery Simpson	75
Forcep Delivery Wrigley	75
Forcep Diathermy Bipolar 11.25cm (4.5") and 3m Lead	36
Forcep Diathermy Bipolar 12.5cm (5") and 3m Lead	101
Forcep Diathermy Bipolar 15cm (6") and 3m Lead	36, 110
Forcep Diathermy Bipolar 20cm (8") and 3m Lead	36

Product	Page
Forcep Dissecting Adson Browne Multi Toothed 12cm (4.8")	85
Forcep Dissecting Adson Fine Non Toothed 12.5cm (5")	35, 81, 85, 123
Forcep Dissecting Adson Fine Toothed 12.5cm (5")	35, 81, 85, 123
Forcep Dissecting Adson Non Toothed 12.5cm (5")	35, 81, 109, 123
Forcep Dissecting Adson Toothed 1:2 Bayonet 11cm (4.5")	35, 81, 109, 123
Forcep Dissecting Adson Toothed 1:2 Bayonet 19cm (7.5")	35, 81, 123
Forcep Dissecting Adson Toothed 12.5cm (5")	35, 81, 123
Forcep Dissecting Adson Toothed Bayonet 19cm (7.5")	81
Forcep Dissecting Bonney Non Toothed 18cm (7.2")	52
Forcep Dissecting Bonney Toothed 18cm (7")	52
Forcep Dissecting Charnley Toothed 14.5cm (5.75")	109
Forcep Dissecting College 15cm (6")	92
Forcep Dissecting Gillies Toothed 15cm (6")	35, 81, 85, 91, 109, 123
Forcep Dissecting Graefe Iris Toothed 1:2 8.75cm (3.5")	100
Forcep Dissecting Iris Non Toothed 10cm (4")	100
Forcep Dissecting Iris Toothed 10cm (4")	100
Forcep Dissecting Iris Toothed Fine 10cm (4")	100
Forcep Dissecting Jean Non Toothed 13.75cm (5.5")	35
Forcep Dissecting Lane Non Toothed 18cm (7")	35, 110
Forcep Dissecting Lane Toothed 1:2 Teeth 12.5cm (5.2")	35
Forcep Dissecting Lane Toothed 1:2 Teeth 15cm (6")	35
Forcep Dissecting Lane Toothed 1:2 Teeth 18cm (7")	35
Forcep Dissector Laparoscopic Maryland Dissector 5 x 170mm (Short)	68
Forcep Dissector Laparoscopic Maryland Dissector 5 x 330mm (Standard)	68
Forcep Dissector Laparoscopic Maryland Dissector 5 x 440mm (Long)	68
Forcep Dissecting McIndoe Non Toothed 15cm (6")	36, 81, 85, 123
Forcep Dissecting McIndoe Toothed 1:2 15cm (6")	36, 81, 85, 123
Forcep Dissecting Moorfields Non Toothed 11.5cm (4.6")	100
Forcep Dissecting Oral Fickling 18cm (7.2")	92
Forcep Dissecting Oral Fickling 1x2 Teeth 18cm (7.2")	92
Forcep Dissecting Ormerod Cup Fine 9cm (3.5")	86
Forcep Dissecting Ramsey Toothed 1:2 Teeth 18cm (7")	36
Forcep Dissecting St Martins 1:2 Teeth and Platform 0.3mm 9cm (3.6")	100
Forcep Dissecting St Martins Non Toothed 9cm (3.6")	100
Forcep Dissecting Struempel Nasal Biopsy 14cm (5.5")	86
Forcep Dissecting T.O.E 12.5cm (5")	36
Forcep Dissecting Tilley Nasal 14cm (5.5")	86
Forcep Dissecting Treves Non Toothed 12.5cm (5")	36

alphabetical full index

Product	Page
Forcep Dissecting Treves Toothed 15cm (6")	36
Forcep Dissecting Treves Toothed Straight 12.5cm (5")	36
Forcep Dissecting Waugh Non Toothed 18cm (7")	36, 81, 86, 123
Forcep Dissecting Waugh Toothed 1:2 20cm (8")	36, 81, 86, 123
Forcep Dressing Cawthorne Crocodile Serrated 10cm (4")	86
Forcep Dressing Denis Browne Tonsil Holding 20cm (8")	86
Forcep Dressing Hartmann Crocodile 14cm (5.5")	86
Forcep Dressing Hartmann Cup 8.5cm (3.3")	86
Forcep Dressing Hartmann Micro Serrated Black 9cm (3.5")	86
Forcep Dressing Lister Sinus 15cm (6")	17, 33, 86, 123
Forcep Dressing Lucs Grasping Oval Jaw 14cm (5.5)	86
Forcep Dressing Tilley Aural 14cm (5.5")	86
Forcep Epilation 8.75cm (3.5")	33, 101
Forcep Epilation 10cm (4")	33, 101
Forcep Extraction Dental 123	92
Forcep Extraction Dental 76N	92
Forcep Fixation Capsulorhexis Forceps Long	101
Forcep Fixation Curved Non Toothed 11cm (4.4")	101
Forcep Fixation Jayle 1:2 Teeth 11.25cm (4.5")	101
Forcep Fixation Serrated Non Toothed Straight 12cm (4.8")	101
Forcep Fixation Straight 1:2 Teeth 11cm (4.4")	101
Forcep Grasper Laparoscopic Grasper With Ratchet Atraumatic Fenestrated Duckbill With Ratchet 5 x 170mm (Short)	68
Forcep Grasper Laparoscopic Grasper With Ratchet Atraumatic Fenestrated Duckbill With Ratchet 5 x 330mm (Standard)	68
Forcep Grasper Laparoscopic Grasper With Ratchet Atraumatic Fenestrated Duckbill With Ratchet 5 x 440mm (Long)	68
Forcep Grasper Laparoscopic Grasper With Ratchet Atraumatic Tapered With Ratchet 5 x 170mm (Short)	68
Forcep Grasper Laparoscopic Grasper With Ratchet Atraumatic Tapered With Ratchet 5 x 330mm (Standard)	68
Forcep Grasper Laparoscopic Grasper With Ratchet Atraumatic Tapered With Ratchet 5 x 440mm (Long)	68
Forcep Grasper Laparoscopic Grasper With Ratchet Babcock With Ratchet 5 x 170mm (Short)	68
Forcep Grasper Laparoscopic Grasper With Ratchet Babcock With Ratchet 5 x 330mm (Standard)	68
Forcep Grasper Laparoscopic Grasper With Ratchet Babcock With Ratchet 5 x 440mm (Long)	68
Forcep Grasper Laparoscopic Grasper With Ratchet Clinch Fenestrated With Ratchet 5 x 170mm (Short)	68
Forcep Grasper Laparoscopic Grasper With Ratchet Clinch Fenestrated With Ratchet 5 x 330mm (Standard)	68
Forcep Grasper Laparoscopic Grasper With Ratchet Clinch Fenestrated With Ratchet 5 x 440mm (Long)	68
Forcep Grasper Laparoscopic Grasper Without Ratchet Atraumatic Fenestrated Duckbill 5 x 170mm (Short)	68

Product	Page
Forcep Grasper Laparoscopic Grasper Without Ratchet Atraumatic Fenestrated Duckbill 5 x 330mm (Standard)	68
Forcep Grasper Laparoscopic Grasper Without Ratchet Atraumatic Fenestrated Duckbill 5 x 440mm (Long)	68
Forcep Grasper Laparoscopic Grasper Without Ratchet Atraumatic Tapered 5 x 170mm (Short)	68
Forcep Grasper Laparoscopic Grasper Without Ratchet Atraumatic Tapered 5 x 330mm (Standard)	68
Forcep Grasper Laparoscopic Grasper Without Ratchet Atraumatic Tapered 5 x 440mm (Long)	68
Forcep Grasper Laparoscopic Grasper Without Ratchet Clinch Fenestrated 5 x 330mm (Short)	68
Forcep Grasper Laparoscopic Grasper Without Ratchet Clinch Fenestrated 5 x 330mm (Standard)	68
Forcep Grasper Laparoscopic Grasper Without Ratchet Clinch Fenestrated 5 x 440mm (Long)	68
Forcep Intubation Magill Adult 25cm (10")	17
Forcep Intubation Magill Child 20cm (8")	17
Forcep Jewellers Angled Smooth No 7 11cm (4.5")	101
Forcep Jewellers Notched Forceps 11cm (4.5")	101
Forcep Jewellers Straight Fine Smooth No 4 11cm (4.5")	101
Forcep Jewellers Straight Medium Smooth No 3 11cm (4.5")	101
Forcep Jewellers Straight Very Fine Smooth No 5 11cm (4.5")	101
Forcep Ligature Lahey 20cm (8")	23
Forcep Ligature Lahey 23cm (9")	23
Forcep Ligature O'Shaughnessy 22cm (8.8")	23
Forcep Sinus Lister 15cm (6")	17, 33, 86, 123
Forcep Splinter Arthurs	36
Forcep Splinter Hunter	36
Forcep Splinter Martins	36
Forcep Rampley Sponge Holding 18cm (7")	34, 52
Forcep Rampley Sponge Holding 24cm (9.5")	34, 52
Forcep Stone Desjardins Grasping Small 23cm (9")	36, 65
Forcep Stone Desjardins Grasping Large 23cm (9")	36, 65
Forcep Suture Moorfields 10cm (4")	101
Forcep Suture Moorfields 12cm (4.8")	101
Forcep Tissue Allis 4:5 Teeth 15cm (6")	35
Forcep Tissue Allis American Pattern 14.5cm (5.75")	123
Forcep Tissue Allis American Pattern 14.5cm (5.75")	35
Forcep Tissue Allis Toothed 2:3 15cm (6")	35, 123
Forcep Tissue Allis Toothed 3:4 15cm (6")	35, 123
Forcep Tissue Allis Toothed 4:5 15cm (6")	123
Forcep Tissue Allis Toothed 4:5 26cm (10")	35
Forcep Tissue Allis Toothed 5:6 20cm (8")	35
Forcep Tissue Babcock 16.5cm (6.5")	35
Forcep Tissue Babcock 20cm (8")	35
Forcep Tissue Bozemann 25cm (10")	52

alphabetical full index

Product	Page
Forcep Tissue DeBakey 18cm (7.2")	24
Forcep Tissue DeBakey 25cm (10")	24
Forcep Tissue DeBakey 30cm (12")	24
Forcep Tissue Duval 23cm (9")	23, 35
Forcep Tissue Green Armytage 20cm (8")	52
Forcep Tissue Judd-Allis 3:4 Teeth 19.5cm (7.75")	35
Forcep Tissue Lanes 17.5cm (7")	35
Forcep Tissue Littlewood Toothed 2:3 18cm (7")	35
Forcep Tissue Maingot Hysterectomy 20cm (8")	52
Forcep Tissue Patterson Grasping 50cm (20")	35
Forcep Tissue Soonawalla 15cm (6")	65
Forcep Tissue Vulsellum Lewer Toothed 1:1 24cm (9.5")	52
Forcep Tissue Vulsellum Teale Uterine Curved Toothed 3:4 18cm (7")	52
Forcep Tissue Vulsellum Teale Uterine Curved Toothed 3:4 23cm (9")	52
Forcep Tissue Vulsellum Teale Uterine Straight Toothed 3:4 23cm (9")	52
Forcep Tying Birks (Tennant) Micro Curved 10cm (4")	101
Forcep Tying Birks (Tennant) Micro Straight 10cm (4")	101
Gag & Blade Boyle Davis Adult and Doughty Blade 10cm (4")	86
Gag & Blade Boyle Davis Adult and Doughty Blade 11.25cm (4.5")	86
Gag & Blade Boyle Davis Adult and Doughty Blade 8.75cm (3.5")	86
Gag & Blade Boyle Davis Child and Doughty Blade 6.25cm (2.5")	86
Gag & Blade Boyle Davis Child and Doughty Blade 7.5cm (3")	86
Gag & Blade Boyle Davis Child and Doughty Blade 8.75cm (3.5")	86
Gag Doyens Adult	86
Gag Doyens Child	86
Handle BP Scalpel BP No 3	33
Handle BP Scalpel BP No 4	33
Handle BP Scalpel BP No 5	33
Handle BP Scalpel BP No 7	33
Handle Mini Blade	33, 86, 99, 114, 123
Handle Scalpel Barron	33, 123
Handle Skin Graft Boards	123
Hook Adson Dissecting Blunt 19.5cm (7.7")	82
Hook Adson Dissecting Sharp 19.5cm (7.7")	82
Hook Cawthorne 14cm (5.5")	87
Hook Dissector Laparoscopic J Blunt High Frequency Electrode 5 x 330mm	68
Hook Dissector Laparoscopic L Blunt High Frequency Electrode 5 x 330mm	68
Hook Dissector Laparoscopic L Cutting High Frequency Electrode 5 x 330mm	68

Product	Page
Hook Emmett Cervical Stay Hook	52
Hook Formby D/E Scoop & Hook	87
Hook Gillies Skin 18cm (7")	38, 88, 102, 124
Hook McIndoe Skin 19cm (7.6")	38, 88, 124
Hook Saunders IUCD Removal 32cm (12.5")	52
Hook St Barts Wax 14cm (5.5")	87
Hook Skin Double (Blunt) 15cm (6")	38, 88, 124
Hook Skin Single (Blunt) 15cm (6")	38, 88, 124
Hook Strabismus Graefe Small 14cm (5.6")	99
Hook Strabismus Graefe Medium 14cm (5.6")	99
Hook Strabismus Graefe Large 14cm (5.6")	99
Hook Strabismus Kilner (Gillies) 14cm (5.6")	99
Hook Strabismus Moorfields 14cm (5.6")	99
Hook Tracheal Double Blunt 18cm (7")	87
Hook Tracheal Single Blunt 18cm (7")	87
Introducers Catheter De Pezzer Curved	65
Laparoscopic Dissectors Maryland Dissector 5 x 170mm (Short)	68
Laparoscopic Dissectors Maryland Dissector 5 x 330mm (Standard)	68
Laparoscopic Dissectors Maryland Dissector 5 x 440mm (Long)	68
Laparoscopic Grasper With Ratchet Atraumatic Fenestrated Duckbill With Ratchet 5 x 170mm (Short)	68
Laparoscopic Grasper With Ratchet Atraumatic Fenestrated Duckbill With Ratchet 5 x 330mm (Standard)	68
Laparoscopic Grasper With Ratchet Atraumatic Fenestrated Duckbill With Ratchet 5 x 440mm (Long)	68
Laparoscopic Grasper With Ratchet Atraumatic Tapered With Ratchet 5 x 170mm (Short)	68
Laparoscopic Grasper With Ratchet Atraumatic Tapered With Ratchet 5 x 330mm (Standard)	68
Laparoscopic Grasper With Ratchet Atraumatic Tapered With Ratchet 5 x 440mm (Long)	68
Laparoscopic Grasper With Ratchet Babcock With Ratchet 5 x 170mm (Short)	68
Laparoscopic Grasper With Ratchet Babcock With Ratchet 5 x 330mm (Standard)	68
Laparoscopic Grasper With Ratchet Babcock With Ratchet 5 x 440mm (Long)	68
Laparoscopic Grasper With Ratchet Clinch Fenestrated With Ratchet 5 x 170mm (Short)	68
Laparoscopic Grasper With Ratchet Clinch Fenestrated With Ratchet 5 x 330mm (Standard)	68
Laparoscopic Grasper With Ratchet Clinch Fenestrated With Ratchet 5 x 440mm (Long)	68
Laparoscopic Grasper Without Ratchet Atraumatic Fenestrated Duckbill 5 x 170mm (Short)	68
Laparoscopic Grasper Without Ratchet Atraumatic Fenestrated Duckbill 5 x 330mm (Standard)	68
Laparoscopic Grasper Without Ratchet Atraumatic Fenestrated Duckbill 5 x 440mm (Long)	68

alphabetical full index

Product	Page
Laparoscopic Grasper Without Ratchet Atraumatic Tapered 5 x 170mm (Short)	68
Laparoscopic Grasper Without Ratchet Atraumatic Tapered 5 x 330mm (Standard)	68
Laparoscopic Grasper Without Ratchet Atraumatic Tapered 5 x 440mm (Long)	68
Laparoscopic Grasper Without Ratchet Clinch Fenestrated 5 x 170mm (Short)	68
Laparoscopic Grasper Without Ratchet Clinch Fenestrated 5 x 330mm (Standard)	68
Laparoscopic Grasper Without Ratchet Clinch Fenestrated 5 x 440mm (Long)	68
Laparoscopic J Hook Blunt High Frequency Electrode 5 x 330mm	68
Laparoscopic L Hook Blunt High Frequency Electrode 5 x 330mm	68
Laparoscopic L Hook Cutting High Frequency Electrode 5 x 330mm	68
Laparoscopic Reposable Handle Without Ratchet	69
Laparoscopic Reposable Scissors Metzenbaum Curved High Insulation 5 x 330mm	69
Laparoscopic Reposable Scissors Metzendaum Curved High Insulation 5 x 170mm (Short)	69
Laparoscopic Reposable Scissors Metzendaum Curved High Insulation 5 x 330mm (Standard)	69
Laparoscopic Reposable Scissors Metzendaum Curved High Insulation 5 x 440mm (Long)	69
Laryngoscope Blade Fibre Optic Macintosh No 1	17
Laryngoscope Blade Fibre Optic Macintosh No 2	17
Laryngoscope Blade Fibre Optic Macintosh No 3	17
Laryngoscope Blade Fibre Optic Macintosh No 4	17
Laryngoscope Blade Fibre Optic Miller No 0	17
Laryngoscope Blade Fibre Optic Miller No 1	17
Laryngoscope Blade Fibre Optic Miller No 2	17
Laryngoscope Blade Fibre Optic Miller No 3	17
Laryngoscope Blade Macintosh No 1 LED	17
Laryngoscope Blade Macintosh No 2 LED	17
Laryngoscope Blade Macintosh No 3 LED	17
Laryngoscope Blade Macintosh No 4 LED	17
Laryngoscope Blade Miller No 0 LED	17
Laryngoscope Blade Miller No 1 LED	17
Laryngoscope Blade Miller No 2 LED	17
Laryngoscope Blade Miller No 3 LED	17
Laryngoscope Handle Fibre Optic	18
Laryngoscope Handle Pen Light	18
Laryngoscope Handle Pen Light Fibre Optic	18
Laryngoscope Handle Standard	18
Mallet Heavy Heath 4.5cm (1.8") (793gms)	110
Mallet Fine Toffee Hammer 19cm (7") (226gms)	110
Marker Scleral Double Ended 3.5mm & 4mm	100

Product	Page
Mirror Laryngeal & Handle 8mm	87
Mirror Laryngeal & Handle 10mm	87
Mirror Laryngeal & Handle 12mm	87
Mirror Laryngeal & Handle 14mm	87
Mirror Laryngeal & Handle 16mm	87
Mirror Laryngeal & Handle 18mm	87
Mirror Laryngeal & Handle 20mm	87
Mirror Laryngeal & Handle 22mm	87
Mirror Laryngeal & Handle 24mm	87
Mirror Laryngeal & Handle 26mm	87
Mirror Post Nasal & Handle 6mm	87
Mirror Post Nasal & Handle 8mm	87
Mirror Post Nasal & Handle 10mm	87
Mirror Dental No 4	36, 93
Mirror Dental No 5	36, 93
Nail Splitter Anvil 13.5cm (5.25")	114
Nail Splitter Thwaites Splitter 13cm (5")	114
Needle Abrams Pleural Biopsy	37
Needle Aneurysm Small 15cm (6")	37
Needle Post Mortem Half (4.4")	37
Needle Holder Converse 12.5cm (5")	24, 37, 124
Needle Holder Combination Foster with R/H Scissors 12.5cm (5")	100
Needle Holder Crile Wood 15cm (6")	24, 37, 110, 124
Needle Holder Crile Wood 20cm (8")	24, 37, 110, 124
Needle Holder DeBakey 18cm (7")	24, 124
Needle Holder Derf 12cm (4.75")	24, 124
Needle Holder Foster with R/H Scissors 12.5cm (5")	24, 37, 124
Needle Holder Gillies 16cm (6.25")	37, 110, 124
Needle Holder Halsey 12.5cm (5")	24, 37, 124
Needle Holder Halsey Tungsten Carbide 15cm (6")	24, 37, 124
Needle Holder Kilner 13.5cm (5.25")	37, 82, 87, 91, 110, 124
Needle Holder Lawrence 15cm (6")	37
Needle Holder Lloyd Davis 25cm (10")	37
Needle Holder Masson 27cm (10.5")	37
Needle Holder Mayo 18cm (7")	37, 110
Needle Holder Mayo 20cm (8")	37
Needle Holder Mayo 22.5cm (9")	37
Needle Holder Mayo Hegar 13cm (5.2")	37
Needle Holder Mayo Hegar 14cm (5.6")	37
Needle Holder Mayo Hegar 16cm (6.4")	37
Needle Holder Mayo Hegar 18cm (7.2")	37, 110
Needle Holder Mayo Hegar Fine 18cm (7")	37

alphabetical full index

Product	Page
Needle Holder Micro Barraquer Curved With Lock 10cm (4")	100
Needle Holder Micro Barraquer Curved Without Lock 10cm (4")	100
Needle Holder Micro Castroviejo Curved Without Lock 14cm (5.5")	100
Needle Holder Micro Castroviejo Curved With Lock 14cm (5.5")	100
Needle Holder Micro Castroviejo Straight Without Lock 14cm (5.5")	100
Needle Holder Micro Castroviejo Straight With Lock 14cm (5.5")	100
Needle Holder Micro Silcock With Lock 14.5cm (5.8")	100
Needle Holder Micro Troutman Curved With Lock 14cm (5.5")	100
Needle Holder Nievert 13cm (5.2")	24, 37, 110, 124
Needle Holder Thomson Walker 20cm (8")	37, 110
Needle Holder Wright 14cm (5.6")	37
Osteotome Nievert Left 20cm (8")	110
Osteotome Nievert Right 20cm (8")	110
Osteotome Nievert Straight 20cm (8")	110
Pin Mayo Instrument 15cm (6")	37
Probes Dental 14W	93
Probes Dental 154	93
Probes Dental 73	93
Probes Dental C.I.T.P.N	93
Probes Dental Caulx	93
Probes Dental Exec 133	93
Probes Dental Flat Plastic	93
Probes Dental H6	93
Probes Dental L9	93
Probes Dental No 152	93
Probes Dental No 5	93
Probes Dental No 6	93
Probes Dental PFI	93
Probes Dental PFI 79	93
Probes Dental Round 8mm	93
Probes Dental Standard	93
Probes Dental Swan Neck	93
Probe Jobson Horne 13cm (5")	87
Probe Jobson Horne Ear 17.5cm (7")	87
Probe Nickel Silver With Eye 12.5cm (5")	37
Probe Nickel Silver With Eye 15cm (6")	37
Probe Wards Double Ended 15cm (6")	37
Probe Lacrimal Double Ended 00/0 Liebreichs	100
Probe Lacrimal Double Ended 0000/000 Liebreichs	100
Probe Lacrimal Double Ended 1/2 Liebreichs	100
Probe Lacrimal Double Ended 3/4 Liebreichs	100

Product	Page
Probe Lacrimal Foster Set of Three 0.6mm, 0.8mm, 1.0mm	100
Proctoscope Gabriel Adult 6.4cm (2.5")	38
Proctoscope Kelly 14cm (5.5")	38
Proctoscope Naunton Morgan Rectal 16cm (6.5")	38
Punch Rongeur Kerrison 18cm Shaft 45° 0.3mm	87, 111
Punch Rongeur Kerrison 18cm Shaft 45° 0.3mm Lite	82
Pusher Ligature Negus Knot Tier 20cm (8")	88
Retractor Minnesota 14cm Long 2cm Wide	93
Retractor Chalazion Desmarres 13mm Oval Plate	102
Retractor Chalazion Desmarres 16mm Oval Plate	102
Retractor Chalazion Desmarres 18mm Oval Plate	102
Retractor Chalazion Desmarres 20mm Oval Plate	102
Retractor Chalazion Desmarres 31mm Oval Plate	102
Retractor Chest Finochietto Large	24
Retractor Chest Finochietto Medium	24
Retractor Chest Finochietto Small	24
Retractor Chest Tudor Edwards Rib Spreader 23cm (9")	24
Retractor Hand Held Allison Lung Retractor 32.5cm (13")	38
Retractor Hand Held Aufricht 16cm (6.25")	88
Retractor Hand Held Czerny 17.5cm (7")	38, 82
Retractor Hand Held Deaver 21.5cm (8.5") Blade 2.5cm (1")	38
Retractor Hand Held Deaver 25cm (10") Blade 2.5cm (1")	38
Retractor Hand Held Deaver 30cm (12") Blade 3.8cm (1.5")	38
Retractor Hand Held Doyen Curved 25cm (12") Blade 10cm (4")	52
Retractor Hand Held Doyen Curved 25cm (12") Blade 8cm (3.5")	52
Retractor Hand Held Gillies Skin Hook 18cm (7")	38, 88, 102, 124
Retractor Hand Held Hook Single (Blunt) 15cm (6")	38, 88, 124
Retractor Hand Held Hook Skin Double (Blunt) 15cm (6")	38, 88, 124
Retractor Hand Held Hurds Pillar 21.5cm (8.5")	88
Retractor Hand Held Kelly 21.5cm (8.5") Blade 3.8cm (1.5")	38
Retractor Hand Held Kilner Alae 8cm (3.25")	88, 91
Retractor Hand Held Kilner Catspaw Double Ended 15cm (6")	38, 88, 91, 111, 124
Retractor Hand Held Lacks Tongue Adult	88
Retractor Hand Held Lacks Tongue Child	88
Retractor Hand Held Langenbeck Large 21.5cm (8.5") Blade 41 x 11mm	38, 82, 88, 111
Retractor Hand Held Langenbeck Medium 21.5cm (8.5") Blade 35 x 11mm	38, 82, 88, 111
Retractor Hand Held Langenbeck Medium Narrow 21.5cm (8.5") Blade 35 x 8mm	38, 82, 88, 111
Retractor Hand Held Langenbeck Medium Wide 21.5cm (8.5") Blade 35 x 15mm	38, 82, 88, 111
Retractor Hand Held Langenbeck Small 21.5cm (8.5") Blade 25 x 6mm	38, 82, 88, 111

alphabetical full index

Product	Page
Retractor Hand Held Langenbeck X-Large Long 21.5cm (8.5") Blade 55 x 11mm	38, 82, 88, 111
Retractor Hand Held Lloyd Davis 29cm (11.5")	38
Retractor Hand Held McIndoes Skin Hook 19cm (7.6")	38, 88, 124
Retractor Hand Held Mollison Pillar 21.5cm (8.5")	88
Retractor Hand Held Morris 24.5cm (9.75") Blade Small 70 x 50mm	38
Retractor Hand Held Morris 24.5cm (9.75") Blade Medium 70 x 65mm	38
Retractor Hand Held Rake 17cm Long 5.5cm Head	38, 111
Retractor Hand Held Waltons Malleable Large 20cm (8") x 50mm (2")	82
Retractor Hand Held Waltons Malleable Medium 20cm (8") x 38mm (1.5")	82
Retractor Hand Held Waltons Malleable Small 20cm (8") x 25mm (1")	82
Retractor Hand Held Watson Orbital Malleable Set of 4	102
Retractor Lid Desmarres Eye (Small) 12mm	102
Retractor Lid Desmarres Eye (Medium) 14mm	102
Retractor Lid Desmarres Eye (Large) 16mm	102
Retractor Lid Jaffe (Pair)	102
Retractor Lid Stallard Small	102
Retractor Lid Stallard Large	102
Retractor Self Retaining Alms Blunt 70mm	111, 124
Retractor Self Retaining Balfour Abdominal	38
Retractor Self Retaining Cone Hemi-Laminectomy 30cm (12")	82
Retractor Self Retaining Gelpi Vulva 16.5cm (6.5")	52
Retractor Self Retaining Jackson Burrows 15.6cm (6.5")	38
Retractor Self Retaining Joll Thyroid 15cm (6")	38
Retractor Self Retaining Kilner Alae 8cm (3.25")	124
Retractor Self Retaining Millin 80 x 25mm	65
Retractor Self Retaining Travers 4:5 Blunt 14cm (5.5")	82
Retractor Self Retaining Travers 4:5 Blunt 20cm (8")	38, 82, 111
Retractor Self Retaining Weitlaner 3:4 Blunt Teeth 14cm (5.5")	39, 82, 88, 111
Retractor Self Retaining Weitlaner 3:4 Sharp Teeth 14cm (5.5")	39, 82, 88, 111
Retractor Self Retaining Weitlaner Small Curved Scalp 14cm (5.5")	82, 111, 124
Retractor Self Retaining West 3:4 Blunt Teeth 14cm (5.5")	39, 88, 111
Saw Gigli Complete with 2 Handles	83, 111
Scissor Enucleation Curved 11cm (4.4")	102
Scissor Iris C.O.F 11.25cm (4.5")	102
Scissor Iris Straight 11cm (4.5")	102
Scissor Iris Straight Micro 11cm (4.5")	102
Scissor McKissock Curved 14.5cm (5.5")	83
Scissor Micro Curved 9cm (3.6")	83
Scissor Stevens Tenotomy Curved Fine 11cm (4.4")	102
Scissor Stevens Tenotomy Straight Fine 11cm (4.4")	89, 102

Product	Page
Scissor Strabismus Curved Blunt 11.5cm (4.6")	40, 102
Scissor Strabismus Straight Blunt 11.5cm (4.6")	40, 102
Scissor Toe Nail Serrated	114
Scissor Westcott Tenotomy Straight Blunt 11cm (4.4")	102
Scissor Westcott Tenotomy Straight Sharp 11cm (4.4")	102
Scissor Angled Potts Angled on Flat 45 Degree 19cm (7.5")	24
Scissor Angled Potts Angled on Side 25 Degree 19cm (7.5")	24
Scissor Angled Potts Angled on Side 45 Degree 19cm (7.5")	24
Scissor Dissecting Aufricht C.O.F 15cm (6")	88
Scissor Dissecting Kilner S/S Curved 11cm (4.5")	88, 124
Scissor Dissecting Kilner S/S Curved 14cm (5.5")	88, 124
Scissor Dissecting Kilner S/S Straight 11cm (4.5")	89, 124
Scissor Dissecting Kilner S/S Straight Fine 11cm (4.5")	89, 124
Scissor Dissecting Lloyd Davies Rectal Curved 25cm (10")	39
Scissor Dissecting Mayo C.O.F 16.5cm (6.5")	39, 53
Scissor Dissecting Mayo C.O.F 19cm (7.5")	39, 53, 112
Scissor Dissecting Mayo C.O.F 22cm (8.5")	53
Scissor Dissecting Mayo C.O.F 23cm (9")	39
Scissor Dissecting Mayo Straight 14cm (5.5")	39, 53
Scissor Dissecting Mayo Straight 15cm (6")	39
Scissor Dissecting Mayo Straight 16.5cm (6.5")	39, 53
Scissor Dissecting Mayo Straight 19cm (7.5")	39, 53, 112
Scissor Dissecting Mayo Straight 22cm (8.5")	53
Scissor Dissecting Mayo Straight 23cm (9")	39
Scissor Dissecting McIndoe C.O.F 12.5cm (5")	39, 124
Scissor Dissecting McIndoe C.O.F 18cm (6")	39
Scissor Dissecting McIndoe C.O.F 18cm (7")	124
Scissor Dissecting McIndoe C.O.F 19cm (7.5")	39, 112, 124
Scissor Dissecting McIndoe C.O.F Tungsten Carbide 23cm (7")	39, 112
Scissor Dissecting Metzenbaum C.O.F 18cm (7")	39
Scissor Dissecting Metzenbaum C.O.F 19cm (7.5")	39, 112
Scissor Dissecting Metzenbaum C.O.F 23cm (9")	39
Scissor Dissecting Metzenbaum C.O.F Tungsten Carbide 18cm (7")	39
Scissor Dissecting Metzenbaum C.O.F Tungsten Carbide 19cm (7.5")	112
Scissor Dissecting Metzenbaum Fine Curved 15cm (6")	39
Scissor Dissecting Metzenbaum Straight 14cm (5.5")	40, 112
Scissor Dissecting Metzenbaum Straight 23cm (9")	40
Scissor Dissecting Sims Uterine Curved 20cm (8")	53
Scissor Dissecting Sims Uterine Straight 20cm (8")	53
Scissor Dissecting Stevens Tenotomy Curved Fine 11cm (4.4")	89, 102
Scissor Dissecting Strabismus Curved Blunt 11.5cm (4.6")	40, 102
Scissor Dissecting Strabismus Straight Blunt 11.5cm (4.6")	40, 102
Scissor Dissecting Wilson Tonsil 20cm (8")	89

alphabetical full index

Product	Page
Scissor Dissecting Woods Tonsil 19cm (7.5")	89
Scissor Dressing B/B Straight 12.5cm (5")	39
Scissor Dressing B/B Straight 20cm (8")	39
Scissor Dressing S/B Straight 12.5cm (5")	39
Scissor Dressing S/S Straight 12.5cm (5")	39
Scissor Dressing S/S Straight 20cm (8")	39
Scissor Dressing Lister Bandage 14cm (5.5")	39, 112
Scissor Dressing Lister Bandage 15cm (6")	39, 112
Scissor Dressing Lister Bandage 18.75cm (7.5")	39, 112
Scissor Episiotomy Bergman 23cm (9")	53
Scissor Episiotomy Braun Stadler 14cm (5.5cm)	53
Scissor Episiotomy Mayo 24cm (9.5")	53
Scissor Episiotomy Vant 24cm (9.5")	53
Scissor Laparoscopic Reposable Handle Without Ratchet	69
Scissor Laparoscopic Reposable Scissors Metzenbaum Curved High Insulation 5 x 330mm	69
Scissor Laparoscopic Reposable Scissors Metzendaum Curved High Insulation 5 x 170mm (Short)	69
Scissor Laparoscopic Reposable Scissors Metzendaum Curved High Insulation 5 x 330mm (Standard)	69
Scissor Laparoscopic Reposable Scissors Metzendaum Curved High Insulation 5 x 440mm (Long)	69
Scissor Micro Castroviejo Corneal Micro 7mm Blade 10cm (4")	103
Scissor Micro Castroviejo Corneal Micro 9mm Blade 10cm (4")	103
Scissor Micro Castroviejo Corneal Micro 12mm Blade 10cm (4")	103
Scissor Micro Vannas Corneal 9mm Blade Curved 8cm (3.2")	103
Scissor Micro Vannas Corneal 9mm Blade Straight 8cm (3.2")	103
Scissor Micro Vannas Corneal 12mm Blade Curved 8cm (3.2")	103
Scissor Micro Vannas Corneal 12mm Blade Straight 8cm (3.2")	103
Scissor Nelson Curved 22.5cm (9")	24
Scissor Nelson Straight 22.5cm (9")	24
Scissor Stitch Littauer 13.5cm (5.25")	89
Scissor Umbilical U.S.A Pattern Straight 10cm (4")	75
Scoop Evacuator	93
Scoop Formby D/E Scoop & Hook	93
Screw Myomectomy 19cm (7")	53
Sizer Orbital Sizer 14mm	103
Sizer Orbital Sizer 16mm	103
Sizer Orbital Sizer 18mm	103
Sizer Orbital Sizer 20mm	103
Sizer Orbital Sizer 22mm	103
Snare Tonsil Eve With Wire	89
Sound Galabin Malleable 29cm (11.5")	53
Sound Horrocks Long 29cm (11.5")	53

Product	Page
Sound Simpson 29cm (11.5")	53
Sound Umbilical 12.5cm (5")	53
Spatula Elliot	103
Specula Aural Tumarkin 4.5cm (1.8")	89
Specula Aural Unispec Oval Set of 5 3.5/4.5/5.5/6.5/7.5mm	89
Specula Aural Unispec Round Set of 5 3/4/5/6/7mm	89
Specula Nasal Thudichum Size 0	89
Specula Nasal Thudichum Size 1	89
Specula Nasal Thudichum Size 2	89
Specula Nasal Thudichum Size 3	89
Specula Nasal Thudichum Size 4	89
Specula Nasal Thudichum Size 5	89
Specula Nasal Thudichum Size 6	89
Specula Nasal Thudichum Size 7	89
Specula Self Retaining Barraquer Open 40 x 26.5mm	103
Specula Self Retaining Barraquer Open 44 x 24mm	103
Specula Self Retaining Kratz Barraquer Open 40 x 33mm Angled	103
Specula Self Retaining Kratz Barraquer Open 40 x 33mm Straight	103
Specula Vaginal Auvard	54
Specula Vaginal Cusco Extra Small	53
Specula Vaginal Cusco Large	54
Specula Vaginal Cusco Large Laser with Smoke Tube Blackened	54
Specula Vaginal Cusco Medium	53
Specula Vaginal Cusco Medium Laser with Smoke Tube Blackened	54
Specula Vaginal Cusco Medium Long	53
Specula Vaginal Cusco Small	53
Specula Vaginal Sims Large Double Ended	54
Specula Vaginal Sims Large Single Ended	54
Specula Vaginal Sims Medium Double Ended	54
Specula Vaginal Sims Medium Single Ended	54
Specula Vaginal Sims Small Double Ended	54
Specula Vaginal Sims Small Single Ended	54
Specula Vaginal Winterton	54
Sucker Frazier 9fg with VC	89
Sucker Frazier 12fg with VC	89
Tonsil Snare Eve With Wire	89

procedure packs

index

Pack	Page
Abdominal Pack	64, 129
Abrams Pleural Biopsy Pack	27, 45
AMD Pack	106
Basic Minor Op Pack	46, 130
Basic Podiatry Pack	118
Biopsy Pack	45
Caesarian Pack	77
Chest Drain Pack	28
Corn Removal Pack	119
Delivery Pack	63
Essential Instrument Delivery Pack	62
Essential Episiotomy Pack	62
IUCD Pack	61
IUCD Pack (With Sound)	61
IUCD Instrument Pack	63
Laparoscopic Cholecystectomy Pack	70
Laparoscopic General Appendix Pack	71
Laparoscopic General Colorectal Pack	72
Laparoscopic General Pack	73
Laparoscopic Gynae Diagnostic Pack	74
Maxillo Facial Pack	127
Minor Op Pack	49, 131
Minor Op Pack	50, 131
Oculoplastic Pack	106
Pacemaker Pack	29
Phaco Pack	105
Podiatry Debridement Pack	118
Rhinoplasty Pack	128
Standard Minor Op Pack	46, 130
Standard Suture Fine Pack	47
Suture Standard Pack	47
Suture Basic Pack	48
Suture Basic Fine Pack	48
Suture Gold Fine Pack	49
Suture Gold Standard Pack	50
Sub Tenon's Pack	105
Vasectomy Pack	67
Zadeks Procedure Pack	119

instrument sets

index

Set	Page
Abdominal Hysterectomy Set	56
Basic Podiatry Set	116
Caesarian Set	76
Cardiovascular Set	26
Chalazion Set	104
D&C Set	59
Dental Examination Set	94
Dental Filling Set	94
Dental Temporary Filling Set	94
Enucleation Strabismus Set	104
Gall Bladder Set	44
General Set	41
General Set Large	42
Gynaecological Basic Set	57
Hysterosalpingogram & Hydrotubation Set	55
I.V Cutdown Set	19
Intestinal Set	43
Laparoscopy Pack	55
Minor Set	43
Nail Surgery Set	116
Nail Surgery Set	117
Perineal Set	55
Podiatry Assistand Set	115
Podiatry Set	115
Prolapse Set	59
Prostaglandin Set	60
Skin Biopsy Set	44, 126
Skin Graft Set	126
Standard Podiatry Set	115
Sub Tenon's Set	104
Suture Set	44, 126
Thoracotomy Basic Set	25
Thyroidectomy Set	95
Tracheostomy Set	19, 95
Tubal Ligation Set	58
Vaginal Suture Set	57
Vaginal Repair & Hysterectomy Set	58
Vasectomy Set	66

numerical full index

Product Code	Page
F200186	27, 45
F200187	28
F200194	29
F300178	77
F300945	64, 129
F300952	70
F300953	71
F300954	72
F300955	73
F300956	74
F500071	105
F500072	105
F500073	106
F500074	106
F700083	127
F700084	128
F810002	54
F810003	54
F810004	53
F810005	53
F810006	53
F810007	53
F810008	54
F810009	54
F810122	63
F810130	62
F810131	62
F881001	19
F882001	25
F882002	26
F883001	61
F883002	61
F883003	51
F883004	45
F883005	46, 130
F883006	46, 130
F883008	47
F883009	47
F883010	48
F883011	48
F883012	67
F883013	49, 131
F883014	50, 131
F883015	50
F883016	49

Product Code	Page
F883018	63
F883019	41
F883020	42
F883021	43
F883022	43
F883023	55
F883024	44
F883025	95
F883026	44, 126
F883027	66
F883028	57
F883029	57
F883030	56
F883031	60
F883032	58
F883033	59
F883034	58
F883035	55
F883036	55
F883037	76
F883038	59
F884001	19, 95
F885002	104
F885003	104
F885004	104
F886001	118
F886002	119
F886003	119
F886004	118
F886005	115
F886007	116
F886010	116
F886014	117
F886015	115
F886016	115
F887001	94
F887002	94
F887003	94
F887005	126
F887006	44, 126
F891001	17
F891002	17
F891003	17
F891004	17
F891005	17

Product Code	Page
F891006	17
F891007	17
F891008	17
F891013	17
F891014	17
F891015	17
F891016	17
F891017	17
F891018	17
F891019	17
F891020	17
F891029	18
F891030	18
F891031	18
F891032	18
F891040	18, 90
F891122	17
F891123	17
F892001	24
F892002	24
F892003	24
F892004	24
F892005	23
F892006	23
F892007	23
F892008	23
F893001	65
F893002	33, 86, 99, 114, 123
F893004	23
F893005	23
F893006	23
F893007	23
F893008	23
F893009	23, 33
F893011	33
F893012	35
F893013	35
F893014	35
F893015	84
F893016	23
F893017	23
F893018	33, 53, 75
F893019	33
F893020	33
F893021	33, 82, 109

numerical full index

Product Code	Page
F893022	85, 111, 123
F893023	85
F893024	33
F893025	34, 81, 123
F893026	34, 81, 123
F893027	34
F893028	34
F893029	34, 110
F893030	34, 110
F893031	34
F893032	23, 34, 81
F893033	23, 34
F893035	23
F893036	34
F893037	34
F893038	34
F893039	34, 110
F893040	34
F893041	34
F893042	34
F893043	34, 110
F893044	34
F893045	34
F893046	34
F893047	34
F893048	36
F893049	36
F893050	36
F893051	65
F893052	36
F893057	35, 81, 123
F893058	35, 81, 123
F893059	35, 81, 109, 123
F893060	35, 81, 109, 123
F893061	36
F893062	36
F893063	24
F893064	35, 81, 85, 91, 109, 123
F893065	100
F893066	100
F893067	35, 110
F893068	35
F893069	35
F893070	35, 110
F893071	36, 81, 85, 123

Product Code	Page
F893073	36
F893074	36
F893077	36
F893080	36
F893081	36
F893082	36, 81, 86, 123
F893083	36, 81, 86, 123
F893084	35
F893085	24
F893086	24
F893087	33, 101
F893088	33, 101
F893089	35
F893090	23
F893091	23, 34
F893092	23, 34
F893093	17, 33, 86, 123
F893094	34, 52
F893095	34, 52
F893096	35, 123
F893097	35, 123
F893098	35, 123
F893099	35, 123
F893100	35
F893101	35
F893102	23, 35
F893103	35
F893104	34
F893105	35
F893106	33, 65
F893107	33
F893108	33
F893109	33, 65
F893110	33, 123
F893111	33
F893112	33
F893114	37
F893115	37
F893116	37
F893117	24, 37, 110, 124
F893118	24, 37, 124
F893119	24, 37, 110, 124
F893120	24, 124
F893121	24, 124
F893122	24, 37, 100, 124

Product Code	Page
F893123	37, 110, 124
F893124	24, 37, 124
F893125	24, 37, 124
F893126	37, 82, 87, 91, 110, 124
F893127	37
F893128	37, 110
F893129	37
F893130	37
F893131	37
F893132	37
F893133	37
F893134	37, 110
F893136	37
F893137	24, 37, 110, 124
F893139	37, 110
F893140	37
F893142	37
F893143	37
F893144	37
F893146	37
F893147	38
F893148	38
F893149	38
F893150	82, 111, 124
F893151	38, 82
F893152	38
F893153	38, 88, 124
F893154	52
F893155	38, 88, 124
F893156	38, 88, 102, 124
F893157	38, 88, 124
F893158	38
F893160	88, 91, 124
F893161	38, 88, 91, 111, 124
F893162	38
F893163	38
F893164	38, 82, 88, 111
F893165	38, 82, 88, 111
F893166	38, 82, 88, 111
F893167	38
F893168	39, 82, 88, 111
F893169	39, 82, 88, 111
F893170	38, 82, 111
F893171	82
F893172	38

numerical full index

Product Code	Page
F893173	39, 88, 111
F893176	39
F893177	39
F893178	39
F893179	39
F893180	39
F893181	89, 124
F893182	88, 124
F893183	89, 124
F893184	88, 124
F893185	39, 112
F893186	39, 112
F893187	39, 112
F893188	89
F893189	39
F893190	39, 53
F893191	39, 53, 112
F893192	39, 53
F893194	39, 53
F893195	39
F893196	39, 53
F893197	39, 53
F893198	39, 124
F893199	39, 124
F893200	39, 112, 124
F893201	40, 112
F893202	39, 112
F893203	39
F893204	83
F893205	24
F893206	24
F893207	24
F893208	24
F893210	114
F893211	39
F893212	75
F893213	39, 53, 112
F893214	89
F893215	39
F893217	54
F893218	40, 54, 75
F893219	40, 54
F893220	51
F893221	51
F893222	51

Product Code	Page
F893223	51
F893224	51
F893225	51
F893227	51
F893228	51
F893229	51
F893230	51
F893231	51
F893232	51
F893233	51
F893234	51
F893235	51
F893236	51
F893237	51
F893238	51
F893239	52
F893240	52
F893241	52
F893242	52
F893243	75
F893244	75
F893245	75
F893246	75
F893247	52
F893248	52
F893251	52
F893252	52
F893253	35
F893254	52
F893255	52
F893256	52
F893257	52
F893258	52
F893259	54
F893260	54
F893261	54
F893262	52
F893263	53
F893264	53
F893265	53
F893266	53
F893267	53
F893268	53
F893269	53
F893270	53

Product Code	Page
F893271	53
F893272	54
F893273	54
F893274	54
F893275	65
F893276	65
F893278	52
F893279	36
F893280	36, 110
F893281	36
F893284	36
F893286	35, 81, 85, 123
F893287	35, 81, 85, 123
F893290	35
F893291	34, 110
F893292	40
F893296	33
F893298	35
F893300	39, 112
F893301	39, 112
F893302	36, 81, 85, 123
F893305	33
F893307	34
F893308	37
F893309	38
F893310	38
F893311	52
F893312	38
F893313	38
F893314	52
F893315	53
F893316	54
F893317	36, 65
F893318	36, 65
F893319	65
F893320	37
F893321	38, 82, 88, 111
F893322	38, 82, 88, 111
F893323	38, 82, 88, 111
F893324	24
F893325	38
F893326	51
F893327	51
F893328	51
F894001	84

numerical full index

Product Code	Page
F894002	84
F894003	84
F894005	92
F894006	92
F894007	92
F894008	92
F894009	92
F894010	84
F894011	84
F894012	84
F894013	84
F894014	84
F894015	84, 87
F894016	84
F894017	84
F894018	84
F894019	84
F894020	84
F894021	84
F894022	88
F894023	88
F894024	17, 85
F894025	17, 85
F894026	85
F894027	92
F894028	92
F894029	92
F894030	92
F894031	92
F894032	92
F894033	92
F894034	85
F894035	85
F894036	85
F894038	92
F894039	92
F894042	36, 93
F894043	36, 93
F894044	93
F894045	93
F894046	93
F894047	93
F894048	93
F894049	93
F894050	93

Product Code	Page
F894051	93
F894052	93
F894053	93
F894054	93
F894055	93
F894056	93
F894057	93
F894058	93
F894059	93
F894060	93
F894061	40, 90, 125
F894062	85
F894063	92
F894065	83, 111
F894066	87, 93
F894067	81
F894068	82, 87, 111
F894069	81
F894070	81
F894071	86
F894072	86
F894073	86
F894074	86
F894075	88
F894076	85
F894077	85
F894078	85
F894079	85
F894080	86
F894081	85
F894082	86
F894083	86
F894084	86
F894085	86
F894086	86
F894087	86
F894088	86
F894089	86
F894090	87
F894091	87
F894092	87
F894094	87
F894095	87
F894096	87
F894097	87

Product Code	Page
F894098	87
F894099	87
F894100	87
F894101	87
F894102	87
F894103	87
F894104	87
F894105	87
F894106	87
F894107	84
F894108	87
F894109	87
F894112	86
F894114	85
F894115	86
F894118	86
F894119	88
F894120	88
F894121	88
F894122	89
F894123	88
F894124	89
F894125	89
F894126	89
F894127	89
F894128	89
F894129	89
F894130	89
F894131	89
F894132	89
F894133	18, 90
F894134	89
F894135	90
F894136	90
F894137	90
F894138	90
F894139	90
F894141	92
F894143	38, 111
F894144	93
F894148	89
F894149	86
F894150	89
F894151	89
F894152	86

numerical full index

Product Code	Page
F894153	89
F894154	89
F894156	92
F894157	92
F894158	81
F894159	81
F894160	82
F894161	83
F894162	82
F894163	82
F894164	82
F894165	82
F894166	82
F894167	82
F895002	99
F895003	99
F895004	99
F895005	99
F895006	81
F895007	99
F895008	99
F895009	99
F895010	99, 102
F895011	99, 102
F895012	99, 102
F895013	99, 102
F895014	100, 101
F895015	100
F895016	101
F895017	101
F895018	101
F895019	100
F895020	99
F895021	101
F895022	101
F895023	101
F895024	100, 101
F895025	101
F895026	101
F895027	101
F895028	100
F895029	101
F895030	99
F895031	99
F895032	100

Product Code	Page
F895033	100
F895034	100
F895035	100
F895036	100
F895037	100
F895038	100
F895039	100
F895040	100
F895041	102
F895042	102
F895043	103
F895044	102
F895045	102
F895046	102
F895047	102
F895048	40, 102
F895049	40, 102
F895050	103
F895051	103
F895052	102
F895053	103
F895054	103
F895055	103
F895056	102
F895057	102
F895058	103
F895061	103
F895062	100
F895064	100
F895066	99
F895067	99, 102
F895068	99
F895069	100
F895070	101
F895071	99
F895072	99
F895073	99
F895074	101
F895075	103
F895076	103
F895077	99
F895078	102
F895079	101
F895080	101
F895081	102

Product Code	Page
F895082	102
F895083	100
F895084	103
F895085	100
F895086	103
F895087	103
F895088	103
F895089	103
F895090	103
F895091	103
F895092	101
F895093	99
F895094	102
F895095	102
F896001	109
F896002	109
F896003	109
F896004	109
F896005	109
F896006	109
F896007	23, 40, 84, 109
F896008	33, 109
F896009	34, 110
F896011	110
F896012	110
F896013	110
F896014	89, 102
F896015	111
F896016	109
F896017	109
F896018	113
F896019	113
F896020	113
F896021	113
F896022	113
F896023	114
F896024	114
F896026	114
F896028	114
F896029	114
F896030	113
F896031	113
F896032	113
F896033	93
F896034	113

numerical full index

Product Code	Page
F896035	113
F896036	113
F896037	113
F896038	113
F896040	113
F896041	113
F896042	113
F896043	113
F896045	109
F896046	124
F896047	109
F896048	109
F896049	109
F896050	109
F896051	109
F896052	110
F896053	110
F896054	111
F896055	111
F896056	109
F896057	109
F896058	111
F896059	109
F897001	101
F897002	101
F897003	111
F897004	123
F898002	68
F898003	68
F898004	69
F898005	69
F898008	69
F898009	69
F898011	68
F898012	68
F898016	69
F898017	68
F898018	68
F898019	68
F898020	68
F898021	68
F898022	68
F898023	68
F898024	68
F898025	68

Product Code	Page
F898026	68
F898027	68
F898028	68
F898029	68
F898030	68
F898031	68
F898032	68
F898033	68
F898034	68
F898035	68
F898036	68
F898037	68
F898038	68
F898039	68
F898040	69
F898041	69
F898042	69
F898043	69
F898044	69
F898045	69
F898046	69
F898047	69
F898048	69
F898100	69
F898101	69



conversion chart

Multiply	By	To Obtain
Millimeters (mm)	0.1	Centimeters (cm)
Millimeters (mm)	0.03937	Inches (in)
Centimeters (cm)	0.3937	Inches (in)
Centimeters (cm)	0.01	Meters (m)
Centimeters (cm)	10	Milimeters (mm)
Inches (in)	2.54	Centimeters (cm)

Millimeters

1 mm 2 mm 3 mm 4 mm 5 mm 6 mm 7 mm 8 mm 9 mm 10 mm

French

2 Fr. 4 Fr. 6 Fr. 8 Fr. 10 Fr. 12 Fr. 14 Fr. 16 Fr. 18 Fr. 20 Fr. 22 Fr.

24 Fr. 26 Fr. 28 Fr. 30 Fr. 32 Fr. 34 Fr. 36 Fr. 38 Fr. 40 Fr.

Inches

1/8" 1/4" 3/8" 1/2" 5/8" 3/4" 7/8" 1"

needle and catheter reference chart

Hypodermic Needle Tubing				French Catheter Sizes			Millimeter/Inch/Centimeter Conversion			
	O.D. In.	I.D. Reg. Wall	I.D. Thin Wall	French Size	O.D. mm	O.D. In.	mm	In.	In.	cm
33	0.008	0.004	0.005	1	0.33	0.013	0.25	0.0098	0.25	0.630
32	0.009	0.004	0.005	2	0.66	0.026	0.50	0.0197	0.50	1.270
31	0.010	0.005	0.006	3	1.00	0.039	0.75	0.0295	0.75	1.900
30	0.012	0.006	0.007	4	1.33	0.053	1	0.0394	1	2.540
29	0.013	0.007	0.008	5	1.66	0.066	2	0.0787	2	5.080
28	0.014	0.007	0.008	6	2.00	0.079	3	0.1181	3	7.620
27	0.016	0.008	0.010	7	2.33	0.092	4	0.1574	4	10.160
26	0.018	0.010	0.012	8	2.66	0.105	5	0.1968	5	12.700
25	0.020	0.010	0.012	9	3.00	0.118	6	0.2362	6	15.240
24	0.022	0.012	0.014	10	3.33	0.131	7	0.2755	7	17.780
23	0.025	0.013	0.015	11	3.67	0.144	8	0.3149	8	20.320
22	0.028	0.016	0.018	12	4.00	0.158	9	0.3543	9	22.860
21	0.032	0.020	0.022	13	4.33	0.170	10	0.3937	10	25.400
20	0.035	0.023	0.025	14	4.67	0.184	11	0.4330	11	27.940
19	0.042	0.027	0.031	15	5.00	0.197	12	0.4724	12	30.480
18	0.050	0.033	0.042	16	5.33	0.210	13	0.5118	13	33.020
17	0.059	0.041	0.046	17	5.67	0.223	14	0.5511	14	35.560
16	0.065	0.047	0.052	18	6.00	0.236	15	0.5905	15	38.100
15	0.072	0.054	0.059	19	6.33	0.249	16	0.6299	16	40.640
14	0.083	0.063	0.071	20	6.67	0.263	17	0.6692	17	43.180
13	0.095	0.071	0.077	21	7.00	0.276	18	0.7086	18	45.720
12	0.109	0.085	0.091	22	7.33	0.288	19	0.7483	19	48.260
11	0.120	0.094	0.100	23	7.67	0.302	20	0.7874	20	50.800
10	0.134	0.106	0.114	24	8.00	0.315	21	0.8267	21	53.340
9	0.148	0.118	0.126	25	8.33	0.328	22	0.8661	22	55.880
8	0.165	0.135	0.143	26	8.67	0.341	23	0.9055	23	58.420
7	0.180	0.150	0.158	27	9.00	0.354	24	0.9448	24	60.960
6	0.203	0.173	0.181	28	9.33	0.367	25	0.9842	25	63.500



Name

Contact Details

Name

Contact Details



 SignatureTM

"Make dealing with Unisurge a standard procedure"



UNISURGE International Limited
Farboud Innovation Park, Formula Drive, Newmarket, CB8 0BF United Kingdom
Telephone: +44 (0) 84 55 21 31 11 Fax: +44 (0) 84 55 21 32 10
Email: info@unisurge.com Website: www.unisurge.com

FOR MORE DETAILS ON OUR COMPREHENSIVE PRODUCT RANGE, UNBEATABLE PRICES AND FREE SAMPLES CALL US NOW ON

084 55 21 31 11