



Your hospital. Your pack. Your way.

See the full range at www.unisurge.com

Leading innovation,
quality and cost efficiency
for hospital supplies



About Unisurge

Established in 1988, Unisurge is a leading manufacturer and supplier of custom procedure packs, dressing packs, medical disposables, and theatre products trusted by healthcare professionals across the UK and beyond.

**2005**

First year trading

2010

Hit £10m turnover

2010

Started building our Newmarket facility

2011

Opened Newmarket Headquarters

2012

Launched uniform, gown and drape range

2014

Opened our Castleford site, our Northern distribution centre

2017

Expanded our Newmarket HQ

2020

Proudly supported the healthcare industry throughout COVID-19

2024

Achieved £50m turnover

2025

Lohmann & Rauscher Group acquired Unisurge International Ltd

Lohmann & Rauscher (L&R), a leading international supplier of medical devices and hygiene products, acquired Unisurge, a UK-based company specialising in surgical supplies and custom procedure trays. This acquisition is part of L&R's strategic expansion efforts to strengthen its market presence in the surgical products sector and enhance its portfolio of high-quality medical solutions. The integration of Unisurge's offerings is expected to provide synergies and innovation opportunities that will benefit healthcare providers and patients alike.

Unisurge also offers a broad range of sterile drapes, gowns and supplementary operating theatre products. The contents of each of our packs have been practitioner specified to meet the requirements of all modern surgical procedures.

Our custom-built Ethylene Oxide sterilisation chambers, some of the largest in Europe, ensure unmatched sterility for all our products, helping to reduce the risk of infections and fluid contamination for patients and medical staff.

Technology and innovation are at the heart of everything we do, alongside strong sustainability practices, including one of the largest solar energy setups in the county.

Supported by significant investment in IT, Unisurge operates a reliable distribution network and maintains its own fleet of delivery vehicles. This allows us to guarantee that sterile packs arrive on time, in perfect condition, and ready for immediate use. Surgical teams can trust that every component they need will be supplied exactly as required.

Our Delivery Fleet

Unisurge's dedicated fleet of vehicles ensures unparalleled control over product distribution, guaranteeing reliable, on-time deliveries. This not only enhances supply chain efficiency but also allows for a truly customer-centric service, tailored to meet precise delivery requirements. With their own fleet, Unisurge eliminates reliance on third-party logistics, ensuring quality, speed, and accountability every step of the way—a game-changer in delivering healthcare solutions where timing is critical.

- Real-time tracking and POD system for complete visibility and rapid issue resolution.
- Full control of the customer journey ensures safe handling, on-time delivery, and minimal damage.
- Modern, low-emission fleet with enhanced branding for visibility.
- Flexible delivery solutions – from a single box to 26 pallets.
- Efficient scheduling and routing for faster response times and a reduced carbon footprint.



Sustainable by Design



76%

of our electricity usage is procured from Green Energy

Electric charge points are installed on site to encourage staff to buy electric vehicles and use Unisurge electricity at a reduced rate.

Our full fleet of company cars are now 100% electric, reflecting our commitment to sustainability



NHS

Evergreen Sustainable NHS Supplier

Alongside this our ETO Sterilisation process is used to heat the warehouse negating the need for any other heating systems.

2040

Is our target to meet net zero emissions



We are on course to meet our new improved Net Zero target, which has been amended to end of 2040 from our original challenge of 2050, this shows our commitment to the environment and protecting our future generations of customers. Our 2025 MSAT score (Modern Slavery Assessment Tool) is an impressive 92%, reinforcing our commitment to all areas of sustainability.

In 2019 we incorporated one of the largest solar installations in West Suffolk on our roof. This is now mirrored at our Castleford site, which together cover 24% of our electricity consumption, with 76% of our remaining usage procured from other forms of Green Energy.

We have installed a full fleet of electric company cars with charging points at both Newmarket and Castleford sites, alongside a 'Cycle to Work' scheme to demonstrate our commitment to staff and the environment.

As a long-term major strategic partner of the National Health Service, we were also invited to undertake the Evergreen Sustainable Supplier Assessment. We proudly completed the assessment successfully and ahead of schedule.

Our commitment to social sustainability includes regular litter picks together with the local police and fire services. Regular activities with the local Newmarket Academy, including competitions, careers days and Q&A sessions and also sponsorship of charity events in the local community.

Your hospital. Your pack. Your way.

Customisable Packs

F1: Anaesthetic Packs


Anaesthetic packs are a critical part of medical procedures, particularly surgeries, providing the instruments and materials needed to administer anaesthesia safely and effectively. They play a key role in managing pain and ensuring patient comfort throughout the procedure.

The contents of each anaesthetic pack are customised to practitioner requirements, and may include items such as fenestrated drapes, manometers, epidural accessories, gowns, and supplementary items like swabs or polyware.

Every pack is clearly labelled with the pack name and destination hospital, often using colour coding as well as a product code and lot number for quick and easy recognition. Designed to support a wide range of procedures, Unisurge anaesthetic packs ensure clinical excellence and patient safety. Delivered via our dedicated fleet of vehicles, they arrive on time and in optimum condition, ready for immediate use.

Customise your procedure packs with a range of items including fenestrated drapes, manometers and epidural accessories

Procedure Pack Example



Your pack.
Your way.



Fully bespoke pack solutions to suit your teams needs



Contact us on
08455 21 31 11
to customise your
packs today



Dedicated team to
help you design your
packs your way



Anaesthetic Packs

Unisurge dedicated
national distribution



Customisable
Packs



"With Unisurge, it's the
service and expertise
we rely on everytime."



On site ETO
sterilisation

Tailor-made
products



Evergreen
Sustainable NHS
Supplier

Precision manufacturing tailored to your product specification

Your pack.
Your way.



Fully bespoke pack
solutions to suit your
teams needs



Contact us on
08455 21 31 11
to customise your
packs today



Dedicated team to
help you design your
packs your way



Procedure Pack Example



Customisable Packs

F2: Cardio Thoracic Packs

Cardio Thoracic Packs are specialised surgical kits designed for use in cardiothoracic procedures, which encompass surgeries involving the heart, lungs, esophagus, and other organs in the thoracic cavity. These packs are pre-assembled to ensure efficiency, sterility, and convenience for surgical teams.

At Unisurge, we understand the pressures surgical teams face. Whether it's a chest aspiration, a complex transplant, a coronary bypass, or a pacemaker insertion, our procedure packs offer a single, tailor-made solution. Cardiac teams can work efficiently and confidently, knowing every theatre pack contains exactly what they specified.

Customise your procedure packs with a range of items
including thoracic catheters, sutures and aortic cannulas

Unisurge dedicated
national distribution



Customisable
Packs



On site ETO
sterilisation

Tailor-made
products



85%

of deliveries are
carried out on our
own vehicles



"Very professional from
start to finish. Flawless
delivery service."

Cardio Thoracic Packs

Your hospital. Your pack. Your way.

Customisable Packs

F3: General Packs

General surgery covers a wide range of procedures for diagnosing, treating, and managing conditions affecting the abdomen, digestive tract, skin, breast, and other soft tissues, often forming the foundation for more specialised surgical disciplines.

A general surgical pack from Unisurge is a customisable, sterile kit designed to support these procedures. Each pack contains essential items such as sterile drapes (often fenestrated), surgical gowns, swabs, disposable instruments (forceps, scissors), suction devices, sutures, needles, polyware, and supplementary items like gloves and dressings, all tailored to practitioner specifications.

Assembled in certified clean rooms, these packs ensure sterility, reduce setup time, and help minimise infection risks, providing surgical teams with confidence and efficiency in the operating theatre.

Procedure Pack Example



Your pack.
Your way.



Fully bespoke pack solutions to suit your teams needs



Contact us on
08455 21 31 11
to customise your
packs today



Dedicated team to
help you design your
packs your way

Customise your procedure packs with a range of items including camera drapes, gloves, wound dressing and suction tubing

General Packs

Unisurge dedicated
national distribution



Customisable
Packs



"We've been so
impressed with the
quality of all our packs."



On site ETO
sterilisation

Greater control
over operations



100⁺

Hospitals supplied
across the UK



Precision manufacturing tailored to your product specification

Your pack.
Your way.



Fully bespoke pack solutions to suit your teams needs



Contact us on
08455 21 31 11
to customise your
packs today



Dedicated team to
help you design your
packs your way



Procedure Pack Example



Customisable Packs

F4: Neurosurgery Packs

Neurosurgery demands the highest level of precision, supporting procedures on the brain, spinal cord, peripheral nerves, and cerebrovascular system. Neurosurgeons are trained to handle complex and delicate procedures that require precision and expertise.

Neurosurgery packs are pre-assembled to practitioner specifications, providing all the necessary tools and sterile materials to support these procedures. They streamline preparation in the operating theatre while maintaining the highest standards of sterility and efficiency.

Customise your procedure packs with a range of items including neuro patties, ligating clips, suture strips and skin marker pens

Unisurge dedicated
national distribution



Customisable
Packs



On site ETO
sterilisation

Modern fleet
of low emission
vehicles



Packs supplied in a
pre-sterile pack



"Having the option to
create bespoke packs
has been a huge help."

Neurosurgery Packs

Your hospital. Your pack. Your way.

Customisable Packs

F5: Ophthalmic Packs

Ophthalmic surgery requires exceptional precision and sterility due to the delicate nature of ocular tissues and the narrow margin for error. Unisurge ophthalmology products are designed to support these demands, providing ophthalmologists with reliable, procedure-specific tools and packs that uphold the highest standards of safety and care.

Using high-quality ophthalmic surgical packs is essential in eye surgery, where precision and sterility are paramount. Each pack is tailored to include sterile tools such as micro-instruments, fenestrated drapes, and irrigation devices, supporting accurate interventions while minimising infection risks. By maintaining a sterile field and providing reliable, procedure-specific equipment, these packs streamline surgical efficiency, reduce complications, and help deliver the best possible outcomes for patients.

Customise your procedure packs with a range of items including cannulas, eye pads and shields, knives and speculums

Procedure Pack Example



Your pack.
Your way.



Fully bespoke pack solutions to suit your teams needs



Contact us on
08455 21 31 11
to customise your
packs today



Dedicated team to
help you design your
packs your way

Ophthalmic Packs

Unisurge dedicated
national distribution



Customisable
Packs



"With Unisurge, it's the
service and expertise
we rely on everytime."



On site ETO
sterilisation

Tailor-made
products



Evergreen
Sustainable NHS
Supplier

Precision manufacturing tailored to your product specification

Your pack.
Your way.



Fully bespoke pack
solutions to suit your
teams needs



Contact us on
08455 21 31 11
to customise your
packs today



Dedicated team to
help you design your
packs your way



Procedure Pack Example



Customisable Packs

F6: Orthopaedic Packs

Orthopaedic procedure packs are designed to bring efficiency and consistency to the operating theatre. Pre-assembled with all essential sterile instruments, implants, and consumables, they streamline preparation and support seamless surgical workflows.

Our orthopaedic packs, part of our Formula range, are tailored to a wide variety of musculoskeletal procedures, from hip and knee replacements to complex hand surgery. Assembled in certified clean rooms, they ensure sterility and compliance with rigorous clinical standards. Each pack can be customised to include procedure-specific components such as sterile drapes, surgical gowns, bone saws, retractors, sutures, cement mixing systems, pulse lavage, stockinettes, and bandages, everything required for precision, safety, and confidence in the theatre.

Customise your procedure packs with a range of items including cement mixing systems, pulse lavage and stockinettes

Unisurge dedicated
national distribution



Customisable
Packs



On site ETO
sterilisation

Tailor-made
products



85%

of deliveries are
carried out on our
own vehicles



"Very professional from
start to finish. Flawless
delivery service."

Orthopaedic Packs

Your hospital. Your pack. Your way.

Customisable Packs


F7: Plastic Packs

Plastic procedure packs are pre-assembled kits containing various medical supplies and instruments required for specific surgical or medical procedures. These packs are designed to improve efficiency and convenience by ensuring that all necessary items are readily available in one package.

Plastic surgery encompasses a range of procedures, from reconstructive operations to cosmetic enhancements, all requiring precision and meticulous technique. Unisurge plastic surgery packs provide tailored, sterile tools, including fine-tipped scissors, specialised retractors, sutures, and drapes, ensuring a controlled environment that minimises infection risk and supports complex procedures. By equipping surgical teams with reliable, procedure-specific instruments, these packs enhance efficiency, reduce complications, and help deliver precise outcomes with minimal scarring.

Customise your procedure packs with a range of items including skin graft blades, nasal bolsters, diathermy pencils and gowns

Procedure Pack Example



Your pack.
Your way.



Fully bespoke pack solutions to suit your teams needs



Contact us on
08455 21 31 11
to customise your
packs today



Dedicated team to
help you design your
packs your way

Plastic Packs

Unisurge dedicated
national distribution



Customisable
Packs



"We've been so
impressed with the
quality of all our packs."



On site ETO
sterilisation

Greater control
over operations



100+

Hospitals supplied
across the UK



Precision manufacturing tailored to your product specification

Your pack.
Your way.



Fully bespoke pack
solutions to suit your
teams needs



Contact us on
08455 21 31 11
to customise your
packs today



Dedicated team to
help you design your
packs your way



Procedure Pack Example



Customisable Packs

F8: Essential Packs

Offering tailored packs to meet the needs of healthcare facilities, these kits provide a streamlined, ready-to-use solution. Each pack contains the instrumentation and essentials required for procedures, ensuring efficiency, safety, and ease of use for medical professionals.

Essential surgical packs are sterile, customisable kits designed for a variety of critical procedures, including general, emergency, and specialty surgeries. They include vital components such as sterile drapes, gowns, instruments, and sutures, all tailored to the specific procedure. Having the right equipment ensures a sterile environment, minimises the risk of surgical site infections, and supports efficient, precise performance in the operating theatre.

Customise your essential packs with a range of items including forceps, swabs and sterile dressings

Unisurge dedicated
national distribution



Customisable
Packs



On site ETO
sterilisation

Modern fleet
of low emission
vehicles



Packs supplied in a
pre-sterile pack




"Having the option to
create bespoke packs
has been a huge help."

Essential Packs

Need help and advice? Talk to our team 08455 21 31 11

UniForm Drapes



We offer a wide range of single-use surgical drapes tailored to the needs of different specialties. Our UniForm drapes are available in two types of base fabric, allowing you to select the style and material that best suits your requirements.

Unisurge offers a comprehensive range of sterile drapes, gowns and supplementary operating theatre products. These include individual drapes, drape packs, incise drapes, camera covers, basin liners and diathermy products. The UniForm draping range offers both specific draping packs as well as individual drapes that provide optimum draping solutions for all surgical procedures at a cost effective price.

The contents of each of our packs have been practitioner specified to meet the requirements of all modern surgical procedures. Our sterile draping packs have been designed to help reduce surgical site infections and fluid contamination whilst offering an efficient means of draping the patient.

We offer lots of variations of drapes, including tri-laminate, SMMS, re-inforced and non re-inforced options.

Sterilised by
Ethylene Oxide



500⁺

Types of drapes to
cover your needs

Unisurge dedicated
national distribution



Precision manufacturing tailored to your product specification

This is just a sample of our product range. For full details of everything we offer, please contact us.

Anaesthetics

Products include:


76 x 76cm with 3 x 5cm adhesive fenestration
 76 x 76cm with 6 x 8cm adhesive fenestration
 76 x 76cm with 6 x 8cm oval fenestration
 76 x 76cm with 15 x 15cm fenestration
 70 x 90cm with 10cm circular fenestration
 110 x 110cm with 55 x 55cm clear panel
 193 x 315cm with 15 x 15cm adhesive fenestration
 193 x 315cm with 15 x 15cm adhesive fenestration
 269 x 158cm with adhesive



Cardio Thoracic

Products include:


Brachial Angiography Drape
 Radial Angiography Drape
 Femoral Angiography Drape
 Paediatric Femoral Angiography Drape
 Pacing Drape Packs
 Cardiovascular Drape
 Perineal Drape
 Perineal Bilateral Under Leg Drape
 Foot Covers
 Huckaback Towels
 Suture Organiser



General

Products include:

Basic Universal Pack
 Re-Inforced Universal Pack
 Re-Inforced Universal Pack With Mayo Cover
 Paediatric Universal Pack
 Lap Incise Drape
 Drape PICC Insertion
 Drape Bariatric
 Spinal Drape
 Utility Drape
 Patient Drape
 Operating Table Cover
 Back Table Cover
 Patient Drape Fenestrated
 Split Sheet



Neurosurgical

Products include:


76 x 76cm with 3 x 5cm adhesive fenestration
 76 x 76cm with 6 x 8cm adhesive fenestration
 76 x 76cm with 6 x 8cm oval fenestration
 76 x 76cm with 10cm circular adhesive fenestration
 76 x 76cm with 15 x 15cm fenestration
 70 x 90cm clear plastic
 90 x 90cm with 50 x 50cm clear panel
 80 x 100cm with 50 x 60cm clear central panel
 Craniotomy Drape
 Split Sheet



Ophthalmics

Products include:


Patient Drape Fenestrated
 Patient Drape
 Universal Ophthalmic Drape
 Full Body Ophthalmic Drape



Orthopaedics

Products include:


U - Bar Pack
 Bar Drape
 Hip Drape
 Split Sheet Impervious
 Knee Drape
 Lower Extremity Drape
 Bilateral Limb Drape
 Upper Extremity Drape
 Universal Extremity Drape
 Knee Arthroscopy Drape with Pouch
 Hand Drape
 Shoulder Sheet
 Spinal Drape
 Impervious Stockinette



Plastics




Products include:

Patient Drape Fenestrated
 ENT Drape
 Head Turban Drape
 Basic ENT Drape Pack
 Basic Universal Pack
 Re-Inforced Universal Pack with Mayo Cover



Need help and advice? Talk to our team 08455 21 31 11

UniForm Surgical Gowns



Unisurge's surgical gowns are made from SMMS fabric, a three-layer non-woven material manufactured from polypropylene. All UniForm surgical gowns fully comply with the requirements of European Standard EN13795.

Our gowns meet all safety requirements for flammability, as polypropylene material naturally recedes from heat sources. They also produce very low levels of lint, helping to reduce the risk of cross-infection in the operating theatre.

Selecting the right surgical gown for each procedure is simple. The colour-coded neck trim clearly identifies the level of protection:

The colour coded neck trim makes it easier to choose the correct level of protection.

White Neck Trim: Standard Surgical Gown – suitable for low-fluid procedures.

Blue Neck Trim: Reinforced Surgical Gown – suitable for medium to high-fluid procedures.

Red Neck Trim: Wet Procedure Gown – suitable for high-fluid cases, e.g., urology.

Low linting SMMS Fabric
& 100% latex free



Full size
range available



Precision manufacturing tailored to your product specification

34

Types of gown to cover your needs



Superior strength, durability and protection



Thermal bonded seams



Non flammable



Description

Description	S	M	L	L - X Long	XL	XL - X Long	XXL	XXL - X Long
Front Length	110	115	123.5	144	133	143	140	150
Back Length	115	120	130	150	140	150	150	160
Sleeve Length	53	56	58.5	58.5	62	62	68.5	68.5
Sleeve Length (Neck to Cuff)	71	76	81	81	89	89	94	94
Cuff Length	9.5	9.5	9.5	9.5	9.5	9.5	9.5	9.5
Neck Length	61	63.5	66	66.04	69	69	76	76
Sweep Width	132	142	155	155	168	168	193	193

Product code and size

Product code and size	Standard SMMS 45g (White Neck Trim)	Standard SMMS 50g (White Neck Trim)	Reinforced Blue (Blue Neck Trim)	Zoned Impervious (Red Neck Trim)	Totally Impervious A-Line (Red Neck Trim)
Petite/Small	F934001	F931001	-	F931017	-
Medium	F934002	F931002	F934010	F931018	-
Large	F934003	F931003	F934011	F931019	-
Large - Extra Long	F934004	F931004	F934012	F931020	F931025
XL	F934005	F931005	F934013	F931021	-
XL - Extra Long	F934006	F931006	F934014	F931022	F931026
XXL	F934007	F931007	F934015	F931023	-
XXL - Extra Long	F934008	F931008	F934016	F931024	-

Need help and advice? Talk to our team 08455 21 31 11

Medical Hollowware



Our surgical bowls have been designed and manufactured by Unisurge International using medical-grade polymer to meet the standards of a medical device.

.....
Made from phthalate-free medical-grade plastic, the bowls are available in 50 ml increments, while our medicine measures/gallipots come in 5 ml increments. The design allows fluid to be poured easily from either corner without compromising shape or durability.



Our Hollowware range consists of: gallipots, bowls, bowl liners, jugs, sieves, which are all available in either sterile or non-sterile form.



Completely recyclable product



15⁺

Types of bowls to cover your needs

Ergonomic design for single-handed use



XL

Lids available for larger sizes

Precision manufacturing tailored to your product specification

Swabs



We offer a wide range of gauze swabs and gauze products, available as pre-packed sterile, double-wrapped, or bulk non-sterile. In addition, we provide a selection of specialist swabs to meet the needs of different procedures.

Our swabs are available in white, blue, and green. They are 100% cotton, meet BP (British Pharmacopoeia) standards, comply with EN14079, and are virtually lint-free.

Our range includes:

Non Woven: Soft, super-absorbent, and low-lint swabs for general use.

Gauze (Plain): 100% absorbent cotton swabs for versatile applications.

Gauze (XRD): 100% absorbent cotton swabs with an interwoven X-ray detectable bar.

Secure Swabs: Sterile, fully absorbent swabs designed specifically for wound care.

Foam Swabs: Including oral polygum swabs for mouth care, foam cubes, and applicators.



Designed to be
virtually lint free



XRD X-ray
detectable

100% absorbent
cotton product



5

Ranges of swabs to
cover your needs

Need help and advice? Talk to our team 08455 21 31 11

Instruments

PC: F810100 - F891122



Our surgical instrument packs contain all the necessary components to provide complete practitioner confidence.

Our pack range includes:

Emergency Packs: Designed for urgent situations, including chest aspiration and suturing.

Outpatients: Covers a wide variety of outpatient procedures, including ENT, ophthalmic, and other specialist areas.

Dental: From basic dental examinations to more complex filling procedures.

Podiatry: From simple nail-cutting sets to more advanced procedures, including Zadeks.

Womens Health: Includes delivery, episiotomy, IUCD, and examination packs.

Minor Op: Covers basic minor surgery, biopsy, and suturing packs.

532

Instruments to cover your needs

Unisurge dedicated national distribution



Precision manufacturing tailored to your product specification

This is just a sample of our product range. For full details of everything we offer, please contact us.

Anaesthetics

Products include:

Fibre-Optic Laryngoscope
Blade Macintosh No.3

Fibre-Optic Laryngoscope
Blade Surgilite Macintosh No.1

Fibre-Optic Laryngoscope
Blade Miller Abbott No.1

Laryngoscope Handle Pen Light
Fibre Optic

Forceps Magills Child

Sucker Yankaeur Rose Tip
28cm (11.2")

Cardiac

Products include:

Cannula Heparin Medium
2.5mm Plain End Sterile

Cannula Heparin Small 2mm
Plain End Sterile

Clamp Debakey Bulldog Large
Straight Sterile

Clamp Debakey Bulldog Medi-
um Straight Sterile

Clamp Debakey Bulldog Small
Curved Sterile

Clamp Debakey Bulldog Small
Straight Sterile

Clamp Dieffenbach Bulldog
Sterile

Clamp Tubing Straight 12.5cm
(5") Sterile

Elevator Periosteal Howarths
Sterile

Forceps Dissecting Debakey
25cm (10")

General

Products include:

Blade Mini Handle Sterile

Clip Towel Backhaus 11cm
(4.4") Sterile

Curette McEwan Sterile

Forceps Artery Dunhill Curved
12.8cm (5.12") Sterile

Forceps Dissecting Iris Toothed
10cm (4") Sterile

Needle Aneurysm Small 15cm
(6") Sterile

Needle Holder Mayo-Hegar
18cm (7.2")

Scissors Mayo COF 19cm (7.5")

Scissors Episiotomy Braun
Stadler 14cm (5.5cm) Sterile

Speculum Sims Medium
Double Ended

Trocár & Cannula Nelson
5.5mm (HRT)

Neurosurgical

Products include:

Dural Elevator Straight Adson
Sterile

Hook Foman Double Sterile

Hudson Brace +14mm
Perforator + 16mm Burr Sterile

Penny Backer Elevator Sterile

Punch Kerrison Rongeur 18cm
Shaft 45' 0.3mm Lite Sterile

Saw Gili Complete + 2 Handles
Sterile

Ophthalmics

Products include:

Curettes Eye Chalazion
Lambert Size D Oval 3.8 x
2.7mm Sterile

Forceps Desmarres 16mm
Wide Sterile

Forceps Dissecting Moorfield
Non Toothed 11.5cm (4.6")
Sterile

Needle Holder Barraquer Micro
No Lock 10cm (4") Sterile

Scissors Castroviejo 12mm
Sterile

Speculae Stallards Lid Large
Sterile

Speculum Eye Self Retaining
Barraquer 44 x 24mm Sterile

Scissors Vannas Corneal 8cm
(3.2") Sterile

Orthopaedics

Products include:

Applicator Tube Gauze Sterile

Cutter Bone Liston 22cm (8.8")
Sterile

Dissector Elevator Watson
Cheyne 190mm Sterile

Nail Splitter Nipper Thwaites
Sterile

Scissors Toe Nail Serrated

Spike Bone Trethowan Sterile

Spoon Volkmann Bone Medium
Size B Sterile

File Nail Medium (Blacks)
13.5cm (5.25") Sterile

Curette Volkmann Medium
Double Ended 14cm (5.5")
Sterile

ENT

Products include:

Speculum Tumarkin 4.5cm
(1.8")

Spoon Evacuator Sterile

Sucker Yankaeur Rose Tip
28cm (11.2")

Suction Magill Sucker Small

Suction Tube Lemperts With
VC & Obturator

Suction Tube Schuknecht

Suction Tube Zoellner with
Stilette 175mm Sterile

Syringe Dental Cartridge Sterile

Trimmer Mitchell Sterile

Trocár Tilley Lichwitz 195mm
8Fg

Tweezer College Sterile

Wards Carver Sterile

Probe Dental Exec133 Sterile

Need help and advice? Talk to our team 08455 21 31 11

Difa® Iodine Drapes

PC: F996648 - F996653



Provides a sterile working field around the operation site and a broad spectrum antimicrobial activity



Available in a variety of sizes to meet your procedural and incision-size needs

The Difa Iodine Incise Drape is a transparent, adhesive polyurethane film that allows the skin to breathe, preventing moisture build-up beneath the drape. It adheres securely to the surrounding skin throughout long surgical procedures, including right up to the wound edge. The adhesive contains PVP-iodine, providing continuous broad-spectrum antimicrobial protection along with superior adhesion to both skin and wound edge.

Key benefits:

Adheres right up to wound edge: Stays in place throughout major operations, including prolonged retraction of wound edges, providing protection from lateral infection.

Transparent non-glare surface: The incision site remains clearly visible without glare.

Thin but extremely strong: Tear-resistant yet can be cleanly cut with a scalpel.


Conformable and elastic: Moulds easily to the operative site. Its elasticity makes it suitable for procedures involving extensive manipulation, such as total hip replacement.

Maintains a sterile operation site: Provides a sterile working surface on which internal organs may be handled safely, minimising contamination risk.

Continuous broad spectrum: The iodine-containing adhesive delivers ongoing broad-spectrum protection.

Precision manufacturing tailored to your product specification

Our EO Sterilisers



Unisurge's ethylene oxide (EO) sterilisers combine cutting-edge technology with rigorous ISO 11135 validation, delivering an exceptional sterility assurance level (SAL 10^{-6}) for maximum safety and reliability. This ensures medical devices are sterile and trustworthy, meeting the highest clinical standards.

.....

We maintain sterility from start to finish, with packaging that complies with ISO 11607 to provide dependable sterile barriers. Our EO sterilisation process is both effective and eco-friendly, reflecting our commitment to sustainability. Each device receives a customised sterilisation approach, tailored to its specific requirements, ensuring sterility throughout its lifecycle while minimising environmental impact.

2 State of the
art EO Sterilisers
on site



Ensures sterility
throughout the
product lifecycle

Eco-friendly
EO sterilisation
processes



Our Services

The right products, in the right hands, at the right price.

At Unisurge, our philosophy is simple: to provide the service and expertise you would expect from a supplier of products critical to the day-to-day running of healthcare departments.

1. A Dedicated Team of Experts

Our support teams are highly trained and knowledgeable, understanding the intricacies of medical supplies. Whether you need product advice, help placing an order, or assistance resolving an issue, our team is just a call or email away.

2. Tailored Solutions for Every Customer

No two clients are the same, and we recognise the importance of personalised care. Our experienced sales team ready to take your order and advise you on delivery times and other Unisurge products - we aim to surpass all the standards of your current service.

3. Reliable and Responsive Communication

We pride ourselves on being accessible and responsive. Clear communication ensures that every question or request is addressed promptly and efficiently.

4. Proactive Problem Solving

We tackle challenges before they escalate. From managing logistics to handling special requests, we ensure your experience remains smooth and stress-free.

5. Unwavering Commitment to Excellence

At Unisurge, customer service is not just a department – it's our philosophy. Every interaction is driven by our commitment to providing you with the best in care, quality, and reliability. We strive to build lasting relationships by consistently putting your needs at the center of everything we do.

6. Deliveries you can rely on

We control the entire customer journey ensuring proper handling of goods and on time deliveries, and over 85% of our national deliveries are carried out on our own low-emission vehicles. Our state of the art tracking and POD system means we can also track deliveries in real time and resolve any issues quickly without compromise.

Fast turnaround
time for orders
and queries



Operating from Newmarket and Castleford, our dedicated fleet guarantees that products arrive on time and in optimal condition



Dedicated
Account Manager
to support you



Our Accreditations



Full Compliance with ISO 13485 and MDR

Our Quality Department ensures every procedure pack is designed, assembled, and sterilised under a certified Quality Management System meeting ISO 13485 standards. All products are fully aligned with the EU MDR (2017/745), guaranteeing complete regulatory compliance and consistent quality.



Rigorous Risk and Safety Management

By applying ISO 14971 risk management principles alongside the MDR's focus on patient safety, our team proactively identifies hazards, implements controls, and ensures the integrity of every pack delivered.



Traceability and Documentation Control

In line with ISO 13485 and MDR requirements, our team maintains complete device history records, ensures robust batch traceability, and uses validated sterilisation processes to guarantee consistency, accountability, and full compliance during audits.



Customer and Patient-Centred Outcomes

Beyond compliance, our team's commitment to continuous improvement ensures every pack is safe, effective, and reliable, giving healthcare professionals confidence and supporting the ultimate goal of patient wellbeing.

For more details on our accreditations and credentials, please visit our website:

unisurge.com/about-us/accreditations



UNISURGE International Limited,
Farboud Innovation Park, Formula Drive,
Newmarket, CB8 0BF, United Kingdom



For more details on our comprehensive product range,
unbeatable prices and free samples call us now on

08455 21 31 11

Contact Us



08455 21 31 11



info@unisurge.com



@unisurge